
Analysis of Chlorzoxazone

In general, a drug ingested is mainly metabolized by cytochrome P450 (CYP) in liver . CYP is an enzyme which converts foreign substances taken into living bodies to the forms easily excreted. Many molecular species of CYP have been identified while 9 species among those (CYP1A2, 2A6, 2B6, 2C8, 2C9, 2C19, 2D6, 2E1, and 3A4) participate in drug metabolism.

An example of the analysis for chlorzoxazone, the index substrate of CYP2E1, and 6-hydroxychlorzoxazone, the metabolite produced from chlorzoxazone, is introduced here. This is the metabolic reaction which can be used as a reference for the HLM (human liver microsome) activity assessment in drug metabolism evaluation.



4

0

4

0

[Sample Preparation Method*]

0.2 mg/mL Human liver microsomes 50 mM Kpi (pH 7.4) 200 μΜ Chlorzoxazone NADPH generating system Total 200 μL

Incubate at 37 °C for 30 min

Stop reaction with 100 µL of acetonitrile

Add 20 μL of 0.1 mM Coumarin as internal standard

Centrifuge at 15,000 rpm for 5 min

Inject 10 µL of supernatant to HPLC

CI Chlorzoxazone
$$C_7H_4CINO_2 = 169.57$$

6-Hydroxychlorzoxazone $C_7H_4CINO_3 = 185.56$

3 DAm 1		6-OH CZX	<u>~</u>	CZX			
Ċ	ָר ' ר	5	10	15 min	20	25	30
	[Chromatogram of Chlorzoxazone (Column: 3 μm)]						

SAMPLE 10 μL of Std. Soln. *	PRESSURE		
PACKING MATERIAL HITACHI LaChrom C18 (5, 3 μm)	TEMPERATURE 35°C		
COLUMN SIZE	SEPARATION METHOD Partition/Adsorption		
4.6 mm I.D. \times 150 mm (P/N : 891-5050) (5 μ m) 4.6 mm I.D. \times 100 mm (P/N : 891-5030) (3 μ m)	DETECTOR UV 295 nm		
ELUENT 50 mmol/L KH_2PO_4 (pH 4.5) / $CH_3CN = 80 / 20$ (v/v)	INSTRUMENTS Chromaster 5110 (Pump), Chromaster 5210 (Autosampler) Chromaster 5310 (Column Oven),		
FLOW RATE 1.0 mL/min	Chromaster 5420 (UV-VIS Detector)		

Sample: Human liver microsome extract (Chlorzoxazon: central muscle relaxant) NOTE This analysis sample was provided by Miki Nakajima, PhD., Division of Pharmaceutical Sciences, Graduate School of Medical Science, Kanazawa University

Bio/Medical Science/Food/Pharmaceutical, Medicine/Pharmaceutical,

Biochemistry, Metabolite, Chlorzoxazone, UV-VIS Spectrometry, LaChrom C18, Chromaster, Partition/Adsorption

High Performance Liquid Chromatograph (HPLC)

Sheet No. LC100021-01