



Analysis of Pb in River Water (Flame Method)



ZA3000

INTRODUCTION : The data obtained by analyzing lead in river water by a ZA3000 series using the flame method is introduced here. ZA3000 series employs the polarized zeeman correction method in the flame method as a background correction, just like a conventional instrument, and thus, very stable baseline can be obtained. As a result, ZA3000 series allows the accurate analysis of a trace amount (0.1 mg/L level) of lead.

INSTRUMENT CONDITIONS		MEASUREMENT PARAMETERS
Element : Pb	Atomizer : STD Burner	Meas. Mode : Working Curve
Instrument : ZA3000	Flame : Air-C ₂ H ₂	Signal Mode : BKG Corrected
Atomization : Flame	Fuel (C ₂ H ₂) : 2.0 L/min	Curve Order : Linear
Wavelength : 283.3 nm	Oxidant (Air) : 160 kPa	Calculation : Integration
Lamp Current : 5.0 mA	15.0 L/min	Time Constant : 2.0 sec
Slit Width : 1.3 nm	Burner Height : 7.5 mm	Calculation Time: 5.0 sec
		Delay Time : 5 sec

NOTE : Lead was added to JSAC 0301, River Water Certified Reference Material, to make concentrations of 0.05 mg/L and 0.1 mg/L.

	CONC(mg/L)	Mean ABS	SD	RSD	REF
STD 1	0.0000	0.00001	0.00004	400.00 %	-0.00034
STD 2	0.1000	0.00057	0.00005	8.77 %	-0.00071
STD 3	0.2000	0.00109	0.00003	2.75 %	-0.00138
STD 4	0.3000	0.00158	0.00001	0.63 %	-0.00194
1	ND	-0.00001	0.00001	- %	-0.00295
1+ Pb 0.05mg/L	0.0492	0.00029	0.00001	3.45 %	-0.00297
1+ Pb 0.1mg/L	0.0960	0.00053	0.00000	0.00 %	-0.00304



KEY WORDS

Environmental Analysis Related, Environmental Water, Clean Water, Environmental Chemistry, River Water, Lead, Pb, Flame, AA, ZA3000, Flame, Environment

Atomic Absorption Photometer
(AA)

Sheet No. AA120005-00