



Analysis of Ni in River Water (Electrothermal Method)

ZA3000




INTRODUCTION: By using the twin injection function, newly installed in the ZA3000 series instruments, nickel in river water was analyzed. By using the specially designed twin cuvette (Pyro D HR), a sample is injected into two different injection ports and thus, a large volume injection is possible. For the measurement, the drying time can be set at the same as that for a conventional cuvette (Pyro C HR). By a conventional instrument, it has been difficult to detect nickel at a level less than 1 µg/L. However, with the twin injection function, the detection without concentration is possible. The analysis result of SLRS-4, a river water certified reference material, was within the range of the certified value, indicating that the accurate analysis is possible.

INSTRUMENT CONDITIONS	MEASUREMENT PARAMETERS	GA AUTOSAMPLER
Element : Ni	Meas. Mode : Working Curve	Sample Volume : 60 µL
Instrument : ZA3000	Signal Mode : BKG Corrected	Addition : Speed : 4
Atomization : GA	Curve Order : Linear	MATRIX MODIFIER
Wavelength : 232.0 nm	Calculation : Peak Height	Matrix Modifier : 100 mg/L Pd+Mg
Lamp Current : 10.0 mA	Time Constant : 0.1 sec	Volume : 20 µL Order : After
Slit Width : 0.2 nm	Temp. Control : ON	
Cuvette : Pyro D HR		

TEMPERATURE PROGRAM					NOTE
Stage	Initial/Final Temperature (°C)	Heating/Keeping (sec)	Gas Flow Rate (mL/min)	Gas	
1 Drying	80 / 140	40 / 0	200	Normal	Pyro D HR, a cuvette specially designed for twin injection was used for the measurement. Sample 1 : SLRS-4 River water Reference Material for Trace Metals
2 Incineration	1000 / 1000	20 / 0	200	Normal	
3 Atomization	2700 / 2700	0 / 3	0	Normal	
4 Cleaning	2800 / 2800	0 / 4	200	Normal	

	CONC (µg/L)	Mean ABS	SD	RSD	REF	ABS
STD 1	0.000	0.0031	0.0005	16.13 %	0.0472	
STD 2	0.250	0.0152	0.0006	3.95 %	0.0446	
STD 3	0.500	0.0277	0.0000	0.00 %	0.0431	
STD 4	1.000	0.0528	0.0005	0.95 %	0.0415	
1	0.613	0.0335	0.0001	0.30 %	0.0422	

Certified value SLRS-4 0.67 ± 0.08 µg/L



KEY WORDS

Environmental Analysis Related, Environmental Water, Clean Water, Environmental Chemistry, River Water, Nickel, Ni, Flameless, Graphite Furnace, AA, ZA3000, GA, Pyro D HR, Environment

Atomic Absorption Photometer (AA)

Sheet No. AA120014-00