



**ZA3000**



### Analysis of As in Drink (Flame Method)

**INTRODUCTION:** The standard for the ingredients in refreshing drinks as specified in the Food Sanitation Act states that arsenic should not be detected. Therefore, for the purpose of quality control, the analysis of arsenic is being performed. ZA3000 series instruments employ the polarized zeeman method as the BKG correction method even for the hydride generation method and it is possible to perform accurate BKG corrections and also to analyze arsenic with an extremely stable baseline. The result of the spike recovery test indicates that arsenic in drinks can be correctly analyzed.

INSTRUMENT CONDITIONS		MEASUREMENT PARAMETERS
Element : As	Atomizer : STD Burner	Meas. Mode : Working Curve
Instrument : ZA3000	Flame : Air-C <sub>2</sub> H <sub>2</sub>	Signal Mode : BKG Corrected
Atomization : Flame	Fuel (C <sub>2</sub> H <sub>2</sub> ) : 1.2 L/min	Curve Order : Linear
Wavelength : 193.7 nm	Oxidant (Air) : 160 kPa	Calculation : Integration
Lamp Current : 12.0 mA		Time Constant : 2.0 sec
Slit Width : 1.3 nm	Burner Height : 10.0 mm	Calculation Time: 5.0 sec
		Delay Time : 5 sec

NOTE : Sample 1 : Deep ocean water, Sample 2: Tea drink

	CONC (μg/L)	Mean ABS	SD	RSD	REF
STD 1	0.00	0.0005	0.0001	20.00 %	0.0175
STD 2	1.00	0.0152	0.0006	3.95 %	0.0196
STD 3	3.00	0.0443	0.0011	2.48 %	0.0216
STD 4	5.00	0.0751	0.0002	0.27 %	0.0282
1	ND	0.0013	0.0001	7.69 %	0.0176
1+ As 1 μg/L	1.07	0.0162	0.0001	0.62 %	0.0220
2	ND	0.0014	0.0001	7.14 %	0.0260
2+ As 1 μg/L	1.12	0.0169	0.0001	0.59 %	0.0296



**KEY WORDS**

Bio/Medical Science/Food/Pharmaceutical,, Food, Food Chemistry, Toxins in Food, Drink, Arsenic, As, Flame, Hydride Generator, AA, ZA3000, HFS-3

Atomic Absorption Photometer (AA)

Sheet No. AA120015-00