



## Analysis of Cr in River Water (Electrothermal Method)


**ZA3000**

**INTRODUCTION:** Chromium is rarely found in natural water. However, as hexavalent chromium is highly toxic, the environmental standard for the chromium is specified (NMT 0.05 mg/L). By using the twin injection function, newly installed to ZA3000 series instruments, total chromium in river water was analyzed. There are two sample injection ports on the twin cuvette (Pyro D HR) and a large volume can be injected while the drying time for the analysis can be set at the same as that for a conventional cuvette (Pyro C HR). The data shown below indicates that chromium in river water at a ng/L level can be detected without concentration.


INSTRUMENT CONDITIONS	MEASUREMENT PARAMETERS	GA AUTOSAMPLER
Element : Cr	Meas. Mode : Working Curve	Sample Volume : 60 $\mu$ L
Instrument : ZA3000	Signal Mode : BKG Corrected	Addition : Speed : 4
Atomization : GA	Curve Order : Linear	<b>MATRIX MODIFIER</b>
Wavelength : 359.3 nm	Calculation : Peak Height	Matrix Modifier : 100 mg/L Pd+Mg
Lamp Current : 7.5 mA	Time Constant : 0.1 sec	Volume : 20 $\mu$ L Order : After
Slit Width : 1.3 nm	Temp. Control : ON	
Cuvette : Pyro D HR		

TEMPERATURE PROGRAM					NOTE
Stage	Initial/Final Temperature ( $^{\circ}$ C)	Heating/Keeping (sec)	Gas Flow Rate (mL/min)	Gas	Pyro D HR, a cuvette specially designed for twin injection, was used for the measurement.  Sample 1: SLRS-4 River Water Reference Material for Trace Metals Sample 2: River water certified reference material (without addition) JSAC 0301-3
1 Drying	80 / 140	40 / 0	200	Normal	
2 Incineration	700 / 700	20 / 0	200	Normal	
3 Atomization	2600 / 2600	0 / 3	0	Normal	
4 Cleaning	2800 / 2800	0 / 4	200	Normal	

	CONC ( $\mu$ g/L)	Mean ABS	SD	RSD	REF	ABS
STD 1	0.000	0.0012	0.0003	25.00 %	0.0280	
STD 2	0.250	0.0269	0.0001	0.37 %	0.0315	
STD 3	0.500	0.0563	0.0021	3.73 %	0.0335	
STD 4	1.000	0.1129	0.0014	1.24 %	0.0215	
1	0.327	0.0369	0.0002	0.54 %	0.0138	
2	0.167	0.0189	0.0004	2.12 %	0.0443	



Certified value of SLRS-4  $0.33 \pm 0.02 \mu\text{g/L}$   
 Certified value of JSAC 0301-3  $0.16 \pm 0.01 \mu\text{g/L}$



**KEY WORDS**

Environmental Analysis Related, Environmental Water, Clean Water, Environmental Chemistry, River Water, Chromium, Cr, Flameless, Graphite Furnace, AA, ZA3000, GA, Pyro D HR, Environment

Atomic Absorption Photometer (AA)

Sheet No. AA120020-00