



Analysis of In in Work Environment Air (Electrothermal Method)

ZA3000


INTRODUCTION: For a workplace where indium tin oxide (ITO) is being handled, a notice including the health hazard preventive measures and technical guidelines for laborers has been issued by the Ministry of Health, Labour and Welfare (Notification No. 1222-2 by the Labour Standards Bureau, Industrial Safety and Health Dept.). The concentration at which the work environment needs to be improved is 0.01 mgIn/m³ (In 10 µg/L in sample solution) and the allowable concentration for the long-term exposure is 3 × 10⁻⁴ mgIn/m³ (In 0.3 µg/L in sample solution). In the polarized zeeman electrothermal atomic absorption spectrometry, the twin injection allows a large volume injection. Thus, it is possible to quantitatively analyze In 0.3 µg/L, the level corresponds to the allowable concentration, without in-furnace concentration. The addition recovery for the extract from the work environment air was also found to be good.

ANALYTICAL CONDITIONS	MEASUREMENT PARAMETERS	GA AUTOSAMPLER
Element : In	Meas. Mode : Working Curve	Sample Volume : 60 µL
Instrument : ZA3000	Signal Mode : BKG Corrected	Addition : Speed : 2
Atomization : GA	Curve Order : Linear	MATRIX MODIFIER
Wavelength : 325.6 nm	Calculation : Peak Height	Matrix Modifier : 100 mg/L Pd+Mg
Lamp Current : 10.0 mA	Time Constant : 0.1 sec	Volume : 20 µL Order : After
Slit Width : 1.3 nm	Temp. Control : ON	
Cuvette : Pyro D HR		

TEMPERATURE PROGRAM					NOTE
Stage	Initial/Final Temperature (°C)	Heating/Keeping (sec)	Gas Flow Rate (mL/min)	Gas	Pyro D HR, a cuvette specially designed for twin injection, was used for the measurement. [Preparation Method] Inhalable dust was collected by the filter collection method at 2 L/min for 20 minutes. The dust on the filter was dissolved in 15 mL of the mixed acid for extraction (water: nitric acid: hydrochloric acid = 4 : 1: 3) and then, the volume was made up to 40 mL.
1 Drying	80 / 110	40 / 0	200	Normal	
	110 / 300	40 / 0	200	Normal	
2 Incineration	1100 / 1100	20 / 0	200	Normal	
3 Atomization	2800 / 2800	0 / 3	0	Normal	
4 Cleaning	2800 / 2800	0 / 4	200	Normal	

	CONC (µg/L)	Mean ABS	SD	RSD	REF	ABS
STD 1	0.00	0.0000	0.0001	- %	0.0184	
STD 2	0.25	0.0033	0.0000	0.00 %	0.0206	
STD 3	0.50	0.0072	0.0000	0.00 %	0.0192	
STD 4	1.00	0.0149	0.0001	0.67 %	0.0206	
1	0.01	0.0000	0.0000	- %	0.0182	
2	0.01	0.0000	0.0000	- %	0.0178	
3	0.30	0.0043	0.0002	4.65 %	0.0180	

Recovery rate (No.3-No.2) : $(0.30-0.01)/0.3 \times 100 = 96.7 \%$



KEY WORDS

Environmental Analysis Related, Air, Environmental Chemistry, Air Contamination, Work Environment, Handling of ITO, etc., Indium, In, Flameless, Graphite Furnace, Twin Injection, AA, ZA3000, GA, Pyro D HR, Environment

Atomic Absorption Photometer (AA)

Sheet No. AA120026-00