



### ZA3000

## Analysis of Pb in Food Additives (Flame Method)

**INTRODUCTION:** The test method for purity in food additives is specified in the Japanese Standards of Food Additives. The Heavy Metal Limit Test in which the color reaction of lead is expressed in a numerical value is employed for the test of heavy metals. The current component specification is 10 – 20 µg/g, calculated as lead. However, to harmonize with the international specification by JECFA, it is predicted that the specification will become applicable only to lead and the value will be lowered. ZA3000 series instruments employ the polarized Zeeman method for BKG corrections even for the flame method. The accurate BKG corrections and stable baseline allow the analysis of 2 µg/g of lead in food additives.

ANALYTICAL CONDITIONS		MEASUREMENT PARAMETERS	
Element : Pb	Atomizer : STD Burner	Meas. Mode : Working Curve	
Instrument : ZA3000	Flame : Air-C <sub>2</sub> H <sub>2</sub>	Signal Mode : BKG Corrected	
Atomization : Flame	Fuel (C <sub>2</sub> H <sub>2</sub> ) : 2.0 L/min	Curve Order : Linear	
Wavelength : 283.3 nm	Oxidant (Air) : 160 kPa	Calculation : Integration	
Lamp Current : 7.5 mA		Time Constant : 5.0 sec	
Slit Width : 1.3 nm	Burner Height : 7.5 mm	Calculation Time : 5.0 sec	
		Delay Time : 5 sec	



NOTE : 1: L-Ascorbic acid, 2: Potassium nitrate

1 g of each sample was weighed out and the volume was made up to 50 mL with purified water.

Addition recovery rate: (Sample 1: 90.5 %), (Sample 2: 92.0 %)

\* The addition of 2 µg/g of Pb to the sample will result in 0.04 mg/L of Pb in the solution.

	CONC (mg/L)	Mean ABS	SD	RSD	REF	ABS
STD 1	0.0000	0.00002	0.00001	50.00 %	-0.00085	
STD 2	0.0500	0.00039	0.00000	0.00 %	-0.00123	
STD 3	0.1000	0.00075	0.00003	4.00 %	-0.00157	
STD 4	0.2000	0.00154	0.00001	0.65 %	-0.00198	
1	ND	0.00001	0.00001	00.00 %	-0.00242	
1+ Pb 2 µg/g	0.0362	0.00029	0.00001	3.45 %	-0.00271	
	0.0362 (mg/L)×50 (mL)/1(g) = 1.81(µg/g)					
2	ND	0.00001	0.00001	00.00 %	-0.00258	
2+ Pb 2 µg/g	0.0368	0.00029	0.00000	0.00 %	-0.00302	
	0.0368 (mg/L)×50 (mL)/1 (g) = 1.84 (µg/g)					



#### KEY WORDS

Bio/Medical Science/Food/Pharmaceutical, Food, Food Chemistry, Food Additive, L-Ascorbic Acid, Potassium Nitrate, Lead, Pb, Flame, AA, ZA3000

Atomic Absorption Photometer (AA)

Sheet No. AA120032-00