

Absorption Spectrum of Green Dopant Materials for OLED Display (Coumarin 6)

INTRODUCTION

Coumarin 6 is one of the green-emitting phosphorescence used for organic EL (Organic Electro Luminescence) emitting layers.

The fluorescence properties of Coumarin 6 were confirmed by using F-7000.


Spectral corrections are necessary to obtain accurate emission properties. By using F-7000 with the substandard light source and R928F photomultiplier, the spectral corrections over a broad wavelength range, from the UV to visible region (200 - 800 nm), are possible.

SAMPLE

Sample : Coumarin 6
 $C_{20}H_{18}N_2O_2S$ mol.wt. 350.44
 (Sigma-Aldrich, Tokyo, Japan)

Solven : Ethanol

Concentration : 500 μ g/L



ANALYSIS CONDITIONS


WAVELENGTH (nm)

Instrument : U-3900H

Scan speed : 300 nm/min

Slit : 2 nm

1. 459 nm



KEY WORDS

Material-Processing Material Related,
 Other Material-Processing Material Related,
 Green Dopant Materials for OLED Display, Coumarin 6, Absorption
 Spectrum, Green Dopant Materials, OLED, UV, U-3900H, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL100007-01


Excitation Spectrum of Green Dopant Materials for OLED Display (Coumarin 6)


INTRODUCTION

Coumarin 6 is one of the green-emitting phosphorescence used for organic EL (Organic Electro Luminescence) emitting layers.

The fluorescence properties of Coumarin 6 were confirmed by using F-7000.

Spectral corrections are necessary to obtain accurate emission properties. By using F-7000 with the substandard light source and R928F photomultiplier, the spectral corrections over a broad wavelength range, from the UV to visible region (200 - 800 nm), are possible.

SAMPLE	ACCESSORY
<p>Sample : Coumarin 6 (Coumarin 6) $C_{20}H_{18}N_2O_2S$ mol.wt. 350.44 (Sigma-Aldrich, Tokyo, Japan)</p> <p>Solvent : Ethanol</p> <p>Concentration : 10 μg/L</p>	<p>Substandard Light Source (P/N : 5J0-0135/5J0-0136)</p>
	
Analysis Conditions	Wavelength (nm)
<p>Instrument : F-7000</p> <p>Fluorescence wavelength : 500 nm Response : Automatic</p> <p>Slit on excitation side : 5 nm Detector : R928F</p> <p>Slit on fluorescence side : 5 nm Photomultiplier voltage : 400 V</p> <p>Scan speed : 1200 nm/min</p>	<p>1. 450 nm</p>



[With Spectral Corrections]

KEY WORDS

Material-Processing Material Related,
 Other Material-Processing Material Related,
 Green Dopant Materials for OLED Display, Coumarin 6, Excitation Spectrum,
 Green Dopant Materials, OLED, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL100007-02


Fluorescence Spectrum of Green Dopant Materials for OLED Display (Coumarin 6)


INTRODUCTION

Coumarin 6 is one of the green-emitting phosphorescence used for organic EL (Organic Electro Luminescence) emitting layers.

The fluorescence properties of Coumarin 6 were confirmed by using F-7000.

Spectral corrections are necessary to obtain accurate emission properties. By using F-7000 with the substandard light source and R928F photomultiplier, the spectral corrections over a broad wavelength range, from the UV to visible region (200 - 800 nm), are possible.

SAMPLE	ACCESSORY
<p>Sample : Coumarin 6 $C_{20}H_{18}N_2O_2S$ mol.wt. 350.44 (Sigma-Aldrich, Tokyo, Japan)</p> <p>Solvent : Ethanol</p> <p>Concentration : 10 μg/L</p>	<p>Substandard Light Source (P/N : 5J0-0135/5J0-0136)</p>
	
ANALYSIS CONDITIONS	WAVELENGTH (nm)
<p>Instrument : F-7000</p> <p>Excitation wavelength : 450 nm</p> <p>Slit on excitation side : 5 nm</p> <p>Slit on fluorescence side : 5 nm</p> <p>Scan speed : 1200 nm/min</p>	<p>1. 500 nm</p>
<p>Response : Automatic</p> <p>Detector : R928F</p> <p>Photomultiplier voltage : 400 V</p>	



[With Spectral Corrections]

KEY WORDS

Material-Processing Material Related,
 Other Material-Processing Material Related,
 Green Dopant Materials for OLED Display, Coumarin 6,
 Fluorescence Spectrum, Green Dopant Materials, OLED, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL100007-03

3D Fluorescence Spectrum of Green Dopant Materials for OLED Display (Coumarin 6)


INTRODUCTION

Coumarin 6 is one of the green-emitting phosphorescence used for organic EL (Organic Electro Luminescence) emitting layers.

The fluorescenc properties of Coumarin 6 were confirmed by using F-7000.

Spectral corrections are necessary to obtain accurate emission properties. By using F-7000 with the substandard light source and R928F photomultiplier, the spectral corrections over a broad wavelength range, from the UV to visible region (200 - 800 nm), are possible.



SAMPLE		ACCESSORY
Sample	: Coumarin 6 C ₂₀ H ₁₈ N ₂ O ₂ S mol.wt. 350.44 (Sigma-Aldrich, Tokyo, Japan)	Substandard Light Source (P/N : 5J0-0135/5J0-0136)
Solvent	: Ethanol	
Concentration	: 10 µg/L	



ANALYSIS CONDITIONS

Instrument	: F-7000	Response	: Automatic	Full scale	: 100
Excitation wavelength	: 450 nm	Detector	: R928F	Contour line interval	: 2
Slit on excitation side	: 5 nm	Photomultiplier	: 400 V		
Slit on fluorescence side	: 5 nm				
Scan speed	: 60000 nm/min				

Excitation Wavelength (nm)



[With Spectral Corrections]

KEY WORDS

Material-Processing Material Related, Other Material-Processing Material Related, Green Dopant Materials for OLED Display, Coumarin 6, 3D Fluorescence Spectrum, 3D, Green Dopant Materials, OLED, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL100007-04

Calibration Curve of Green Dopant Materials for OLED Display (Coumarin 6)


INTRODUCTION

Coumarin 6 is one of the green-emitting phosphorescence used for organic EL (Organic Electro Luminescence) emitting layers.

The fluorescence properties of Coumarin 6 were confirmed by using F-7000.



Spectral corrections are necessary to obtain accurate emission properties. By using F-7000 with the substandard light source and R928F photomultiplier, the spectral corrections over a broad wavelength range, from the UV to visible region (200 - 800 nm), are possible.

SAMPLE	ACCESSORY
<p>Sample : Coumarin 6 $C_{20}H_{18}N_2O_2S$ mol.wt. 350.44 (Sigma-Aldrich, Tokyo, Japan)</p> <p>Solvent : Ethanol</p> <p>Concentration : 0 to 100 $\mu\text{g/L}$</p>	<p>Substandard Light Source (P/N : 5J0-0135/5J0-0136)</p>



ANALYSIS CONDITIONS

Instrument	: F-7000		
Excitation wavelength	: 450 nm	Response	: Automatic Calibration curve
Slit on excitation side	: 5 nm	Detector	: R928F Fluorescence wavelength : 500 nm
Slit on fluorescence side	: 5 nm	Photomultiplier	: 400 V
Scan speed	: 1200 nm/min		



[With Spectral Corrections]

KEY WORDS

Material-Processing Material Related,
 Other Material-Processing Material Related,
 Green Dopant Materials for OLED Display, Coumarin 6, Calibration Curve,
 Green Dopant Materials, OLED, Working Curve, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL100007-05