




Analysis of Contaminants in Aspirin Tablet using 2-wavelength Measurement

The UV/UV-VIS detector of Chromaster allows the two-wavelength simultaneous analysis by real-time wavelength switching. An analysis example of the aspirin decomposition products is introduced here. In the US pharmacopoeia and Japanese pharmacopoeia (dissolution test), the analytical wavelength for aspirin is specified as 270 nm. However, by also including the measurement at the wavelength of 205 nm, the presence of decomposition product such as salicylic acid or additives can be confirmed.

Chromaster allows a short data loading time of 400 ms when switching the wavelength and thus, peaks can be accurately detected without missing the data points, contributing to the improved reliability of the analysis values.




Aspirin
C₉H₈O₄ = 180.16



Salicylic Acid
C₇H₆O₃ = 138.12

Enlarged Chromatogram of Aspirin (270 nm)



[Sample Preparation *]

- Sample 1 tablet
- ← Purified water 1000 mL
- Sonication 3 min
- Filtration Pore Size 0.45 μm
- ↓
- Sample for injection

SAMPLE	20 μL of Sample Soln. *	PRESSURE	
PACKING MATERIAL	HITACHI LaChrom C18 (5 μm)	TEMPERATURE	40 °C
COLUMN SIZE	4.6 mm I.D. × 150 mm (P/N : 891-5050)	SEPARATION METHOD	Partition/Adsorption
ELUENT	20 mmol/L CH ₃ COONH ₄ / CH ₃ OH = 80 / 20 (v/v)	DETECTOR	UV 205, 270 nm
FLOW RATE		1.0 mL/min	INSTRUMENTS
KEY WORDS	Bio/Medical Science/Food/Pharmaceutical, Medicine/Pharmaceutical, Medical Science/Pharmaceutical Science, Aspirin, Salicylic Acid, 2-wavelength Detection, UV-VIS Spectrometry, Health, Aspirin, Acetyl Salicylic Acid, Salicylic Acid, 2-wavelength Measurement, Chromaster, LaChrom C18, Partition/Adsorption		High Performance Liquid Chromatograph (HPLC)
			Sheet No. LC100026-01