

## Measurement of Prism

### INTRODUCTION

Prism is used as a part of various products such as video cameras, projectors, and liquid crystal televisions. By installing a prism measurement accessory into U-4100 spectrophotometer, the transmission spectrum and reflectance spectrum of a prism can be measured. This time, two types of prisms, a beam splitter and a polarized beam splitter, were analyzed. It was found that the beam splitter split the transmittance and reflectance of the incident light at 1:1 regardless of the polarization condition while the transmittance and reflectance by the polarizing beam splitter depended on the polarizing element.

### SAMPLE

Sample : Beam splitter

#### INSTRUMENT CONDITIONS

Instrument : U-4100 spectrophotometer (solid sample measurement system)

Analysis wavelength range : 485 - 550 nm

Sampling interval : 1 nm

[UV/VIS]

Scan speed : 300 nm/min



Slit : 2 nm

#### OTHER NECESSARIES


Polarizer holder  
(P/N : 132-0325)

Glan Taylor polarizing prism  
MGTYB 20

Prism measurement  
accessory  
(P/N : 134-0110)



Transmission Spectrum of Beam Splitter



Reflectance Spectrum of Beam Splitter

#### KEY WORDS

Material/Processing Material Related, Glass/Ceramics, Industrial Chemistry, Prism, Transmittance, Reflectance, Spectrophotometer, Beam Splitter, U-4100, Material

Spectrophotometer (UV)

Sheet No. UV100019-01

## Measurement of Prism

### INTRODUCTION

Prism is used as a part of various products such as video cameras, projectors, and liquid crystal televisions. By installing a prism measurement accessory into U-4100 spectrophotometer, the transmission spectrum and reflectance spectrum of a prism can be measured. This time, two types of prisms, a beam splitter and a polarized beam splitter, were analyzed. It was found that the beam splitter split the transmittance and reflectance of the incident light at 1:1 regardless of the polarization condition while the transmittance and reflectance by the polarizing beam splitter depended on the polarizing element.

### SAMPLE

Sample : Polarizing beam splitter

#### INSTRUMENT CONDITIONS

Instrument : U-4100 spectrophotometer (solid sample measurement system)

Analysis wavelength range : 550 - 800 nm

Sampling interval : 1 nm

[UV/VIS]

Scan speed : 300 nm/min



Slit : 2 nm

#### OTHER NECESSARIES


Polarizer holder  
(P/N : 132-0325)

Glan Taylor polarizing prism  
MGTYB 20

Prism measurement  
accessory (P/N : 134-0110)



Transmission Spectrum of Polarizing Beam Splitter



Reflectance Spectrum of Polarizing Beam Splitter

#### KEY WORDS

Material/Processing Material Related, Glass/Ceramics, Industrial Chemistry, Prism, Transmittance, Reflectance, Spectrophotometer, Polarizing Beam Splitter, U-4100, Material

Spectrophotometer (UV)

Sheet No. UV100019-02