

Reflectance Spectrum of SiO₂ Film by Small Size 5° Specular Reflectance Accessory

INTRODUCTION

The small size 5° specular reflectance accessory allows the measurement of a reflectance spectrum for a micro sample of about $\phi 5$ mm at an incident angle of 5° for the measurement part. This time, the reflectance spectrum of the SiO₂ film on a silicon substrate used for semi-conductors, etc. was measured. As a result, a reflectance spectrum with an interference fringe was obtained. By using this accessory, the low noise levels of micro samples can be easily measured.

SAMPLE

Sample : SiO₂ film of silicon substrate
10 mm(H) × 10 mm(W) × 1 mm(T)

INSTRUMENT CONDITIONS

Instrument : U-4100 spectrophotometer (solid sample measurement system)

Measurement wavelength region : 240 - 2500 nm

Sampling interval : 1 nm

[UV/VIS]

Scan speed : 300 nm/min

Slit : 8 nm

[NIR]

Scan speed : 750 nm/min

Slit : Automatic control

PbS sensitivity : 2

ACCESSORY

Small size 5° specular reflectance accessory
(P/N : 134-0103)




Figure 1 : Reflectance Spectrum of SiO₂ Film on Silicon Substrate

KEY WORDS

Electronics/Semiconductor Related, Semiconductor,

Silicon Substrate, Reflectance Spectrum, Absolute Reflectance,

Small Size 5° Specular Reflectance Accessory, Spectrophotometer, U-4100

Spectrophotometer (UV)

Sheet No. UV110004-01

Measurement of Film Thickness of SiO₂ Film (Observation of Scanning Electron Microscope)

INTRODUCTION

Based on the obtained spectrum, the film thickness was calculated by using the film thickness calculation, the function enabled by the optional application measurement package. By entering the refractive index and incident angle, the film thickness can be easily calculated. As a result, the film thickness was found to be 1 μm. In addition, by using S-3400N scanning electron microscope which allows the clear confirmation of the microscopic structure of a sample, the cross section of the sample was observed. The film thickness of 1 μm was confirmed and thus, the validity of the measurement by the spectrophotometer was confirmed. By measuring the reflectance spectrum, its film thickness can be calculated easily.

SAMPLE

Sample : SiO₂ film on silicon substrate
10 mm(H) × 10 mm(W) × 1 mm(T)

OBSERVATION CONDITIONS

Instrument : Hitachi S-3400N Scanning Electron Microscope
Observation condition : high vacuum reflection electron image
Accelerated voltage : 5 kV
Observation magnification : × 20,000




Figure 2 : Reflectance Spectrum of SiO₂ Film on Silicon Substrate
(Enlarged Spectrum of Figure 1)

* Refer to UV110004-01 for the instrument conditions.




Figure 3 : SEM Image of Sample Cross Section

KEY WORDS

Electronics/Semiconductor Related, Semiconductor, Silicon Substrate, Reflectance Spectrum, Absolute Reflectance, Small Size 5° Specular Reflectance Accessory, Spectrophotometer, Hitachi Scanning Electron Microscope, U-4100, S-3400N

Spectrophotometer (UV)

Sheet No. UV110004-02