




UH5300


Analysis of PO₄³⁻ in Drainage Water with UH5300

ANALYSIS METHOD	PREPARATION METHOD FOR AMMONIUM MOLYBDATE-ASCORBIC ACID MIXED SOLUTION
<p>Sample 25 mL ← Reagent addition (*1) 2 mL</p> <p>Stir ← Let stand for 25 min</p> <p>Solution for measurement</p> <p>Temperature: room temperature</p> <p>*1 Ammonium molybdate-ascorbic acid mixed solution</p>	<p>Solution 1: Ammonium molybdate solution Weigh 6 g of Hexaammonium heptamolybdate tetrahydrate and 0.24 g of Bis[(+)-tartrato]diantimonate(III) dipotassium trihydrate, add water to dissolve, and make the volume to 300 mL. ← Sulfuric acid (2+1) 120 mL Mix ← Ammonium amidosulfate 5 g Mix and add water to make the volume to 500 mL</p> <p>Solution 2: L(+)-Ascorbic acid solution Weigh 7.2 g of L(+)-ascorbic acid and add water to make the volume to 100 mL</p> <p>Mix the solution 1 and solution 2 at the ratio of 5:1 to prepare the ammonium molybdate-ascorbic acid mixed solution (Prepare at the time of use)</p>

ANALYSIS RESULT



Calibration Curve of Phosphoric Acid



Absorption Spectrum of Phosphoric Acid

Addition Recovery Test for Drainage Water

Drainage water	Drainage water + 0.5 mg/L	Recovery rate
ND	0.48 ± 0.003	95.2 ± 0.7 %

ND: Not detected, n = 3

KEY WORDS

Environmental Analysis Related, Drainage Water, Environmental Chemistry, Environment, PO₄³⁻, Phosphate Ion, Absorption Spectrum, Calibration Curve, Coloring Reagent, Phosphoric Acid, UV, UH5300, U-5100, U-2900

Spectrophotometer
(UV)

Sheet No. UV120004-02