

## Investigating the limit of determination and uncertainty in ultra-thin Pd film measurements

### 1. Overview

This application brief reports on our investigation into quantitative limits and uncertainty when measuring ultra-thin Pd films 0.1  $\mu\text{m}$  or less in thickness using the SEA5000 manufactured by Hitachi High-Tech Science.

### 2. Quantitative Limits

#### (a) Definition

Limit of determination is defined by the following expression,

$$DL = 10\sigma_{BG}/m$$

where  $m$  is the slope of the calibration curve and  $\sigma_{BG}$  is the background standard deviation.

#### (b) Calibration Creation

Because the relationship of Pd film thickness and Pd X-ray intensity in an ultra-thin Pd film of 0.1  $\mu\text{m}$  or less can be approximated in a single equation, a linear calibration curve is created with a base and Pd 0.11  $\mu\text{m}$ . Figure 1 shows the calibration curve.




Figure 1 Calibration Curve

(c) Measurement conditions are shown in table 1.

**Table 1 Measurement Conditions**

|                        | Condition 1        | Condition 2        |
|------------------------|--------------------|--------------------|
| Collimator Size (mm)   | 0.1                | 1.8                |
| Voltage (kV)           | 50                 | 15                 |
| Current (uA)           | 1500               | 76                 |
| Analytical Line        | Pd-K $\alpha$ line | Pd-L $\alpha$ line |
| Sample Chamber         | Air                | Vacuum             |
| Measurement Time (sec) | 300                | 300                |

(d) Limit of determination

Limit of determination is shown in table 2. Here, background standard deviation is the actual measured value obtained from 10 measurements.

**Table 2 Limit of Determination**

|                     | Condition 1 | Condition 2 |
|---------------------|-------------|-------------|
| <i>m</i>            | 36.927      | 455.46      |
| $\sigma_{BG}$ (cps) | 0.054       | 0.254       |
| DL (um)             | 0.015       | 0.0056      |

### 3. Uncertainty

In this section we will consider uncertainty when measuring ultra-thin Pd film thickness by the fluorescent X-ray method.

(a) Method

Data for a Pd film thickness of 0.1 um or less can be output with a two-point calibration curve using the base and Pd 0.11 um. Pd film thickness  $t_0$  is given in equation (1) below.

$$(1) \quad t_0 = \frac{I_0 - I_1}{I_2 - I_1} \times t_1$$

$I_0$  = X-ray intensity of unknown sample  
 $I_1$  = Background intensity of base measurement  
 $I_2$  = X-ray intensity of Pd 0.11 um measurement  
 $t_1$  = Pd 0.11 um

(b) Measurement by fluorescent X-ray method

Table 3 shows the measured intensity of 10 measurements of 0.11 um and base, as well as a 0.05 um sample.

**Table 3 Measurement Conditions**

|                                      | Condition 1 |        |         | Condition 2 |          |          |
|--------------------------------------|-------------|--------|---------|-------------|----------|----------|
|                                      | 0.11μm      | 0μm    | unknown | 0.11μm      | 0μm      | unknown  |
| Average                              | 5.485       | 1.423  | 3.158   | 60.1465     | 10.0462  | 33.2315  |
| Standard Deviation                   | 0.180       | 0.054  | 0.135   | 0.598576    | 0.253706 | 0.241256 |
| Standard Deviation/ n <sup>1/2</sup> | 0.0570      | 0.0169 | 0.0427  | 0.189286    | 0.080229 | 0.076292 |
| Unknown thickness                    | 0.0470      |        |         | 0.0509      |          |          |

From here the standard uncertainty in each item I<sub>0</sub>, I<sub>1</sub>, I<sub>2</sub> is given. **The standard uncertainty of t<sub>1</sub>, from the guaranteed range of the standard sample being 10%, considering rectangular distribution, the value dividing 10% by 3<sup>1/2</sup> of the displayed value is the standard uncertainty.** Because all items in Equation (1) are not independent, the following formula must be applied in order to find the uncertainty of t<sub>0</sub> found in equation (1).

$$u(y(p, q, \dots)) = \left[ \left( \frac{\partial y}{\partial p} \right)^2 \times \{U(p)\}^2 + \left( \frac{\partial y}{\partial q} \right)^2 \times \{U(q)\}^2 + \dots \right]^{1/2}$$

**The contribution of each variable in this equation is shown by the value squared of partial differentiation by the variable to the square of each standard deviation.** Standard uncertainty can be calculated using a spreadsheet as in Tables 4 and 5 if all items are given. From the standard uncertainty obtained, if the expansion standard uncertainty is calculated as comprising of two coefficients, then, as shown below, there will be no significant difference in either party within the range of uncertainty.

**Condition 1:**                      **0.047μm ± 0.0060μm**  
**Condition 2:**                      **0.051μm ± 0.0059μm**

**Table 4 Uncertainty calculations under condition 1**

|                      |          | u(I0)       | u(I1)    | u(I2)    | u(t1)    |
|----------------------|----------|-------------|----------|----------|----------|
|                      |          | 0.042735374 | 0.016926 | 0.056966 | 0.006351 |
| I0                   | 3.1577   | 3.200435374 | 3.1577   | 3.1577   | 3.1577   |
| I1                   | 1.4233   | 1.4233      | 1.440226 | 1.4233   | 1.4233   |
| I2                   | 5.4853   | 5.4853      | 5.4853   | 5.542266 | 5.4853   |
| t1                   | 0.11     | 0.11        | 0.11     | 0.11     | 0.116351 |
| t0                   | 0.046968 | 0.048125281 | 0.046704 | 0.046318 | 0.04968  |
|                      |          | 0.001157285 | -0.00026 | -0.00065 | 0.002712 |
| Standard Uncertainty | 0.003031 | 1.33931E-06 | 6.96E-08 | 4.22E-07 | 7.35E-06 |

**Table 5 Uncertainty calculations under condition 2**

|                      |          | u(I0)       | u(I1)    | u(I2)    | u(t1)    |
|----------------------|----------|-------------|----------|----------|----------|
|                      |          | 0.076291954 | 0.080229 | 0.189286 | 0.006351 |
|                      |          |             |          |          |          |
| I0                   | 33.2315  | 33.30779195 | 33.23    | 33.23    | 33.23    |
| I1                   | 10.0462  | 10.0462     | 10.12643 | 10.0462  | 10.0462  |
| I2                   | 60.1465  | 60.1465     | 60.1465  | 60.33579 | 60.1465  |
| t1                   | 0.11     | 0.11        | 0.11     | 0.11     | 0.116351 |
|                      |          |             |          |          |          |
| t0                   | 0.050906 | 0.05107305  | 0.050807 | 0.050711 | 0.053841 |
|                      |          | 0.000167506 | -9.8E-05 | -0.00019 | 0.002936 |
| Standard Uncertainty | 0.002948 | 2.80584E-08 | 9.62E-09 | 3.8E-08  | 8.62E-06 |