

Chromaster

-Automation of Pretreatment by Column-Switching-

The analysis of biological specimens and samples containing large quantities of impurities requires pretreatment, such as removal of the impurities and the condensation of the sample. Pretreatment, exerting a significant impact on the results of the analysis and being a critical element, to a large extent depends on the experience and skill levels of the individual performing the analysis. Therefore, reliable pretreatment operations are of paramount importance.


Among the processes that are widely employed in pretreatment are solid-phase extraction and column-switching. The Chromaster system permits the automation of pretreatment through the use of column-switching combined with the use of optional valves.

The following is an example of analysis and evaluation of drugs found in blood serum.

[Analysis of phenytoin in blood by column-switching]

Sample: Standard phenytoin reagent added to human serum
(the sample was centrifuged, and the supernatant was used as an injection specimen)

Purpose: To study the automation of pretreatment for phenytoin analysis



Given that phenytoin is toxic from 10 to 20 µg/mL (effective concentration in blood), and 20 µg/mL or higher in in-blood concentration, it has been verified that the method employed in this procedure can accurately determine the amount of phenytoin in blood.

[Configuration of the system employed]


Chromaster 5110 pump x 2
Chromaster 5210 AutoSampler
Chromaster 5310 Column Oven
(with a 6-way 2-position valve)
Chromaster 5420 UV-VIS Detector

LC Conditions: Analysis column

Mobile phase
(A) 50 mM KH₂PO₄-K₂HPO₄(pH6.9)/CH₃CN=95/5
(B) CH₃CN
*(A)/(B)= 65/35
Flow rate 1.0 mL/min
Column LaChrom C18 (5 µm) 4.6 mm I.D. × 150 mm
Temperature 40°C
Detection 210 nm

LC Conditions: Pretreatment column

Mobile phase
(A) 50 mM KH₂PO₄-K₂HPO₄(pH6.9)/CH₃CN=95/5
(B) CH₃CN
*(A)/(B)= 95/5
Flow rate 1.0 mL/min
Injection vol. 20 µL
Column Mspak PK-4A (4.0 mm I.D. × 10 mm Shodex¹)
Temperature 40°C



¹This is a pretreatment column for column-switching analysis; it is designed to retain low-molecular weight components only, without retaining protein.


Chromaster

-Automation of Pretreatment by Column-Switching-

[Evaluation of a column-switching pretreatment method]

The analysis of phenytoin (anti-epilepsy agent) in blood requires the removal of protein and other components. In view of this fact, we studied a method of exclusively introducing the target component, phenytoin, into the analysis column by removing protein through the use of a pretreatment column.

Steps leading from the automation of pretreatment to the analysis




Elutes the protein in the serum and at the same time exclusively traps the drug.

Only the drug trapped in the pretreatment column is eluted.

During the analysis, the flow path is switched and simultaneously the pretreatment column is cleaned and stabilized in preparation for another analysis run.

[Evaluation of the conditions for the pretreatment column (1)]

Purpose: To verify the conditions under which protein is eluted



LC Conditions: Pretreatment column

Mobile phase
 (A) 50 mM KH_2PO_4 - K_2HPO_4 (pH6.9)/ CH_3CN =95/5
 (B) CH_3CN
 *(A)/(B)= 95/5
 Flow rate 1.0 mL/min
 Injection vol. 20 μL
 Column Mspak PK-4A (4.0 mm I.D. \times 10 mm Shodex
 Temperature 40°C

Results: It was verified that proteins are not retained in the column.


Chromaster

-Automation of Pretreatment by Column-Switching-

[Evaluation of a column-switching pretreatment method (2)]

[To verify the conditions under which the target component (phenytoin) is introduced into an analysis column.]

The evaluation was conducted under the following conditions:
phenytoin standard, 10 mg/L (injection volume: 20 μ L); detection: 254 nm.



LC Conditions: Pretreatment column

Mobile phase

(A) 50 mM KH_2PO_4 - K_2HPO_4 (pH6.9)(B) CH_3CN

Flow rate 1.0 mL/min

Column Mspak PK-4A (4.0 mm I.D. \times 10 mm)

Temperature 40°C



In this manner, it can be verified that the phenytoin can be eluted by increasing the concentration of organic solvent.




It can be verified that the phenytoin does not elute under these conditions.

Evaluation of analysis column conditions

To verify the conditions under which the target component (phenytoin) is retained in the analysis column and separated.

The evaluation was conducted under the following conditions: phenytoin standard, 10 mg/L (injection volume: 20 μ L); detection: 254 nm.



LC Conditions: Analysis column

Mobile phase

(A) 50 mM KH_2PO_4 - K_2HPO_4 (pH6.9)(B) CH_3CN

Flow rate 1.0 mL/min

Column LaChrom C18 (5 μ m) 4.6 mm I.D. \times 150 mm

Temperature 40°C

Results:

Phenytoin was introduced from the pretreatment column to the analysis column.

The analysis condition was designated as A/B=60/40.

Chromaster



-Automation of Pretreatment by Column-Switching-

[Evaluation of valve-switching time]

(1) Introduce serum into a pretreatment column; verify the amount of time required for protein removal (protein detected at 280 nm).

- For these tests, the valve-switching time was set to 5 minutes.



(2) The concentration of the mobile phase was determined based upon the results of the previous pretreatment column condition evaluation.



[Results of analysis by column-switching]

Sample: A standard phenytoin sample (10 µg/mL) added to human serum

* Supernatant injected after centrifuge separation



Calibration line for phenytoin added to serum (1 - 50 µg/mL)

The rate of recovery of phenytoin by the automation of pretreatment was 85.3%. At 10 µg/mL, the reproducibility of elution time was CV: 0.04%, and the reproducibility of the area under the curve was CV: 1.91%. Given that the toxic region of phenytoin is 10 to 20 µg/mL in effective concentration in blood, and 20 µg/mL or higher in in-blood concentration, it has been verified that the method employed in this procedure can accurately determine the amount of phenytoin in blood.

NOTE: These data are an example of measurement; the individual values cannot be guaranteed.

The system is for research use only, and is not intended for any animal or human therapeutic or diagnostic use.