

Chromaster

-Simultaneous analysis of four agricultural chemical components-

Iprodione, asulam, thiophanate-methyl, and siduron are four agricultural chemicals that are designated as water quality management target items. Simultaneous analysis of these chemicals typically requires the use of LC/MS. The accuracy of this analysis must be within 20% of the coefficient of variation, and is determined by analyzing samples that are 1/100 of the concentration of target values (0.003 mg/L or 0.002 mg/L for asulam). If this sample is concentrated by a factor of 500, the resulting concentration (1.5 mg/L or 1.0 mg/L for asulam) is within the measurement capability of a high-sensitivity DAD system (See Technical Data No. 187).

The text below describes examples of analyses performed.

Analysis of four agricultural chemical components

[Component and Target value]

Name of chemical	Component	Molecular formula	Application	Target value [mg/L]
Iprodione	3-(3, 5-dichlorophenyl)-N-isopropyl-2, 4-dioxoimidazolidine-1-carboxamide	C ₁₃ H ₁₃ Cl ₂ N ₃ O ₃	Germicide	0.3
Asulam	methyl sulfanilylcarbamate	C ₈ H ₁₀ N ₂ O ₄ S	Herbicide	0.2
Thiophanate-methyl	dimethyl 4, 4' -(o-phenylene) bis (3-thioallophanate)	C ₁₂ H ₁₄ N ₄ O ₄ S ₂	Germicide	0.3
Siduron	1-(2-methylcyclohexyl)-3-phenylurea	C ₁₄ H ₂₀ N ₂ O	Herbicide	0.3

Sample preparation: Standard solution: Prepared with acetonitrile so that the concentration will be 0.2 to 20 mg/L.

[Solid phase extraction procedure]

SPE Cartridge NOBIAS RP-SG1*

Conditioning

- ← Acetonitrile 20 mL
- ← Diluted nitric acid (pH3.5) 20 mL

Sample loading

- ← 500 mL
 - Mixture of pesticides
 - Sodium ascorbate 20 mg/L
 - 10% EDTA 2 mL
 - pH3.5 (Nitric acid)

Washing

- ← Pure water 10 mL

Drying

- ← Blow the N₂ gas into the cartridge

Elution

- ← Acetonitrile 3 mL

Concentration

- ← Blow the N₂ gas into the eluate and fill up to 1 mL by eluent

Measurement (HPLC)

* NOBIAS RP-SG1


A reverse-phase polymer-based solid-phase packed column.

Chromaster

-Simultaneous analysis of four agricultural chemical components-

[Example of analysis of a standard sample]

Sample concentration: 1.0 mg/L per component



[System configuration]

Chromaster 5110 Pump
 Chromaster 5210 Autosampler
 Chromaster 5310 Column Oven
 Chromaster 5430 DAD
 Empower 2 Data Processing System

[LC conditions]

Column	: HITACHI LaChrom C18 (3 μ m) 4.6 mm I.D. \times 150 mm
Eluent	: 50 mM potassium dihydrogen phosphate (pH 3.0)/ acetonitrile = 45 / 55
Flow rate	: 1.0 mL/min
Temperature	: 40°C
Detection	: 270 nm Asulam 230 nm thiophanate-methyl, siduron, iprodione
Inj.vol	: 10 μ L

[Linearity]

At a 0.2 to 20 mg/L concentration, the sample produces a favorable straight line.

[Asulam] R^2 : 0.9998, [Thiophanate methyl] R^2 : 0.9999

[Siduron(1)] R^2 : 0.9999, [Siduron(2)], R^2 : 0.9999, [Iprodione] R^2 : 0.9999