

Yoshikazu Nishiyama¹, Masayuki Iwasa¹, Brian Goolsby², Kenichi Shibata³
¹Hitachi High-Tech Science Corp., ²Hitachi High Technologies America, Inc.,
³Hitachi High-Technologies Science America, Inc.

New thermal analyzers were developed that can observe samples during measurements. These systems can obtain optical images and collect thermal data including DSC, TG, DTA, DMA. Furthermore, the optical images are linked in the time-scale with thermal data, so we can properly correlate the phenomena occurring during thermal analysis.

We measured polyethylene terephthalate (PET) undergoing tensile stretching with this system. During glass transition, the observed physical transformation of PET was well correlated with changes in the DSC and DMA signal. Furthermore, we found that PET expanded unexpectedly during decomposition in air. When we used aluminum pans for measuring this decomposition, the aluminum pan was deformed by PET expansion, causing a dummy peak on the DTA signal. We report the results of various polymer samples analyzed by the new system.

Experimental

Material	PET			PE	Rubber (NBR)
System	TG/DTA	DSC	DMA (1 Hz oscillation)	TG/DTA	TG/DTA
Sample Detail	5 mg	10 mg	20 mm L x 10 mm W x 0.1 mm T	5 mg	5 mg
Heating Program	10°C/min	10°C/min	2°C/min	20°C/min(N ₂) to 230°C and hold at 230°C	20°C/min(N ₂) , 10°C/min(Air)
Purge Gas	Air or N ₂	N ₂	N ₂	N ₂ →O ₂ (change after temp. stability)	N ₂ →Air (change after drop to 300°C)

Notes: In DSC and TG/DTA, PET sample is from a drink bottle. For DMA, PET film is a consumer product. PE sample is cut from a water pipe. NBR is a consumer product.

The USB-linked camera can be configured to take pictures based on various parameters (For example: every 1°C temperature change, every second, over a range of times, etc.). Furthermore, image data are linked to thermal analysis data, so data and images can be overlaid using the analysis software.






Fig.1 Thermal analyzers with optical observation units (red arrows) by Hitachi High-Tech Science Corp.

Results


TG/DTA of PET




DSC of PET




DMA of PET




OIT measurement of PE



While sample A oxidation starts before sample B, the TGA curve indicates that sample A decreases less than sample B. This result means the decomposition reaction of sample B is more drastic. In pictures #4 to #7, we can see sample B collapses and flows widely, while sample A is stable. So, from the TG/DTA curve and pictures, we can determine that sample B had a more drastic decomposition reaction than sample A.



TG/DTA of NBR



Conclusion

Using thermal instruments with optical observation, we can collect much more information about sample behavior during analysis.

This can help explain both the expected phenomena polymers undergo as well as anomalous or unexpected behaviors. This technology could be an important tool in education, failure analysis, and experiment design, to name a few.