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IPC TM-650 Test Methods Manual 2.4.25D

Glass Transition Temperature and Cure Factor by DSC 11/17

Determination of T_g The midpoint temperature T_m (°C) as described in Figure 1 reported as the T_g . T_g is the point on the thermal curve corresponding to 1/2 the heat flow difference between the extrapolated onset and extrapolated end. If suitable computer software is available, the automatic calculation of the glass transition temperature is allowable provided the value calculation is either the midpoint or the steepest deflection and not the onset temperature. See Figure 1.

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2.4.24.1 Time to Delamination (TMA Method)

2.4.24.3 Glass Transition Temperature of Organic Films

2.4.24.5 Glass Transition Temperature and Thermal Expansion of Materials Used In High Density Interconnection (HDI) and Microvias

2.4.24C Glass Transition Temperature and Z-Axis Thermal Expansion

Glass Transition Temperature – T_g Construct a tangent line to the curve above and below the transition in the curve. The temperature where these tangents intersect is the TMA determined T_g for the material. (中略) For consistency, it is recommended that the TMA computer analysis software be used for this calculation (see Figure 2).

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Depending on sample preparation, two heating cycles may be required to obtain accurate T_g and CTE above T_g . If the sample shows unexpected shrinkage above T_g (see Figure 1), the two heat test method is required. If the sample does not show anomalous behavior, only one heat cycle (the second heat cycle at 5°C/min) is required.

- First heat: The first heat cycle of the specimen shall be run at 10°C/min.
- Second heat (reportable data heat cycle): The second heat cycle of the specimen shall be run at 5°C/min.

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
Figure 1 TMA Expansion Curves: First Heat Cycle and Second Heat Cycle

The time to delamination is determined as the **time from the onset of the isotherm to failure**. Failure is any event or deviation of the data plot **where the thickness is shown to have irreversibly changed**. The scan in Figure 1 is typical of an epoxy material at 260°C [500°F] isothermal temperature. On occasion, some materials will delaminate before the isotherm is reached. In this case, the temperature at the time of failure shall be recorded.

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2.4.24.2 Glass Transition Temperature of Organic Films

2.4.24.4 Glass Transition Temperature and Thermal Expansion of Materials Used In High Density Interconnection (HDI) and Microvias

Glass Transition Temperature (T_g) Construct a tangent line to the curve below the transition temperature in the modulus curve. Construct a tangent to the storage modulus curve at or near the inflection point approximately midway through the step change in the transition. The temperature where these tangents intersect is the reported T_g for the material. For consistency it is recommended that the DMA computer analysis software be used for this calculation. See Figure 1 for an example of this tangent intersection method.

Storage Modulus (E') The sample storage modulus (E') shall be calculated at room temperature (22°C) and reported in units of Pa (N/m²). For consistency it is recommended that the DMA computer analysis software be used for this geometry specific calculation.

Alternative thermal transitions may be reported as the transition peak temperature in the sample loss modulus (T_i) or tan δ plots (T_i) (see Figure 2 and Figure 3).


Figure 1 DMA Modulus Plot


Figure 2 DMA Tan Delta Plot

Note: T_i is the transition peak temperature.




Figure 3 DMA Loss Modulus Plot

Note: T_i is the transition peak temperature.

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
Reference


1. IPC-TM-650 TEST METHODS MANUAL 2.4.24.1,

2. The Institute for Interconnecting and Packaging Electronic Circuits, (1994)

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The 8th International and the 10th Japan-China Joint Symposium on Calorimetry and Thermal Analysis (CATS-2017), November 24, 2017, Fukuoka University

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
Method A Flexural bending geometry – thick specimens (>0.5 mm): Specimens shall be approximately 8 mm to 12 mm wide, 20 mm to 40 mm long, and 1 mm to 2 mm thick. The thickness shall be a minimum of 0.5 mm; for thicknesses <0.50 mm, use Method B. An aspect ratio of length/thickness = 10/1 or greater should be maintained. Exact specimen dimensions should be determined by the apparatus used.

Method B Thin film tension geometry – thin specimens (<0.50 mm): Specimens shall be approximately 15 mm to 20 mm long and 2 mm wide. The minimum thickness is determined by the strength of the material; it should not break during testing. Exact specimen dimensions may be determined by the apparatus used.

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