

## Resolving Individual Magnetic Domains in Single-Crystal $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> via High-Sensitivity *in vacuo* Magnetic Force Microscopy

SHEET No. 011

Instrument: Hitachi AFM5300E with High-Vacuum and Environmental Control

### Introduction

Iron(III) oxide (Fe<sub>2</sub>O<sub>3</sub>) has four known phases:  $\alpha$ -,  $\beta$ -,  $\gamma$ -, and  $\epsilon$ -. While  $\alpha$ - and  $\gamma$ -Fe<sub>2</sub>O<sub>3</sub> exist naturally in minerals, the other two phases ( $\beta$ - and  $\epsilon$ -) are typically generated through a synthetic route. In 2004, Prof. Shinichi Ohkoshi, *et al.* from the University of Tokyo reported the novel observation of an extremely large coercive field (20 kOe) comprised of nanometer-sized  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> that was generated via nanoparticle synthesis. This discovery led to extensive investigations on the magnetic properties and domain structures of  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub>. This application data sheet demonstrates Hitachi's high-sensitivity magnetic force microscopy (MFM) studies of  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> with single-crystal rod morphology.<sup>1)</sup>




Figure 1: Crystal structure and unique magnetic characteristics of  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub>. (a) The  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> unit cell viewed along the b-axis (left) and a-axis (right). (b) Schematic image of a bar magnet based on a mesoscopic  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> rod. (c) Magnetization versus external field plot of the  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> rod measured at room temperature along a single in-plane direction.

## Results

Figure 2 shows high-sensitivity MFM imaging of a  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> single-crystal ferrite bar magnet in vacuum. Several of these bar magnets are fixed on a common substrate and each bar magnet morphology is that of a single-crystal  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> with its own spontaneous magnetization. The N- and S-pole regions of individual rods as well as overlapping rod complexes were clearly observed with MFM.




Figure 2: (a) Hitachi high vacuum AFM5300E. (b) Topography and (c) MFM images of an  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> single-crystal ferrite bar magnet. (d) Topography-MFM overlay image. (e) Magnetization arrangement model.

## More Applications of $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub>

Single-crystal  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub> exhibits a resonance frequency in response to the magneto-optical effect at the terahertz level (10<sup>12</sup> Hz). Therefore, one example of utilizing this property could be to serve as a high-frequency microwave-absorbing material in advanced highway safety systems.

Additionally, single-crystal  $\epsilon$ -Fe<sub>2</sub>O<sub>3</sub>, when attached to an AFM cantilever, can be utilized as a high-performance MFM probe due to its high coercivity (H<sub>c</sub>).

Author: T. Yamaoka, Analytical Application Engineering Section Tokyo 2, Hitachi High-Tech Science Corporation