

Quantified Elastic Modulus Mapping of Rubber Blend Under Temperature Control Using Hitachi Proprietary SIS-QuantiMech




Fig.1 Temperature-dependent elastic modulus mapping of rubber blend revealed by SIS-QuantiMech (Model: AFM5300E)

Hitachi SIS-QuantiMech offers quantified mechanical-property mapping based on the collection of force curve via Hitachi propriety Sampling Intelligent Scan (SIS) at single pixel level and subsequent derivation of the modulus of elasticity that is governed by AFM tip sharpness, spring constant, and the deflection sensitivity of the AFM cantilever as well as selection of the most appropriate model to fit the data.

Fig. 1 shows the elastic modulus mapping of a rubber blend composed of both natural rubber (NR) and butadiene rubber (BR) at three different temperatures. Fig. 2 plots the peak value of the elastic modulus distribution obtained in Fig. 1 and measurements of the storage modulus by a dynamic viscoelasticity device (DMA).




Fig.2 Result of elastic modulus of rubber blend by DMA and AFM (SIS-QuantiMech)
Models: AFM5300E and DMA7100

The derived modulus of elasticity measured by the AFM at the low temperature of -120 °C was about 1.9 GPa, which is in agreement with the storage modulus measured by the DMA. At the room temperature (25 °C), the modulus of elasticity is significantly lower at round 3 MPa. Both elastic modulus mappings are with very little contrast variation due to the fact that the two rubber components are either all in a glassy state at -120 °C or all in a rubbery state at 25 °C. On the other hand, binary contrasts were exhibited in the elastic modulus mapping at -90 °C. As the T_g of BR is lower than -90 °C and T_g of NR is higher than -90 °C, regions of the dark contrast with a calculated modulus of elasticity of 0.24 GPa and those areas with a brighter contrast with the calculated modulus of elasticity of 1.9 GPa can be assigned as the BR and NR domains, respectively.

In this case study, it is demonstrated that SIS-QuantiMech can render high-resolution and localized measurements of elastic modulus that are not achievable by macro analysis techniques such as DMA. It also proves that SIS-QuantiMech in conjunction with temperature control is an effective approach for compositional mapping of polymer complex, especially for those containing components with a very low T_g. In this regard, Hitachi's high-vacuum AFM5300E provides superior and wide-ranging temperature control from -120 °C to 800 °C and thus allows for the component differentiation of NR/BR mixtures by sample cooling at -90 °C.

Authors: M. Iwasa, Hitachi High-Tech Science Corporation
M. Iyoki, Hitachi High-Technologies Corporation

Measurement of Elastic Modulus of Rubber Blend under Cooling Using SIS-QuantiMech



Recommended configuration	Recommended configuration
Environment Control Unit AFM5300E	
▪ Evacuation System 1 for (TMP+RP)	
▪ Heat/Cool Sample Stage	
▪ Temperature Controller	
▪ Cantilever: SI-DF3P2	
Probe Station AFM5000 II	
▪ SIS-ACCESS/SIS-QuantiMech	
Dynamic Mechanical Analyzer DMA7100	



AFM5300E



DMA7100