

On-Site Non-Destructive Analysis of Colorants in Cultural Properties



Susumu Shimoyama

PhD*1

Color Materials Science Laboratory, Den Material Company Ltd.

Professor Emeritus, Kibi International University

Table of Contents

1. Infrared Photographs
 2. Microscopy Images
 3. Visible/Near-infrared Reflection Spectra
 4. X-ray Fluorescence Spectra
 5. Three-dimensional Fluorescence (3DF) Spectra
 6. Conclusions
-

Cultural properties are precious assets that furnish living, breathing witness to history and heritage. For this reason, scientific studies of cultural properties must rely exclusively on non-destructive analytical techniques, avoiding any chemical or physical processing steps that may risk damaging precious artifacts. In this article we focus in particular on the *colorants* used in cultural properties such as Ryukyuan *bingata* garments, lacquerware, ceramic pottery, *Nihonga* paintings, oil paintings, and *ukiyo-e* woodblock prints. We describe an experimental technique—including both analytical instruments for on-site measurements and methods for processing the information they provide—for analyzing such colorants non-destructively. This technique, which has been used extensively for practical field work, involves a multi-stage procedure for non-destructive analysis of colorants in cultural properties. We begin by taking infrared photographs to determine the infrared absorption properties of the colorants. Then we select points at which to make measurements and capture microscopy images at these points to yield enlarged views of the coloration states. Next, we obtain visible/near-infrared reflection spectra at different measurement points and use the shape of these spectra to analyze the colorants. Finally, we measure X-ray fluorescence (XRF) spectra and attempt to use the resulting elemental analysis data to identify pigments; when pigments cannot be identified in this way, we measure three-dimensional fluorescence (3DF) spectra and use the results to identify dyes.

*1 Professor Emeritus, Kibi International University and Director, Den Material Company Ltd. Color Materials Science Laboratory, (3-4-19 Kadotabunkamachi, Naka-Ku, Okayama City, Okayama 703-8273 Japan). PhD in Science and Engineering, Iwaki Meisei University Graduate School of Science and Engineering, Department of Science and Engineering. Author of "Research on Non-destructive Methods for Identifying Dyes Used in Ancient Dyed Fabrics and Ukiyo-e Prints from 3D Fluorescence Spectra" (1998).

1. Infrared Photographs

1-1. Instruments

Our infrared photography setup is shown in Figure 1. We use a PENTAX 645D IR camera (effective pixel count: approximately 40 million pixels) equipped with an smcPENTAX-FA645 75 mm F2.8 lens. As light sources we use two Canon SPEEDLITE 420EX flash units with stroboscopic capabilities. On both the front surface of the camera lens and the light-emitting surfaces of the flash units we install IR-86 filters (a sharp-cutoff filter from Fuji Film that transmits infrared wavelengths of 860 nm or longer). We use a wireless remote-flash controller (NEEWER), triggered by the camera with receivers on both flash units, to synchronize the flash units with the shutter. This setup allows us to capture clear images with no need to operate in darkened surroundings: The flash units irradiate the sample surface with infrared light at wavelengths longer than 860 nm, and light reflected from the sample is received by the camera and recorded to yield digital images. Note that we removed the band-pass filter that is usually installed inside the camera to eliminate ultraviolet and infrared wavelengths and transmit only visible light.




Fig. 1 Setup for infrared photography.

1-2. Information derived from infrared photography¹⁾

The longer wavelengths of infrared light compared to visible light ensure that infrared light is not reflected from the surfaces of material bodies—as is visible light—but instead penetrates into the body interior. For example, if the colorants in a *bingata* garment transmit infrared light without absorption, the light penetrates to deeper layers of the garment before eventually being reflected from the surface of fibers such as silk, cotton, or hemp; note that the reflected light passes through the colorant a second time as it exits the garment. Thus, when a colorant with high infrared transmittance is irradiated with infrared light for a photograph, the light is reflected at fiber surfaces and passes through the colorant a second time as it exits the sample to be captured and recorded by the camera, yielding white regions in captured images. In contrast, if the sample contains colorants that strongly *absorb* infrared radiation, incident infrared light is absorbed by the colorants before it can penetrate deeper into the material; in this case the light does not reach the fiber surfaces and is thus not reflected back toward the camera, yielding dark regions in captured images. Intermediate between these extremes are colorants with moderate infrared transmittance; in this case, only the transmitted portion

Table 1 Behavior of colorants (pigments, metals, and dyes) in infrared photographs.

Pigment (alternative name) or metal		Appearance in infrared photographs
White	Chalk	☐ White
	White lead (silver white)	☐ White
	Chinese white (zinc white)	☐ White
Red	Vermilion	☐ White
	Red lead (minium)	☐ White
	Bengala (red ochre)	■ Grey
Yellow	Orpiment	☐ White
	Chrome yellow	☐ White
	Yellow ochre	■ Grey
Green	Malachite	■ Black
	Viridian	☐ White
	Emerald green	■ Black
Blue	Prussian blue	■ Black
	Azurite	■ Black
	Smalt	■ Grey
Black	Sumi (carbon black)	■ Black
Metals	Gold	☐ White (surface reflections)
	Brass (copper/zinc alloy)	
	Silver	
	Tin	
Dye (coloring matter)		Appearance in infrared photographs
Red	Safflower (carthamin)	☐ White
	Lac (laccic acid)	
	Japan madder (pseudo-purpurin)	
	Sappanwood (braziliin)	
Yellow	Gardenia (crocin)	☐ White
	Turmeric (curcumin)	
	Amur cork (berberine)	
	Gamboge yellow (gambogic acid)	
Blue	Indigo	☐ White
Purple	Purple gromwell (shikonin)	☐ White
Black	Tannin black	■ Grey - ■ black

of the incident radiation proceeds deeper into the sample to be reflected at fiber surfaces, yielding a lower-intensity reflected signal that produces grey regions in infrared photographs. Also, beneath high-transmittance colorants there may exist written characters or line drawings marked in *sumi ink* (or charcoal); although these are not visible by eye, they appear as dark regions in infrared photographs. In *bingata* garments and similar samples, ultra-fine ink lines drawn by hand on infrared-transmitting colorants appear as dark lines in infrared photographs. This is due to the fact that the ink used, which consists primarily of carbon, exhibits high infrared absorbance. Colorants (pigments, metals, and dyes) appearing as white, black, or grey regions of infrared photographs are listed in Table 1.

Analyzing infrared photographs of model oil painting

When sketches beneath oil paintings are drawn in charcoal (carbon), it is well known that these sketches may be seen in infrared photographs. In fact, this phenomenon is not limited to charcoal, but also applies to infrared-absorbing Prussian blue, as we confirmed by using a HORIBA XGT (X-ray Guide Tube)-5000 X-ray analytical microscope to carry out an elemental-mapping analysis of "Daubigny's Garden"—one of the best-known oil paintings of van Gogh's final years, held by the Hiroshima Museum of Art (dimensions: 530 mm tall × 1030 mm wide)—to determine whether or not a black cat, said to be drawn by van Gogh himself, was present behind the grass (that is, underneath the paint used to depict the grass) in the lower-left portion of the painting (Figure 2)^{2,3}. To this day, the black cat remains hidden by brush strokes added by an unknown third party.

In this case, we cannot conduct experiments on the actual cultural property itself. Thus, when a preliminary investigation is required we perform experiments on models that we fabricate ourselves. For the analysis of Figure 2, we used a model of van Gogh's painting to test whether or not our elemental-mapping analysis could successfully identify a black cat hidden beneath the paint layers. This model was constructed by painter Hiroshi Yoshida, who prepared a facsimile—using the same paints van Gogh is thought to have used—of the region of the painting believed to be hiding the black-cat drawing. Yoshida began by blending Prussian blue with chrome yellow to yield "chrome green," then added more Prussian blue to yield a dark blue—nearly black—hue with which he painted a black cat (Figure 3A). Next, using colors thought to be used by van Gogh—primarily zinc white with some silver white, and with chrome yellow as a yellow hue and viridian and emerald green as green hues—Yoshida painted over the black cat to replicate the appearance of the grass in the van Gogh original, superposing multiple grains of grassy textures on top of each other until the black cat was hidden (Figure 3B). Finally, on top of this Yoshida added additional strokes of red ochre—as the unknown third party is believed to have used—to complete the preparation of our model painting⁴.

Figure 3C shows an infrared photograph of the model painting. Two of the paints used to prepare the facsimile—Prussian blue and emerald green—are infrared-absorbing and thus appear as black regions in the photograph. Of the other paints, zinc white, silver white, chrome yellow, and viridian appear as white regions, while red ochre appears as grey. Although infrared-absorbing emerald green is present at points scattered throughout the grass in the painting, the black cat appears to be drawn primarily in Prussian blue. This is why the cat is visible in the infrared photograph despite being hidden beneath paint layers depicting the grass. This demonstrates the ability of infrared photography to reveal features drawn in charcoal or other infrared-absorbing substances hidden beneath paint layers.

Unfortunately, our infrared photographs of the original "Daubigny's Garden" did *not* reveal a black cat; we attribute this to the thickness of the paint layers being greater than we had anticipated. However, an elemental-mapping analysis based on scanning X-ray fluorescence (XRF) measurements—performed on-site (in a temporary laboratory set up at the museum) using an X-ray analytical microscope (Figure 4)—revealed, in the mapping image for elemental chrome, drawings outlining the head, neck, front paws, torso, and tail of a cat, thus confirming that van Gogh did in fact draw a black cat beneath his masterpiece.





Fig. 2 Location of black cat hidden beneath van Gogh's "Daubigny's Garden".




A. The black cat and the surrounding grass are painted first.

B. Additional brushstrokes are then added until the black cat is hidden, completing the model painting.

C: An infrared photograph of the model painting (B) clearly reveals the black cat.

Fig. 3 Model painting prepared to reproduce portion of "Daubigny's Garden" (A, B) and infrared photograph of this model painting (C).



※The sample stage is translated horizontally while measuring X-ray fluorescence emissions induced by X-ray beam 100 μm irradiation.

Fig. 4 Original "Daubigny's Garden" mounted in X-ray analytical microscope for scanning X-ray fluorescence measurements.

Analyzing infrared photographs of *bingata* garments

Figure 5 shows standard and infrared photographs of a *bingata garment**². With the exception of the black tannin-based mordant used as a black dye, most of the dyes, metals, and pigments used in this garment—including the *chalk* (Oyster Shell White) and *white lead* used as white pigments, *vermilion* and *red lead* used as red pigments, and *orpiment* used as a yellow pigment—are infrared-transmitting and appear as white regions in the infrared photograph. The garment also contains some colorants that are infrared-absorbing and appear as black regions in the infrared photograph—including *malachite* and *emerald green* for green coloring, *Prussian blue* and *azurite* for blue coloring, and *sumi* ink for black coloring—as well as some colorants that are partially infrared-absorbing and appear as grey regions in the infrared photograph, including *Bengala*, which contains iron oxide, for red coloring, *yellow ochre* for yellow coloring, and *smalt* (smalt mineral glass) for blue coloring. This demonstrates how white, black, and grey features in infrared photographs can be used to narrow down the list of candidates for the colorants that may be present in a specimen. This information, in combination with information obtained by other non-destructive analytical methods, can often be used to identify these colorants precisely.

For example, comparing the color photograph on the left with the infrared photograph on the right in Figure 5, we see that the base yellow color of the garment corresponds to white regions in the infrared image, while the blue and green tail feathers of the Chinese phoenix correspond to black regions. From this we infer, first, that the yellow base color of the garment is due either to *orpiment* (As_2S_3) or to *yellow dye*, as these two colorants appear white in infrared photographs. If XRF spectral measurements, discussed below, detect the presence of arsenic (As), we conclude the garment contains orpiment. On the other hand, if the XRF spectra do not suffice to settle the question, we measure 3DF spectra and identify the dyes present from the unique fluorescence fingerprints of their molecular structures.

Next, because the blue portions of tail feathers in the garment appear as black regions in the infrared photograph, we infer that the blue colorant is either *Prussian blue* ($KFe[Fe(CN)_6]_3 \cdot nH_2O$) or *azurite* ($Cu_3(CO_3)_2(OH)_2$). If XRF spectral measurements detect iron (Fe), we conclude that it is the former; if XRF measurements detect copper (Cu), we conclude that it is the latter.

Finally, because the green portions of tail feathers in the garment appear as black regions in the infrared photograph, we hypothesize that the green colorant is *malachite* (basic copper carbonate $CuCO_3 \cdot Cu(OH)_2$). This will be confirmed if a subsequent XRF detects the presence of Cu.

Actually, at this point we pause to mention one insight gleaned from the many *bingata* specimens we have non-destructively analyzed thus far: in general, *azurite* and *malachite* are never used to yield color features in *bingata* garments. The reason for this is that, when a brush is used to apply a pigment to fabric (fibers), large pigment particles fail to become thoroughly embedded among the cloth fibers, and the fabric does not retain the intended color. Thus, the colorants used to brush patterns onto *bingata* garments must be pigments containing *fine* particles and yielding vivid colors. However, although both *azurite* and *malachite* produce deep and vivid colors when used as *large* particles, when used as *small* particles they yield only mild "whitish" colors typically known as "whitish-blue" and "whitish-green." Consequently, green features in *bingata* garments are typically produced by combining *Prussian blue*—which yields a vivid blue color even as small particles—with *orpiment*, which yields a vivid yellow color. In the case of Figure 5, the presence of Prussian blue in green regions causes these regions to appear black in infrared photographs. This is confirmed by the XRF spectral measurements discussed below, in which green regions of garments are found to contain both Fe (in Prussian blue) and As (contained in orpiment), but *not* Cu (in azurite and malachite).

*² Decorative patterns are added to *Bingata* garments by cutting out a piece of paper in the shape of, e.g., a peony, an iris, or a bird, placing this paper atop the garment, using an applicator tool known as a hera to apply a paste, and allowing the paste to dry. Then, a hake (or pressing brush) is used to press colorant into remaining regions of the pattern in which paste was not applied, after which colorant is painted over the region or a color shading is applied. This is known as hueing. After this step, the paste is rinsed away with water. Paste is then applied to the pattern subjected to hueing to hide it, and it is dried. Finally, the hake is used to apply the base color to fabric regions to which no paste was applied.



Fig. 5 Standard color photograph (left) and infrared photograph (right) of *bingata* garment.

2. Microscopy Images

2-1. Instruments

To capture microscopy images, we use a PENTAX WG-3 camera (lens focal length: 4.5-18.0 mm, aperture range F2.0(W) - F4.9(T), effective pixel count: approximately 16 million pixels) equipped with a micro-stand ring. As shown in Figure 6, for microscopy images of *bingata* garments we place a protective sheet—from which is cut out a circular hole the size of the camera lens—on the garment specimen, then place the camera on this sheet with its lens aligned in the hole. (The weight of the camera, including batteries and SD card, is 247 grams). Then, referring to the camera monitor, we move the camera until the image is centered on one of the predetermined measurement points, set the optical magnification to digital microscopy mode and capture images via remote control. The available optical magnifications are 1.2 \times , 2.5 \times , and a maximum of 4.0 \times , and the image data captured for each magnification can be viewed on a PC.



Fig. 6 Setup for capturing microscopy images.

2-2. Information derived from microscopy images

Microscopy images of *bingata* garments present enlarged views of the colorants used for the base color and for accent colors in the fiber structures selected as measurement points. If these colorants are *dyes*, then the images allow us to observe transparent dyes for chromatic colors penetrating into the interior of fibers to achieve coloration, and also to observe the structural composition of the fabric. On the other hand, if the colorants are *pigments*, then we observe non-transparent colorants affixed to the surfaces of fibers; in this case, there will be defect regions, through which the base fibers will appear as white regions in microscopy images. For example, Figure 7 shows a microscopy image (optical magnification 4.0 \times) of red flower petals within a peony pattern formed on a *bingata* specimen. The image reveals non-

transparent red colorants affixed to fiber surfaces; small defect regions, in which the colorants have peeled away or detached, appear white. From this we conclude that flower petals in the peony pattern were formed using a red pigment. Note that these red areas will appear as white regions in infrared photographs; visible/near-infrared reflection spectra, discussed below, reveal spectral shapes similar to that for *vermilion* in Figure 9A, while XRF spectra, also discussed below, indicate the presence of mercury (Hg). From these findings we conclude that the red coloration seen in the microscopy image on Figure 7 is due to particles of the pigment *vermilion* (HgS) affixed to fabric fibers.




Fig. 7 Microscopy image of red flower petals in peony pattern formed on *bingata* garment.

3. Visible/Near-infrared Reflection Spectra

3-1. Instruments⁶⁾

Our setup for measuring visible/near-infrared reflection spectra is shown in Figure 8. We use a two-branch fiber-optic cable, (R400-7-VIS-NIR, Ocean Optics USA), connected at one end to a light source (LS-1 tungsten halogen lamp, Ocean Optics USA) and at the other end to a miniature multichannel spectrometer (USB4000, Ocean Optics USA). For example, to measure visible/near-infrared reflection spectra for a *bingata* specimen, we place a protective sheet—from which is cut out a circular hole—on top of the specimen, positioned with its hole centered on a preselected measurement point, and then place a light-blocking box (390 g), into which is inserted the end of the optical fiber, on top of the sheet. Next, we look through the view portal of the light-blocking box and make fine adjustments to ensure that light from the optical fiber is irradiating the vicinity of the preselected measurement point; we maintain a distance of approximately 3 mm from the end of the cable to the specimen surface, yielding a spot diameter of 1 mm irradiated by light with wavelengths ranging from the visible to the near-infrared region (380-1000 nm). Light reflected from the specimen is then carried by the other branch of the optical fiber to the spectrometer, where it is subdivided by wavelength to yield a visible/near-infrared reflection spectrum—that is, a curve plotting spectral reflectance vs. wavelength from the visible to the near-infrared region—which is displayed on a PC. Note that this setup can also be used without the light-blocking box: with the optical-fiber tip positioned approximately 3 mm from the preselected measurement point, reflected light is immediately captured and indicated on the PC display, allowing non-destructive measurement of visible/near-infrared reflection spectra. We have used this approach to measure large numbers of woodblock prints in a short period of time, as discussed below in the section *Identifying blue colorants used in ukiyo-e woodblock prints*.




Fig. 8 Setup for measuring visible/near-infrared reflection spectra.

3-2. Information derived from visible/near-infrared reflection spectra

The visible/near-infrared reflection spectrum of a body plots the reflectance of the body (in units of percent) on the vertical axis vs. wavelength (in units of nanometers) on the horizontal axis. This yields a curve indicating how the body's reflectance varies with wavelength over the visible (380-700 nm) and near-infrared (700-1000 nm) regions. The visible-light portion of this curve contains information on the hue, brightness, and color intensity of colorants in the specimen, while the full curve spanning the visible and near-infrared regions contains information on the physical properties of the specimen. It is often possible to distinguish distinct colorants from the unique shapes of their spectra over this broad wavelength range, even if the colorant substances appear to be identical in color when viewed with the naked eye. Visible/near-infrared reflection spectra for some common colorants are shown in Figures 9-11.

Looking first at the spectrum of the pigment *vermilion* in Figure 9A, we see that the reflectance is low for short wavelengths, but begins to rise rapidly around 600 nm—a visible wavelength that the human eye recognizes as *red*—and remains high as the wavelength increases into the near-infrared region. Similarly, in the spectra of the iron-oxide-based pigments *Bengala* and *yellow ochre* in Figure 9B, we see characteristic absorption features near 655 nm in the visible region, and near 870 nm in the near-infrared region. Next, in the spectrum in Figure 10C for the pigment *smalt* we see a characteristic peak indicating a slight increase in reflectance near 565 nm in the visible region. In the spectrum in Figure 11B for the dye *lac* we see the reflectance increase rapidly near the 600 nm wavelength of *red*; and although this behavior is similar to that seen in the spectrum of the pigment *vermilion* (Figure 9A), the spectrum in Figure 11B is distinguished by a characteristic peak indicating a slight increase in reflectance near 545 nm in the visible region. Finally, in the spectrum in Figure 11C for the dye *indigo*, we see a characteristic spectral shape including a reflectance peak near the visible wavelength of 420 nm—interpreted by the human eye as *blue*—and a rapid increase in reflectance as we enter the near-infrared region beyond 700 nm. As these examples demonstrate, many varieties of colorants—including the red pigment *Bengala*, the yellow pigment *yellow ochre*, the blue pigment *smalt*, and the red dye *lac*—can be identified from the characteristic shapes of their reflection spectra.




Fig. 9 Visible/near-infrared reflection spectra of vermilion (A), iron-oxide-based pigments Bengala and yellow ochre (B), and malachite (C).




Fig. 10 Visible/near-infrared reflection spectra of Prussian blue (A), azurite (B), and smalt (C).




Fig. 11 Visible/near-infrared reflection spectra for sumi ink (A), lac (B), and indigo (C).

Identifying blue colorants used in ukiyo-e woodblock prints

The *ukiyo-e* woodblock prints produced during Japan's Edo period used a variety of substances to achieve blue coloration. *Commelinin* was the blue pigment of choice throughout the Meiwa and Kansei periods (1765-1800), but the years between the late Kansei period and the end of the Bunka period (1817) witnessed the discovery of *indigo*, which became increasingly popular after the start of the Bunsei period in 1818 and was the most common choice by the mid-Bunsei era, around the year 1824. The transition from *indigo* to *Prussian blue* began in 1830, the first year of the Tenpo era, and proceeded rapidly: by Tenpo 2 (1831) Prussian blue was used in fully 86% of all *ukiyo-e* prints, and by Tenpo 3 (1832) the figure was nearly 100%^{6,7)}. To investigate this historical evolution in the blue colorants of *ukiyo-e* prints, we analyzed a large set of *ukiyo-e* specimens, focusing in particular on *shibai-e* and *yakusha-e* prints^{*3} made to promote *kabuki* dramas and actors. Where necessary, we combined visible/near-infrared reflection spectra with other types of analyses discussed below—including XRF spectra to assist in identifying *Prussian blue* and 3DF spectra to assist in identifying *indigo*—but in many cases we found that visible/near-infrared reflection spectra—which are quick and easy to measure in our non-destructive setup—alone sufficed to distinguish *commelinin*, *indigo*, and *Prussian blue*. This is illustrated by the visible/near-infrared reflection spectrum of *commelinin* shown in Figure 12, which we readily observe to differ sufficiently from both the *Prussian blue* spectrum in Figure 10A and the *indigo* spectrum in Figure 11C to allow rapid identification of the blue colorant used in a given *ukiyo-e* specimen from the shape of measured reflection spectra.




Fig. 12 Visible/near-infrared reflection spectrum of *commelinin*.

Emergence of landscape images as subgenre of ukiyo-e prints

In addition to the insights described above, our use of visible/near-infrared reflection spectra as a tool for non-destructive analysis also achieves an interesting feat of historical detective work: We have established that Katsushika Hokusai's famous woodblock-print series "Thirty-Six Views of Mount Fuji," published around the first year of the Tenpo era, marks the first appearance of landscape images as a particular subgenre of *ukiyo-e* prints. As shown in Figure 13, each of the 36 prints in this series, labeled "Front views of Mount Fuji," makes use of both dark and light shades of blue to depict the sky and the sea in ways that lend the viewer a sense of closeness—or of distance. This is achieved via the method of *bokashi* printing⁸⁾, and by measuring visible/near-infrared reflection spectra for the sky and sea regions of each of the 36 prints we have confirmed that the spectrum of *Prussian blue* ($\text{KFe}[\text{Fe}(\text{CN})_6]_3 \cdot n\text{H}_2\text{O}$), shown in Figure 10A, is present in all cases, while XRF spectral measurements, discussed below, confirm the presence of iron Fe ⁹⁾. The attractive *bokashi* appearance evoked by this genre of landscape images is impossible to achieve without the ultra-miniature particles and bright blue coloration of *Prussian blue*, and thus we may say that the use of *Prussian blue* was a key enabling factor allowing the genre of landscape images to be introduced into the world of *ukiyo-e* woodblocks.

In the production process for *ukiyo-e* prints, the first step was to use a printing block known as a "key-block" (*omohan*⁹⁾) to print the title of the work and the signature of the illustrator, as well as pictures consisting of line drawings, and a sequence of "color-blocks" (*irohan*)—sculpted separately for each color in the image—was applied, each superposed on top of the last, to yield multicolored prints. In typical *ukiyo-e* processes, the *omohan* was used with *sumi* ink, but the

*3 *Shibai-e* and *yakusha-e* were the most frequently published types of *ukiyo-e* prints. These prints were rarely modified or reprinted after their initial publication, and the year of their publication can be determined from *kabuki* records or other historical resources. Indeed, as *kabuki* performances in this era occurred at two-month intervals, it is even possible to pinpoint the month of publication. As these prints were published as advertisements announcing upcoming performances, their production predates those performances.

"Thirty-Six Views" deviated from this pattern by using an *omohan* with blue *indigo*. Consequently, our measurements of visible/near-infrared reflection spectra for the title and signature regions and the outlines in the *Thirty-Six Views* series exhibit the characteristic shape of the *indigo* spectrum (Figure 11C)⁶.



Fig. 13 Katsushika Hokusai's "Thirty-Six Views of Mount Fuji (Front Views of Mount Fuji)".

Relationship between visible/near-infrared reflection spectra and infrared photographs

In regions of specimens containing pigments with high infrared transmittance, line-drawing sketches and other features that lie beneath surface layers are easily detected in infrared photographs¹⁾. To this end, the *reflectance* values plotted in the visible/near-infrared reflection spectrum of a colorant may simply be reinterpreted in terms of *transmittance* values^{*4}. This observation allows us to understand how a colorant appears in infrared photographs—specifically, the extent to which the colorant produces white regions or black regions in infrared photographs—based on reflectance data measured for the colorant in the near-infrared region. For example, if the *reflectance* of a colorant exceeds 60% in the infrared region, we may infer that the colorant will also exhibit high *transmittance* above 60%; consequently, most infrared radiation will pass through the colorant into the specimen, reflect from fiber surfaces or other features lying deeper within the specimen, then pass through the colorant a second time as it travels back out of the specimen to be captured by the camera, producing white regions in infrared photographs. In contrast, a colorant whose infrared *reflectance* is below 20% will also exhibit low *transmittance* below 20%; when capturing infrared photographs, such colorants will absorb most of the incident infrared light before it can travel deeper into the specimen and be reflected from fiber surfaces, and the absence of reflected light captured by the camera will produce black regions in infrared photographs. For intermediate colorants with near-infrared reflectance values around 40%, incident infrared radiation will be partially absorbed and partially transmitted; the transmitted portion will be reflected from fiber

*4 The software package we use to measure visible/near-infrared reflection spectra, from Ocean Optics USA, uses the same formulas to determine reflectance and transmittance. This allows us to interpret *reflectance* and *transmittance* interchangeably.

surfaces or other features beneath the colorants and will travel back out of the material—passing through the colorant a second time—to be captured by the camera, but the relatively low intensity of this reflected light will produce grey regions in infrared photographs.

Applying these observations to the colorants considered above, *vermilion* (Figure 9A), *lac* (Figure 11B), and *indigo* (Figure 11C) all exhibit a near-infrared reflectance of more than 60%, and thus appear as white regions in infrared photographs. In contrast, *malachite* (Figure 9C), *Prussian blue* (Figure 10A), *azurite* (Figure 10B), and *sumi ink* (Figure 11A) all exhibit a near-infrared reflectance of less than 20%, and thus appear as black regions in infrared photographs. Finally, the iron-oxide-based pigments *Bengala* and *yellow ochre* (Figure 9B), together with *smalt* (Figure 10C), have near-infrared reflectance values around 40%, and thus appear as grey regions in infrared photographs. This indicates how correlations between these two types of data—infrared photographs and visible/near-infrared reflection spectra—may be used to check the accuracy of analytical results.

Using visible/near-infrared reflection spectra to confirm the use of lac as a red colorant in bingata garments

One of the dyes commonly used in *bingata* garments is the red colorant *lac* (laccic acid). Figure 14 shows visible/near-infrared reflection spectra measured for red flower-petal regions in an iris pattern on a *bingata* specimen. In this spectrum we see a characteristic peak near 545 nm in the visible region, indicating a slight increase in reflectance, with the reflectance rising for wavelengths longer than 600 nm and remaining high into the near-infrared region. Note that the characteristic reflection peak near 545 nm matches the measured spectrum for *lac* shown in Figure 11B. However, the difference between the maximum and minimum measured reflectance values in Figure 14 is smaller than the corresponding gap in Figure 11B, suggesting that the colors for the specimen in Figure 14 are less intense (less brilliant) than those in the standard specimen^{*5} whose principal ingredient was *lac* and whose spectrum is shown in Figure 11B. The reason of this discrepancy is explained below in our discussion of XRF spectral measurements.

The visible/near-infrared reflection spectra that we measured for the red flower petals in the iris pattern also revealed the presence of lead (Pb), indicating that *white lead* ($2\text{PbCO}_3 \cdot \text{Pb(OH)}_2$) was used as a white colorant together with *lac* as a red colorant. Blending with *white lead* has the effect of increasing the brightness of the *lac*—that is, making it appear whiter—while decreasing the color intensity (brilliance), thus reducing the gap between the maximum and minimum reflectance values in Figure 14. Nonetheless, the appearance in the measured reflection spectra of a weak peak near 545 nm clearly indicates that the substance in question is in fact *lac*. As discussed below, the fluorescence fingerprint for *lac* obtained from 3DF spectra confirm that that red flower-petal regions of the iris pattern contain a blend of *white lead* and *lac* that has the effect of modifying the observed shade of red.




Fig. 14 Visible/near-infrared reflection spectrum measured for red region of *bingata* specimen.

*5 The visible/near-infrared reflection spectrum for lac in Figure 11B was measured for a standard specimen of Edo-era rouge-dyed cotton provided by Kunisuke Ueda of Uematsu Japanese Art Supply (Shibuya, Tokyo, Japan).

4. X-ray Fluorescence Spectra

4-1. Instruments¹⁰⁾

Figure 15 shows the setup we use for on-site measurement of XRF spectra. This apparatus does not use a fluorescent X-ray tube as an X-ray source; instead, X-rays are generated by an X-ray-emitting radioactive isotope. More specifically, the X-ray source (made by AET Technologies) consists of pellets of ²⁴¹Am arranged in a ring and hermetically sealed in a ceramic container, yielding a fully enclosed annular source of X-rays with an activity of 1.85 MBq. To measure XRF spectra, as shown in Figure 15, we position this source (labeled A in the figure) at a distance of approximately 5 mm from a preselected measurement point on a specimen, thus irradiating the specimen with X-rays emitted by ²⁴¹Am. (These X-rays have energies of 13.95 keV, 17.74 keV, and 59.54 keV.) The source produces a circumferential pattern of X-ray irradiation that induces XRF from elements present within a 12-mm-diameter circular region of the specimen centered on the measurement point. The energies of these fluorescent X-rays are measured by a semiconductor detector (Amptek XR-100CR Si-PIN, Be window 0.3 mil=7.62 μm, energy resolution 180-205 eV, labeled C in Figure 15). The signal from this detector is amplified by a preamp (Amptek PX2T/CR-type, labeled D), separated by a miniature multichannel pulse-height analyzer (Amptek PMCA-8000A, labeled E), and converted into XRF spectral data by a software package installed on a PC (F). Because this is an energy-dispersive apparatus, it is capable of simultaneous multi-element analysis. Shielding material (lead, labeled B) enclosing the upper surface, the perimeter, and the inner cavity region of the X-ray source (A) prevents X-rays emitted by the source from entering the window of the detector (C). The tip of the detector is affixed to the upper surface of the source across the shielding material.




Fig. 15 Setup for X-ray fluorescence measurements.




Fig. 16 Measuring XRF spectra for *bingata* specimen.



Fig. 17 Measuring X-ray fluorescence spectra for cultural property deemed a Japanese national treasure: square-base candle stand with green *Ryoku-Yu* glaze.




Fig. 18 Conceptual diagram illustrating typical X-ray fluorescence spectra.

To measure XRF spectra of *bingata* specimens using this apparatus, we begin by mounting the X-ray source A / detector C assembly on an adjustable arm designed for use as a stationary camera mount (Figure 16). We use the arm to position the source/detector assembly above the *bingata* specimen at a distance of approximately 5 mm from the preselected measurement point; the assembly will remain in this position throughout the non-contact measurement. The specimen is then irradiated by X-rays emitted from the ^{241}Am source, inducing elements in the vicinity of the preselected measurement point to emit fluorescent X-rays, whose energies are measured by the detector C.

Similarly, to measure XRF spectra for a piece of pottery or similar specimen, we fix the position of the source/detector assembly at a distance of approximately 5 mm from a flat portion of the specimen, as shown in Figure 17; again, X-rays emitted by the ^{241}Am source irradiate the specimen and induce the emission of fluorescent X-rays whose energies are measured by the detector. As the X-ray irradiation time increases, the measured XRF spectra increasingly come to resemble the spectrum shown in the conceptual diagram in Figure 18. The length of time required to obtain accurate spectra for a given element depends on how much of the element is present in the specimen. The larger the volume of the element contained in the specimen, the more rapidly the measured spectral curve will take shape^{*6}, and thus the shorter the measurement time required to yield clear spectra. Typical measurement times [characterized by *Live Time* (LT), the time required to accumulate an effective number of X-ray counts] depend on the volume of the given element contained in the specimen, but are on the order of 180–600 s (3–10 min) if the element in question is abundant in the specimen, or 1200–1800 s (20–30 min) if the element is scarce in the specimen. For specimens containing particularly low volumes of a given element, the time required to measure clear XRF spectra may exceed 3600 s (1 h). These long measurement times are a consequence of the weak X-ray source in our apparatus: the ^{241}Am source has an activity of just 1.85 MBq. However, this source was purchased through the Japan Radioisotope Association, with certification from AET Technologies USA, in 1999, prior to the 2005 enactment of the "Act on the Regulation of Radioisotopes, etc." that is currently in effect in Japan. Because the activity of the source does not exceed the limit (3.7 MBq) set by the previous regulatory framework, at present the source can be used with no need to notify, or receive permission from, relevant oversight bodies—assuming, of course, that it is used and managed properly. Moreover, because the source is hermetically sealed and has an intensity below 3.7 MBq, there is no need to worry about exceeding the regulated maximum dose of external radioactivity (1.3 mSv over 3 months); consequently, we can use the source for non-destructive on-site analyses at museums and other facilities without needing to establish controlled X-ray areas. Also, although the activity of this source was 1.85 MBq in the year it was constructed (1998), this activity has decreased over the ensuing 25 years. Specifically, based on the 432-year half-life of ^{241}Am , we calculate that the activity of the source has declined by 5.8% since its construction, placing its current activity at 1.74 MBq. We consider this to be more than adequate for use in X-ray fluorescence analysis. However, if we ever purchase a new X-ray source, we will need to receive permission to use it in accordance with current laws and regulations.

4-2. Information derived from XRF spectra

As shown in the conceptual diagram in Figure 18, XRF spectra are obtained by plotting the intensity (counts) on the vertical axis against the energy (keV) on the horizontal axis. These curves describe the energy and intensity of fluorescent X-rays emitted by the constituent elements of the specimen. The XRF *energy* is indicated by the horizontal-axis position of peaks in the spectrum, while the *intensity* is indicated by the peak height.

XRF spectra form doublets and triplets

The X-rays emitted by the ^{241}Am source, which have energies of 13.95 keV, 17.74 keV, and 59.54 keV, induce different types of fluorescent X-ray emissions from different elements. Elements in column 4 of the periodic table—including calcium (Ca), titanium (Ti), chromium (Cr), manganese (Mn), iron (Fe), cobalt (Co), nickel (Ni), copper (Cu), zinc (Zn), and arsenic (As)—emit fluorescent X-rays with the characteristic energy of the $K\alpha$ line. For column-5 elements such as silver (Ag), cadmium (Cd), and tin (Sn), fluorescent X-ray emissions have the characteristic energy of the $K\beta$ line.

^{*6} When an X-ray tube is used as an X-ray source, measurement times may be shortened by increasing the tube voltage (kV) to increase the energy of emitted X-rays or by increasing the beam current (μA) to increase the intensity (number) of emitted X-rays, although the details depend on the elemental composition of the measurement specimen.

Meanwhile, column-6 elements—such as barium (Ba), tungsten (W), platinum (Pt), gold (Au), mercury (Hg), and lead (Pb)—emit fluorescent X-rays with the energies of the $L\alpha$ line, the $L\beta$ line, and the $L\gamma$ line. Consequently, fluorescent X-ray emissions from a single element may be organized into groups: the $K\alpha$ -line and $K\beta$ -line spectra in Figure 18 form a *doublet*, while the $L\alpha$ -line, $L\beta$ -line, and $L\gamma$ -line spectra form a *triplet*. The energies of fluorescent X-ray emissions from a given element are characteristic of that element, and reference manuals for X-ray analysis contain detailed tables listing characteristic $K\alpha$ -line and $K\beta$ -line energies, or characteristic $L\alpha$ -line, $L\beta$ -line, and $L\gamma$ -line energies, for the various elements of the periodic table, ordered by atomic number¹¹⁾. Using these tables, measured XRF spectra can be compared to tabulated data: if the horizontal-axis position of a peak in a measured spectrum matches the characteristic fluorescent X-ray energy tabulated for a given element, it can be concluded that this element is present in the sample. Table 2 lists the energies of various fluorescent X-ray emissions from the constituent elements of common pigments and metals. If the measured energy of a fluorescent X-ray emission from a specimen (that is, the energy read off from the horizontal-axis position in the XRF spectrum) matches the energy tabulated for a given element, the pigments or metals having that element as a constituent can be identified.

Also, for a specimen irradiated by X-rays of a given intensity, the intensity of fluorescent X-ray emissions from an element is proportional to the volume (number of atoms) of that element contained in the specimen. We can make use of this observation as follows: We prepare standard specimens containing different known concentrations of a particular element, measure the XRF intensity for each of these specimens, then plot the intensity for each specimen versus the concentration of the given element in that specimen. The resulting plot, known as a *calibration curve*, can then be used to quantify the concentration (volume) of the given element in an unknown specimen based on the XRF intensity for that specimen.

Table 2 Energies of fluorescent X-rays emitted by constituent elements of pigments and metals

Color	Name	Constituent element (chemical formula)	Energy of fluorescent X-ray emission from pigment or metal		
			α line	β line	γ line
White	Chalk	CaCO_3	Ca $K\alpha$ 3.691 keV	Ca $K\beta$ 4.013 keV	—
	White lead	$2\text{PbCO}_3 \cdot \text{Pb(OH)}_2$	Pb $L\alpha$ 10.550 keV	Pb $L\beta$ 12.612 keV	Pb $L\gamma$ 14.764 keV
Red	Vermilion	HgS	Hg $L\alpha$ 9.987 keV	Hg $L\beta$ 11.821 keV	Hg $L\gamma$ 13.830 keV
	Red lead	Pb_3O_4	Pb $L\alpha$ 10.550 keV	Pb $L\beta$ 12.612 keV	Pb $L\gamma$ 14.764 keV
	Bengala	Fe_2O_3	Fe $K\alpha$ 6.400 keV	Fe $K\beta$ 7.058 keV	—
Yellow	Orpiment	As_2S_3	As $K\alpha$ 10.532 keV	As $K\beta$ 11.726 keV	—
	Yellow ochre	$\text{Fe}_2\text{O}_3 \cdot n\text{H}_2\text{O}$	Fe $K\alpha$ 6.400 keV	Fe $K\beta$ 7.058 keV	—
Green	Malachite	$\text{CuCO}_3 \cdot \text{Cu(OH)}_2$	Cu $K\alpha$ 8.041 keV	Cu $K\beta$ 8.905 keV	—
	Emerald green	$\text{Cu}(\text{C}_2\text{H}_3\text{O}_2)_2 \cdot 3\text{Cu}(\text{AsO}_2)_2$	Cu $K\alpha$ 8.041 keV As $K\alpha$ 10.532 keV	Cu $K\beta$ 8.905 keV As $K\beta$ 11.726 keV	— —
Blue	Prussian blue	$\text{KFe}[\text{Fe}(\text{CN})_6]_3 \cdot n\text{H}_2\text{O}$	Fe $K\alpha$ 6.400 keV	Fe $K\beta$ 7.058 keV	—
	Azurite	$\text{Cu}_3(\text{CO}_3)_2(\text{OH})_2$	Cu $K\alpha$ 8.041 keV	Cu $K\beta$ 8.905 keV	—
	Smalt	$[\text{Co}, \text{Ni}]\text{As}_{3-2}, \text{Fe} + \text{SiO}_2, \text{K}_2\text{O}, \text{CaO}$	Fe $K\alpha$ 6.400 keV	Fe $K\beta$ 7.058 keV	—
			Co $K\alpha$ 6.925 keV	Co $K\beta$ 7.649 keV	—
		Ni $K\alpha$ 7.472 keV	Ni $K\beta$ 8.265 keV	—	
		As $K\alpha$ 10.532 keV	As $K\beta$ 11.726 keV	—	
Black	Sumi	C	C $K\alpha$ 0.277 keV	XRF is undetectable	
Gold	Gold (foil/slurry)	Au	Au $L\alpha$ 9.712 keV	Au $L\beta$ 11.440 keV	Au $L\gamma$ 13.382 keV
	Brass (foil/shavings)	Cu:Zn = 70-90:30-10 wt%	Cu $K\alpha$ 8.041 keV Zn $K\alpha$ 8.631 keV	Cu $K\beta$ 8.905 keV Zn $K\beta$ 9.572 keV	— —
Silver white	Silver (foil/slurry)	Ag	Ag $K\alpha$ 22.105 keV	Ag $K\beta$ 24.942 keV	—
	Tin	Sn	Sn $K\alpha$ 25.295 keV	Sn $K\beta$ 28.486 keV	—

(Intensity ratio) $K\alpha/K\beta$: approximately 10:1; $L\alpha/L\beta$: between 10:6 and 10:5; $L\beta/L\gamma$: approximately 10:1

Intensity ratio: The relative intensity of $K\alpha$ and $K\beta$ lines—and of $L\alpha$, $L\beta$, and $L\gamma$ lines

$K\alpha$ -line and $K\beta$ -line emissions, which form a doublet, differ in intensity by roughly a factor of 10, i.e., the intensity ratio is approximately $K\alpha:K\beta = 10:1$. For $L\alpha$ -line and $L\beta$ -line emissions, which are part of a triplet, the intensity ratio is typically adjusted by varying the voltage of the X-ray tube used to generate X-rays to a value somewhere between $L\alpha:L\beta = 10:10$ and $L\alpha:L\beta = 10:5$ ¹²⁾. However, this procedure is not available in our apparatus, which uses a ²⁴¹Am source (activity 1.85 MBq) in place of an X-ray tube. Consequently, while the $K\alpha$ -line/ $K\beta$ -line intensity ratio in our experiments assumes its usual value of roughly 10:1, the $L\alpha$ -line/ $L\beta$ -line intensity ratio is roughly 10:6 or 10:5. Also, if $L\alpha$ -line and $L\beta$ -line emissions are observed, then $L\gamma$ -line emissions are also observed, although with low intensity: the $L\beta$ -line/ $L\gamma$ -line intensity ratio is roughly 10:1.

As an example, Figures 19a and 19b show XRF spectra measured by our apparatus for the green pigment *emerald green* ($\text{Cu}(\text{C}_2\text{H}_3\text{O}_2)_2 \cdot 3\text{Cu}(\text{AsO}_2)_2$) and for the white pigment *white lead* ($2\text{PbCO}_3 \cdot \text{Pb}(\text{OH})_2$). Because our apparatus is an energy-dispersive fluorescent X-ray analyzer based on a semiconductor detector, it can simultaneously detect fluorescent X-ray emissions from multiple elements. Emerald green contains Cu and As as constituent elements, and both elements emit fluorescent X-rays, whose energies we detect. As shown in Figure 19a, from Cu we obtain spectrum A, corresponding to the 8.04 keV Cu- $K\alpha$ line, and spectrum B, corresponding to the 8.92 keV Cu- $K\beta$ line. For As we obtain spectrum C, corresponding to the 10.54 keV As- $K\alpha$ line, and spectrum D, corresponding to the 11.75 keV As- $K\beta$ line. In both cases, the $K\alpha$ -line/ $K\beta$ -line intensity ratios are approximately 10:1. On the other hand, white lead contains Pb as a constituent element; as shown in Figure 19b, we obtain spectrum A, corresponding to the 10.55 keV Pb- $L\alpha$ line, spectrum B, corresponding to the 12.64 keV Pb- $L\beta$ line, and spectrum C, corresponding to the 14.76 keV Pb- $L\gamma$ line. The $L\alpha$ -line/ $L\beta$ -line intensity ratio is roughly 10:6 or 10:5, while the $L\beta$ -line/ $L\gamma$ -line intensity ratio is roughly 10:1.

From the above discussion, the intensity ratio for $K\alpha$ -line and $K\beta$ -line emissions, which form a doublet, differs from that for $L\alpha$ -line and $L\beta$ -line emissions, which form a triplet. This observation may be used to determine whether peaks in measured spectra correspond to $K\alpha$ -line and $K\beta$ -line emissions or to $L\alpha$ -line and $L\beta$ -line emissions*⁷. Note that, in XRF spectra, $K\beta$ -line energies are higher than $K\alpha$ -line energies, so that the horizontal-axis positions of $K\beta$ -line peaks lie to the right of the positions of $K\alpha$ -line peaks. Similarly, $L\beta$ -line energies are higher than $L\alpha$ -line energies (so $L\beta$ -line peaks lie to the right of $L\alpha$ -line peaks) and $L\gamma$ -line energies are higher than $L\beta$ -line energies (so $L\gamma$ -line peaks lie to the right of $L\beta$ -line peaks). This observation is useful for, e.g., deciphering whether features appearing in measured XRF spectra are attributed to $K\alpha$ and $K\beta$ lines for a given element or to $L\alpha$ and $L\beta$ lines for that element.

Overlapping features in XRF spectra

Our apparatus is an energy-dispersive X-ray fluorescence analyzer based on a semiconductor detector with an energy resolution of 180-205 eV; this means it is capable of simultaneously detecting the energies of fluorescent X-ray emissions from multiple elements. However, if the energies of X-rays emitted from distinct elements are nearly equal—specifically, if they differ by less than our detector resolution of 205 eV = 0.205 keV—then the two spectra overlap in our setup. For example, as noted above, the green pigment *emerald green* includes both Cu and As among its primary constituent elements, and thus fluorescent X-rays emitted simultaneously from both of these elements appear in measured XRF spectra. If *emerald green* was the only colorant used in a specimen, then, as shown in Figure 19a, the $K\alpha$ -line and $K\beta$ -line emissions from Cu and As are not close in energy, so there would be no overlap and we would have no difficulty measuring both spectra simultaneously. However, in practice emerald green is used together with other colorants; for example, the white pigment *white lead* is sometimes added to increase brightness and yield a brighter hue. For a specimen in which the white pigment *white lead* is combined with dark-green *emerald green* to yield a brighter shade of green, XRF measurements simultaneously detect the energies of X-ray emissions from *emerald green* (Figure 19a) and from *white lead* (Figure 19b), yielding the spectrum shown in Figure 19c, in which the As- $K\alpha$ line emission from *emerald green* overlaps with the Pb- $L\alpha$ line emission from *white lead*. Note that this phenomenon is not limited to

*⁷ The elemental analysis software installed on the PC attached to our X-ray fluorescence analyzer can automatically detect peaks in XRF spectra and determine the corresponding elements from their energies. Nonetheless, an understanding of how $K\alpha$ / $K\beta$ -line intensity ratios differ from $L\alpha$ / $L\beta$ -line intensity ratios is important not only for confirming the accuracy of results reported by this software but also for interpreting the results of XRF spectra published in the literature.

the specific case of *emerald green* mixed with *white lead*, but would arise for any specimen in which As and Pb are simultaneously present, due to the similar energies of the As-K α line (10.532 keV) and the Pb-L α line (10.550 keV).

Similarly, the yellow pigment *orpiment* (As₂S₃) is sometimes combined with the red pigment *vermilion* (HgS) to yield orange coloration. When these two colorants are simultaneously present, we have K α -line (10.532 keV) and K β -line (11.726 keV) emissions from As, a primary constituent of orpiment, together with L α -line (9.987 keV) and L β -line (11.821 keV) emissions from Hg, a primary constituent of vermilion; this results in an overlap between the As-K β line (11.726 keV) and the Hg-L β line (11.821 keV). The same difficulty would arise even if *orpiment* and *vermilion* were present in separate but neighboring regions of a specimen.

As an example, consider the *bingata* garment pattern shown in Figure 20. For this specimen, XRF spectra measured in yellow regions of the garment contained three spectral peaks: peak A, near 9.98 keV; peak B, near 10.54 keV; and peak C, near 11.73 keV. A first observation is that the intensities for peaks B and C differ by a factor of approximately 10, i.e., the B:C intensity ratio is 10:1. This suggests that the higher-intensity peak B is a K α -line peak, while the lower-intensity peak C is a K β -line peak forming a doublet with peak B. Next, referring to Table 2, we ask: What pigment would emit both a K α line corresponding to an energy of 10.54 keV and a K β line corresponding to an energy of 11.73 keV? As the measured specimen region is primarily yellow in color, we deduce that peak B corresponds to the K α line (10.532 keV) for As—a principal constituent of *orpiment*—while peak C corresponds to the K β line (11.726 keV) for As. Based on these findings, we conclude that the yellow pigment here is *orpiment*.

The only remaining mystery is the origin of peak A at 9.98 keV. Again referring to Table 2, we ask: Which element emits fluorescent X-rays corresponding to this energy? A quick search reveals the answer: Hg, whose L α line energy is 9.987 keV, and whose L β line—which forms a doublet with the L α line—lies at 11.821 keV, and thus overlaps in our measured spectrum with peak C, i.e., with the As-K β line at 11.726 keV.

As this example demonstrates, a thorough understanding of the phenomenon of overlapping spectra is essential for correctly interpreting spectral measurements—and for deciphering the mystery of "Where did this peak come from?!"




Fig. 19a XRF spectrum of *emerald green*.




Fig. 19b XRF spectrum of *white lead*.




Fig. 19c Measured XRF spectrum for a specimen containing both *emerald green* and *white lead*, illustrating overlapping spectral features.




Fig. 20 XRF spectrum (right) measured for the yellow-colored region of the specimen shown on the left, which contains the yellow pigment *orpiment* and lies adjacent to a red-colored specimen region containing the red pigment *vermilion*.

Analyzing measured XRF spectra for the pigment *smalt*¹⁰⁾

In some cases, an understanding of the historical trajectory that led to the introduction of a new pigment can serve to facilitate the challenge of simultaneously analyzing multiple constituent elements in specimens. As one example, Figure 21 shows a cultural-property specimen for which our measured XRF spectrum successfully identified the presence of the classical blue pigment *smalt*. The specimen in this case is an artwork entitled *Images of Rashomon* (dimensions: 134.4 cm tall × 102.3 cm wide) that was painted on a wooden prayer board and left as an offering to the Asakura Shrine in Fukui Prefecture in Japan in the second year of the Tenna era (1682). The painting depicts the samurai warrior Watanabe no Tsuna, and we measured XRF spectra for a blue region painted on the dome of the warrior's helmet.

Smalt, the oldest man-made Co-based pigment, is made by adding Co ore ($[\text{Co},\text{Ni}]\text{As}_{3,2},\text{Fe}$) to molten silicate glass ($\text{SiO}_2,\text{K}_2\text{O},\text{CaO}$) to yield a blue-colored glass, which is then ground into a powder^{13,14)}. The Co content of the ore is the source of the blue coloration, but the ore also contains Fe, Ni, and As. Armed with this knowledge, we can determine the origins of various elements whose presence is indicated by the XRF spectrum in Figure 21. Ca originates from silicate glass, while Fe, Co, Ni, and As originate from Co ore. Although there is a paint named *smalt* that is commercially available today, it is made with Co as the sole source of blue coloration and is distinct from the classical pigment known as *smalt*. The spectral features marked “scatter” in Figure 21 are due to X-rays emitted from the ²⁴¹Am source (with energies of 13.95 keV and 17.74 keV) that undergo scattering processes and are erroneously captured by the detector; they represent an instrumental artifact and do not indicate peaks associated with elements present in the sample.

Analyzing XRF spectra is a technique that works well for inorganic substances such as pigments and metals, but cannot be used for organic substances such as dyes. The reason for this is that the fluorescent X-rays emitted from C, a component of *sumi* ink, have extremely low energies—just 0.277 keV—and cannot be detected by any X-ray fluorescence analyzer currently available. This also means that XRF analysis cannot be used to study dyes or other organic substances structured primarily around C bonds.




Fig. 21 Left: *Images of Rashomon*, an artwork presented as an offering to a shrine in Japan's Edo period [Tenna 2 (1682)]. Right: XRF spectrum measured for a blue-pigment-filled subregion of the dome of the warrior's helmet.

Combining XRF spectra and visible/near-infrared reflection spectra to unravel the mystery of green *tsuikin*¹⁵⁾

Tsuikin is a method for transferring pictorial patterns onto lacquerware. In this technique, a pigment is uniformly blended with lacquer to yield a paste known as *tsuikin-mochi*, which is stretched into a thin sheet, cut into shapes appropriate for the pattern, and affixed to lacquerware items. The cultural property shown in the lower left panel of Figure 22—a tray made from black-lacquered wood that is recognized as a Japanese national treasure—was decorated by *tsuikin*. The other panels in Figure 22 show the results of various non-destructive analyses of this specimen: an XRF spectrum measured for a yellow-colored region of the *tsuikin* decoration (upper left), an XRF spectrum measured for a green-colored *tsuikin* region (upper right), and a visible/near-infrared reflection spectrum measured for the same green-colored *tsuikin* region (lower right).

Looking first at the XRF spectrum for the yellow *tsuikin* region (upper left), we see two peaks—peaks A and B—at energies of 10.54 keV and 11.77 keV, respectively, with the intensity of peak A being approximately ten times greater than that of peak B (intensity ratio A:B = 10:1). This suggests that peaks A and B respectively indicate K α and K β lines for some element. Consulting Table 2, we find that this element is As, whose K α and K β lines have energies of 10.532 keV and 11.726 keV, respectively. Hence we conclude that *orpiment* (As₂S₃) was used to produce the yellow color of this specimen region. In other words, the yellow-colored pictorial pattern on the tray was produced by blending *orpiment* uniformly with lacquer to yield a *tsuikin-mochi* paste.

Looking next at the XRF spectrum for the green *tsuikin* region (upper right), we see two peaks (peaks B and C) with intensities differing by roughly a factor of ten (intensity ratio B:C = 10:1), again suggesting K α and K β lines for some element. Matching the positions of the peaks to data in Table 2, we conclude that peaks B and C are the K α and K β lines for the yellow pigment *orpiment*, respectively. This spectrum also exhibits a peak (peak A) at 6.39 keV, which we identify as the Fe-K α line (energy: 6.400 keV). We attribute the presence of Fe in this region to the black lacquer¹⁶⁾ used to coat the tray, which was prepared by adding Fe powder or pellets to liquid lacquer; the Fe-K α line seen here is not due to Fe present in the yellow *tsuikin* region. This interpretation is also supported by the fact that the Fe-K α and Fe-K β lines both appear prominently in the XRF spectrum measured for a region of the specimen containing only black lacquer, with no *tsuikin* colorations.

Thus, we now face this question: Which blue colorant was used in combination with yellow-colored *orpiment* to yield the green coloration? It has been suggested that *malachite* may have been used in some cases to yield green lacquer¹⁶⁾. However, we see here neither the Cu-K α (8.041 keV) or Cu-K β (8.905 keV) lines that we would expect to find for copper-containing malachite (CuCO₃ · Cu(OH)₂). According to Ref. 16, "Looking at historical records from long-past eras, we find statements indicating that the *indigo* used was the same *indigo* used to make *indigo* dye". But Ref. 17 says that "although *indigo* existed at that time, it is unclear whether or not it was used in mixtures with *orpiment*." To clarify this murky situation, we switch to a different method of analysis: we measure visible/near-infrared reflection spectra for the same green *tsuikin* region, yielding the reflection spectrum shown at the lower right in Figure 22. Because this spectrum is measured for a *green*-colored region, it exhibits no peaks in the 400-500 nm wavelength interval corresponding to blue and violet. On the other hand, the shape of the spectral curve—with the reflectance increasing in the near-infrared region beyond 700 nm—resembles the spectrum of *indigo* in Figure 11C. The *indigo* in *indigo* dye, when stored in the clay or ceramic jars commonly used as storage containers, tends to produce water-soluble leuco salts (green-colored), which are oxidized in air to yield a non-soluble crystalline blue substance used as a pigment¹⁸⁾. Based on these findings, we conclude that green *tsuikin* regions are produced by blending the yellow colorant *orpiment* with the blue colorant *indigo* to yield a green colorant.

A cultural property that does not involve *tsuikin* but is nonetheless relevant here is a *Gold-inlaid green-lacquered wood tray* held by the Okinawa Churashima Foundation. A *green* colorant has been applied to the four corners of this tray, and the XRF spectra and visible/near-infrared reflection spectra measured for these green regions are similar to those in Figure 22. Thus we conclude that, for this specimen as well, *orpiment* was blended with *indigo* to yield a *green lacquer* that was applied to the specimen¹⁹⁾.




Fig. 22 Lower left: A tray recognized as a Japanese national treasure made from black-lacquered wood and decorated by *tsuikin*. Upper left: XRF spectrum measured for a yellow region of the *tsuikin* decoration. Upper right: XRF spectrum measured for a green region of the *tsuikin* decoration. Lower right: Visible/near-infrared reflection spectrum measured for the same green region of the *tsuikin* decoration.

X-ray fluorescence analysis of (glazed) earthenware pottery²⁰⁾

The Japanese national treasure in Figure 17, a square-base candle stand with green *Ryoku-Yu* glaze, is one example of *Ryoku-Yu* glazed pottery. According to Ref. 21, *Ryoku-Yu* glaze—also known as *Ao-Oribe-Yu* ("green glaze")—is prepared by "adding copper oxide (Cu_2O) to transparent wood-ash glaze or lime glaze, then firing; in some cases, impurities such as chalcantite (copper sulfate $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$), shavings of brass [an alloy of copper (Cu) and zinc (Zn)], or shavings of copper may be added to the copper component to add color". Figure 23 shows the XRF spectrum measured for the skin layer of the *Ryoku-Yu* glazed candle stand. We identify six peaks in this spectrum: peak A (3.69 keV), peak B (6.41 keV), peak C (8.00 keV), peak D (8.67 keV), peak E (8.88 keV), and peak F (9.54 keV). The highest-intensity peak is peak C, which corresponds to the Cu-K α line (8.041 keV) and reflects the presence of Cu as the principal colorant in the fired glaze. Peak E corresponds to the Cu-K β line (8.905 keV). The next highest-intensity peak is peak D, at an energy of 8.67 keV, which is roughly ten times more intense than peak F (9.54 keV) (intensity ratio D:F = 10:1). This suggests that peaks D and F correspond to the Zn-K α (8.631 keV) and Zn-K β (9.572 keV) lines. Thus we conclude that Cu and Zn coexist in the portions of this green glaze that add coloring upon firing, and that the glaze consists of brass shavings. Finally, peaks A (3.69 keV) and B (6.41 keV) correspond to the Ca-K α line (3.691 keV) and the Fe-K α line (6.400 keV); Ca and Fe are present in the clay materials and lime-ash glazes (calcium carbonate CaCO_3) used as raw materials.




Fig. 23 XRF spectrum measured for the *Ryoku-Yu* glaze used in the Japanese national treasure in Figure 17.

5. Three-dimensional Fluorescence (3DF) Spectra

5-1. Instruments

At present, we use two similar analytical instruments to measure 3DF spectra: For on-site analyses at museums or other remote locations, we use a Hitachi F-2500 fluorescence spectrophotometer (weight: 35 kg) equipped with an optical fiber, while for analyses in our home laboratory we use a Hitachi F-7000 fluorescence spectrophotometer (weight: 41 kg) permanently installed in a fixed location and similarly equipped with an optical fiber. In some cases, the cultural property we wish to analyze cannot be moved, and in some cases we must analyze not only the cultural property but also various aspects of the environment in which it is housed, such as the vessels in which it is stored. In such cases, we transport our F-2500 system to the remote site—either by hiring a shipping company capable of handling delicate instruments, or by hand-delivering the system ourselves in a private vehicle. In the latter case the spectrophotometer and the optical fiber are packaged in a custom-designed protective metal carrying case. Needless to say, among the accessory items transported with the system for use at remote sites are samples of rhodamine B and other reagents, as well as secondary standard light sources, to facilitate the steps required to prepare the instrument for use at the remote site, including measuring instrument response functions (and applying spectral corrections) and measuring instrument properties in the long-wavelength fluorescence region. Figure 24 shows a simplified schematic diagram of the optical fibers we use for the F-2500 and F-7000; thick lines in this diagram represent the shapes of fibers. Initially, we adopted a fiber configuration in which the excitation side used 48 fiber cores bundled into a single cable, while the fluorescence side used 32 fiber cores arranged in a double-bundle configuration designed to wrap around the perimeter of the excitation cable. However, we now use a different configuration in which both the excitation-side fibers (total of 48 cores) and the fluorescence-side fibers (total of 48 cores) are dispersed and randomly bundled. To evaluate the performance of these two fiber configurations, we measured 3DF spectra (as shown in the center of Figure 26) under identical measurement conditions and compared the two fluorescent intensity values; the results demonstrated that the measured intensity was approximately 2.5 times greater for the random-bundle configuration than for the double-bundle configuration—even with an L39 filter installed on the fluorescence side to eliminate secondary light producing scattering artifacts in 3DF spectra—thereby allowing the detector to capture more of the fluorescent light emitted by measurement specimens²²). However, with the tip of the optical fiber positioned at a distance of approximately 5 mm from the specimen, the spot size—that is, the diameter of the region irradiated by the fiber—is approximately 3 mm for the double-bundled configuration, but expands to approximately 5 mm for the random configuration.

To measure 3DF spectra at remote sites using the F-2500 with this randomly-bundled optical fiber, we use a setup like that shown in Figure 25. We begin by covering the measurement sample (a *bingata* garment in the case of Figure 25) with a protective sheet from which we have cut out a circular hole approximately 30 mm in diameter. The sheet is positioned such that the center of this hole lies at a preselected measurement point on the specimen. Then we place the light-blocking box, into which the tip of the optical fiber is inserted, on top of the protective sheet (the mass of the box is 390 g), open its view portal, and look through this portal while adjusting the positioning to ensure that the tip of the optical fiber lies above the measurement point at a distance of approximately 5 mm from the specimen surface. Upon completing this preparation we begin the measurement. When adjusting the position of the light-blocking box after inserting the optical-fiber tip, it is desirable to pass visible light (e.g., 550 nm) through the fiber.

When we first began to pursue the research program described in this article, our analytical instrument was a Hitachi F-4010 fluorescence spectrophotometer. Natural dyes derived from animal or plant matter, whose molecular structure involves a significant prevalence of covalent double bonds²³, exhibit characteristic excitation and fluorescence wavelengths, and thus it should be possible to use *fluorescence fingerprints*—contour-plots of 3D fluorescent spectra—to directly identify the dyes responsible for lending coloration to dyed fabrics. This was the reasoning that motivated the launch of our research program in the early 1990s²⁴. Later, we replaced the F-4010 with our current F-4500 system, developed the optical-fiber approach described above, and expanded the focus of our research from dyed fabrics to *ukiyo-e* prints^{25,26}. It was at this juncture that we began developing our procedure for making on-site measurements with the F-2500, and since then we have performed non-destructive analyses of *ukiyo-e* prints at museums throughout Japan and at the Museum of Fine Arts, Boston, USA²⁷. At present, our research is focused on using Hitachi's F-7000

system for non-destructive analysis of lacquers used to make lacquerware.

From the moment we began this research program in the 1990s, and continuing through to the present day, there has been a continual expansion of the capabilities of fluorescence spectrophotometers. Among the advances that we find most valuable are the dramatic reduction in measurement times (facilitated by higher scan speeds), the significantly longer lifetimes of xenon-lamp light sources, and the proliferation of data-analysis functionality. We are also considering Hitachi's F-2700 compact fluorescence spectrophotometer (weight: 41 kg) as a possible replacement for the F-2500 system we currently use for on-site measurements.




Fig. 24 Optical fiber configuration.



Fig. 25 Measuring 3DF spectra of *bingata* specimen.

5-2. Information derived from 3DF spectra

For a specimen of the dye *safflower* (coloring matter: carthamin), Figure 26 shows the measured 3DF spectrum and its contour-plot representation or *fluorescence fingerprint*. 3DF spectra are measured by successively irradiating the specimen at a sequence of excitation wavelengths and measuring the fluorescence spectrum at each excitation wavelength. To perform such a measurement, we specify values for various parameters—the excitation wavelength range (typically 250–600 nm), the wavelength interval between successive excitation wavelengths (typically 5 nm), and the wavelength range over which to measure fluorescence spectra for each excitation wavelength (typically 300–700 nm)—then simply initiate the measurement. If we watch in real time as the measurement proceeds, we see a gradual sequence of variations in the fluorescence spectra measured for each excitation wavelength. For typical values of the scan parameters—excitation wavelength ranging from 250 to 600 nm in 5 nm steps—a total of 71 separate fluorescence spectra are measured (Figure 26, left). After the full measurement is complete, these 71 spectral curves are combined to yield a three-dimensional plot (Figure 26, center) depicting the variation of fluorescence intensity with respect to both excitation wavelength (Ex) and fluorescence wavelength (Em). This plot is known as a *three-dimensional fluorescence (3DF) spectrum*²⁸⁾. Note that this plot is reminiscent of bird's-eye view of a topological map depicting a mountain range; in the same way that a three-dimensional topological map can be converted into a two-dimensional contour map—with the altitude information conveyed by three-dimensional features in the former now encoded by contour lines in the latter—we can convert the 3DF spectral plot (Figure 26, center) into a contour-plot (Figure 26, right). The result, known as a *fluorescence fingerprint*, furnishes a unique signature to the measured specimen, which may be used, for example, to identify the dyes used to color cloth fibers or to color *ukiyo-e* prints. The contour-plot shown at the right in Figure 26 is a fluorescence fingerprint measured for the pigment *carthamin* in the dye *safflower*, whose molecular structure is shown in the lower central portions of the figure. Because this fluorescence fingerprint is derived from the particular molecular structure of this dye, any other dye—or any other molecular structure—would have a distinct fingerprint.

In forensic science, human fingerprints are used to identify specific individuals. For example, in Japan, fingerprints left at crime scenes (latent fingerprints) are compared against two databases maintained by the Fingerprint Center of

Japan's National Police Agency: one database stores the fingerprints of more than 8 million convicted criminals, and a second database contains billions of latent fingerprints left at past crime scenes. Fingerprint comparisons are automated by the Automated Fingerprint Identification System (FAIS), with the final determination as to whether or not a given fingerprint matches a fingerprint in the database being made by an expert forensic analyst^{29,30}.

The use of human fingerprints to identify crime suspects furnishes a partial analog to the use of fluorescence fingerprints to identify dyes in cultural properties. The fluorescence fingerprints measured for a cultural property, which correspond to latent fingerprints left at a crime scene, are compared against fluorescence fingerprints measured for laboratory-prepared standard specimens; if a match is found, we may conclude that the dye used in the cultural property is identical to the dye used to prepare the standard specimen. To construct a reference database of fluorescence fingerprints for known substances, we perform measurements of silk fabric specimens dyed with common natural dyes and of *washi* paper specimens colorized by pressing natural coloring matters into the paper.

Figure 27 shows fluorescence fingerprints measured for specimens of silk fabric colorized by single dyes and for *washi* paper colorized by single dyes. When a desired colorization effect cannot be achieved using any one single dye, multiple dyes may be combined; for example, in ancient Japan the colors *scarlet* and *green* were each realized by over-dyeing using two different dyes³¹. Figure 28 shows a fluorescence fingerprint measured for a standard silk-fabric specimen dyed scarlet by the yellow dye *Amur cork* (*berberine*), and then over-dyed with the red dye *safflower*. Figure 29 shows a fluorescence fingerprint measured for a standard silk-fabric specimen dyed green by the blue dye *indigo*, and then over-dyed with the yellow dye *Amur cork*. Figures 28 and 29 also show enlarged subregions of the full fluorescent-footprint images. The low measured values of the fluorescence intensity in these cases have the effect of obscuring contour-line peaks in the full fluorescence fingerprint plots, so we present these enlargements with lower intensities to ensure that contour-line peaks are readily discernible.

As is clear from a comparison of these fluorescence fingerprints, the fluorescence fingerprint measured for a given dye has a unique shape that is characteristic of the dye, with a unique contour-line peak defined by the maximum excitation wavelength (λ_{Ex}) and the maximum fluorescence wavelength (λ_{Em}), which are also characteristic of the dye. In fluorescence fingerprints for specimens dyed scarlet or green by over-dyeing using two different dyes, we see two distinct contour-line peaks, one produced by each of the dyes. The dyed silk fabrics measured to yield fluorescence fingerprints for our reference database of known colorants were provided by our research collaborator Yasuko Shimoyama, an expert in the field of dyes, while the colorized *washi* prints measured to yield fluorescence fingerprints for our database were provided by the late Shinya Katsuhara (Painter name: Inuki Tachihara). Note that fluorescence fingerprints (contour-line patterns) measured for fabrics dyed with a given dye should be independent of the type of fabric—which may be silk, hemp, or another fiber—and should be identical to fluorescence fingerprints measured for *washi* specimens colorized by the same dye.




Fig. 26 Measured 3DF spectrum for safflower dye (carthamin) and corresponding contour-plot (or fluorescence fingerprint).




Fig. 27 Fluorescence fingerprints for standard specimens consisting of silk fabrics colored by single dyes or *washi* paper colored by single dyes.




Fig. 28 Fluorescence fingerprints for standard specimens of silk fabrics dyed scarlet by *Amur cork*, and over-dyed with *safflower*.




Fig. 29 Fluorescence fingerprints for standard specimens of silk fabrics dyed green by *indigo*, and over-dyed with *Amur cork*.

Identifying the colorants used to dye fabrics in past eras

Figure 30 shows a brocade, discovered in 1994 in Lhasa, Tibet, that dates from the 16th or 17th century (from the late Muromachi era to the Azuchi-Momoyama era of Japanese history) and is created from fabric made in China. The decorative motifs woven into this brocade include: a dragon, which was thought to symbolize the Emperor of China, woven from gold thread and incorporating five talons; a cloud, woven from blue and green thread; and a flame woven from scarlet thread. Susumu Shirai, then the chief engineer of Tatsumura Textile Co., Ltd., named this object the "brocade with dragon and cloud patterns on light-brown ground". With dimensions of 290 mm (tall) by 280 mm (wide), the brocade is too large to fit into the specimen chamber of our spectrophotometer; instead, we used our optical-fiber probe setup to make measurements of specific subregions. First, Figure 32 shows the fluorescence fingerprint obtained from measurements (at point X) of the blue thread used to weave the cloud motif shown on the left in Figure 31. This fluorescence fingerprint matches that obtained for the *indigo*-dyed silk fabric in Figure 27G. Next, Figure 33 shows the fluorescence fingerprint obtained from measurements (at point Y) of the green fabric used to weave the cloud motif. This fluorescence fingerprint matches that obtained for the silk fabric dyed with both *indigo* and *Amur cork* in Figure

29. Y1 and Y2 indicate contour-map peaks due to Amur cork dye, while Y3 indicates a contour-map peak due to *indigo* dye. Finally, Figure 34 shows the fluorescence fingerprint obtained from measurements (at point Z) of the scarlet thread used to weave the flame motif. This fluorescence fingerprint and the peak contour line match the measured data shown in Figure 28 for silk fabric dyed with both safflower and Amur cork. Contour-line peaks Z1 and Z2 are due to safflower dye, while peaks Z3, Z4, and Z5 are due to Amur cork dye. Based on these findings, we conclude that the blue thread woven into the brocade is colored by *indigo* alone, the green thread is colored by both *indigo* and Amur cork, and the scarlet thread is colored by both safflower and Amur cork³²⁻³⁴.




Fig. 30 A brocade dating from the 16th or 17th century and known as the "brocade with dragon and cloud patterns on light-brown ground". The measurement positions are indicated on the back side of the brocade (right).




Fig. 31 Measurement points for brocade specimen in Figure 30.




Fig. 32 Fluorescence fingerprint obtained from measurements of blue thread in brocade (measurement point X).




Fig. 33 Fluorescence fingerprint obtained from measurements of green thread in brocade (measurement point Y).




Fig. 34 Fluorescence fingerprint obtained from measurements of scarlet thread in brocade (measurement point Z).

Analyzing green colorants used in ukiyo-e prints

By analyzing measured 3DF spectra and XRF spectra, we showed that a mixture of *indigo* dye and orpiment (As_2O_3) pigment was used as a green colorant in a *ukiyo-e* print³⁵. The *ukiyo-e* print in question was published around 1821 (Bunsei 4, in the late Edo period) as part of a three-print series by Gototei Kunisada titled "Eight views of Edo: Evening Snow at Mokuboji". Figure 35 shows the fluorescence fingerprint and XRF spectra obtained for a measurement point lying within a green region of the print.

We first note that the fluorescence fingerprint obtained for this green region resembles the fluorescence fingerprint for *indigo* shown in Figure 27G. However, the maximum-excitation wavelength (λ_{Ex}), indicated by the peak contour line, is shifted slightly in the long-wavelength direction compared to that for *indigo*. The presence of elemental As is also

detected in this green-colored region, from which we conclude that the yellow pigment *orpiment* was used. To proceed, we prepared a standard specimen consisting of a mixture of blue-colored *indigo* and yellow-colored orpiment applied to *washi* paper, then measured 3DF and XRF spectra for this specimen. As expected, the XRF spectrum for the standard specimen reveals the presence of As, and we obtained the fluorescence fingerprint shown in Figure 36. This fluorescence fingerprint matches that obtained for the green region of the *ukiyo-e* specimen in Figure 35, and the positions of the peak contour lines ($\lambda_{Ex}/\lambda_{Em}$) match as well. Based on this experimental confirmation, we conclude that the green colorant used at that time, not only in this *ukiyo-e* print but more generally, was a mixture of *indigo* and orpiment.

As a side note, this example raises two points of general interest. First, as we see from the fluorescence fingerprint for the green colorant produced by blending *indigo* and orpiment, features in fluorescence fingerprints due to dyes may appear distorted for specimens that also contain pigments. Second, the presence of paramagnetic ions such as Fe, Ni, Cr, Cu, or Co can dramatically reduce the fluorescence intensity emitted by fluorescent molecules (due to extinction effects), and thus extra care is required in such cases^{22,36}.




Fig. 35 3DF and XRF spectra measured for *ukiyo-e* print in series "Eight views of Edo: Evening Snow at Mokuboji".




Fig. 36 3DF spectrum measured for standard specimen of green colorant produced by blending orpiment with *indigo*.

Fluorescence fingerprints for *washi* paper used in *ukiyo-e* prints

The production of *ukiyo-e* prints is a multi-stage collaboration between an artist, a sculptor, and a printer. The process begins when the artist gives the sculptor a rough draft of the image to be depicted by the print. The sculptor then prepares several items: a master printing key-block (or *omohan*) with lettering for the title of the image and the name of the artist, a schematic diagram (an outline sketch) of the image based on the artist's rough draft, and a colorization template (*kyougouzuri*) used to communicate color information between artist and sculptor. More specifically, the

kyougouzuri is a version of the schematic diagram on which the artist adds instructions specifying the colors to be used in various subregions of the image. Based on these instructions, the sculptor then sculpts a set of *color blocks*—one for each color in the image—that will be used to print the colorized image, one color at a time. Next, the color blocks and the master key-block are delivered to the printer, who first uses the latter to form an outline image on *washi* paper before proceeding to production with the former. Consequently, the colorants used to produce the various colors in the image are stamped onto the uppermost layer of the *washi*. If *dyes* are stamped onto the *washi* then we can measure their 3DF spectra, while if *pigments* are stamped onto the *washi* then we can measure fluorescence emissions from the base layer of the *washi* (cellulose) itself and use this data as a fluorescence fingerprint.

Figure 37 shows XRF spectra and fluorescence fingerprints measured for a *ukiyo-e* print by Gototei Kunisada, renamed Utagawa Toyokuni the 3rd, entitled "Eight Views of Floating World: Night Rain at Kioroshi River". This *ukiyo-e* print was published in 1855 (Ansei 2), after Katsushika Hokusai first used Prussian blue to establish landscape images as a subgenre of *ukiyo-e*. Thus, by the time this print was published the use of Prussian blue in *ukiyo-e* prints had become standard.

The upper and lower sets of data plots in Figure 37 respectively indicate results for measurement point A—at which no colorants were stamped—and for measurement point B, at which blue colorant was stamped. Comparing the two datasets, we see that the fluorescence fingerprint for measurement point B—which lies in a blue-colored region of the specimen—matches the fluorescence fingerprint for measurement point A, which lies in an uncolored region. This demonstrates that the fluorescence fingerprint for measurement point B is in fact a fluorescence fingerprint for the *washi* (cellulose) itself, and that the blue colorant present at point B is a pigment that does not emit fluorescence at ultraviolet or visible wavelengths. Thus, based on the presence of Fe revealed by the measured XRF spectrum for measurement point B, we conclude that the blue pigment used here is Prussian blue. This example demonstrates that it cannot be immediately concluded that dyes were used simply because the fluorescence fingerprint was obtained from that region; one must also consider the fluorescence fingerprint for the background (empty regions or regions at which no colorant is present), and if the fingerprints match one must assume that a pigment has been stamped into the specimen and proceed to XRF analysis.




Fig. 37 Fluorescence fingerprints and XRF spectra measured for *white (washi)* and *blue* regions of *ukiyo-e* print "Eight Views of Floating World: Night Rain at Kioroshi River".

Learning new functional capabilities from fluorescence fingerprints

The study of bygone eras can yield valuable new insights and perspectives. Here we discuss a novel twist on this idea of *learning from the past*. Through our studies of colorants used for cultural properties in previous eras, we discovered a new type of optical functionality, which was successfully exploited to develop new commercial products of use in the modern world³⁷⁾.

The fluorescence fingerprint for safflower shown in Figure 26 indicates absorption and excitation of 550 nm light with emission of 605 nm light. 550 nm light is seen by humans as green, while 605 nm light is seen as yellowish-red (or

reddish orange). Interpreting this phenomenon as a demonstration of the functional optical capabilities of safflower, we investigated the possibility of exploiting this phenomenon to productive ends in the modern world. In contrast to LED light sources, which are now starting to become ubiquitous, the light sources used most often in previous eras included daylight fluorescent lamps and white fluorescent lamps. As any makeup artist can attest, even the most beautifully accentuated face can appear sallow under a fluorescent lamp. The reason for this, as shown in Figure 38, is that these fluorescent lamps emit blue-tinted and green-tinted light with high intensity, but emit almost no red-tinted light. This is the origin of the disappointing appearance of cosmetics when viewed under fluorescent lamps. But this diagnosis suggests a potential remedy: If one could find a substance with a particular optical functionality—specifically, the ability to absorb and eliminate blue-tinted and green-tinted light while generating red-tinted light—then this substance could be added to cosmetic products to mitigate the undesirable phenomenon described above. To test this idea we dyed cellulose with safflower, ground the cellulose into a fine powder, and added it to foundation and other cosmetic powders. This led to the emergence of new commercial products released in 1994 by the Japanese cosmetic firms Shiseido and Chifure^{38,39}.

On a separate but related note, in Europe, the red colorant cochineal (carminic acid) has been one of the most widely used colorants ever since the birth of lipstick. Cochineal does not fluoresce even when irradiated by ultraviolet light. However, blending cochineal with safflower has the effect of shifting the maximum fluorescence wavelength of safflower from 605 nm to longer wavelengths, and values in the range 625-635 nm ($\lambda_{Em}/\lambda_{Ex}$ ranging from 625/545 to 635/400 nm) were found to produce red-tinted light. Based on this discovery, in 2005 the French cosmetic firm Chanel introduced a lipstick named "AKA" for sale exclusively in Japan (Figure 39).




Fig. 38 Coloration of light emitted from fluorescent light sources. Solid curve: Daylight fluorescent lamp. Dashed curve: White fluorescent lamp.



Fig. 39 Chanel's AKA lipstick.

6. Conclusions

In this article we presented a collection of measurement techniques we have developed for non-destructive analysis of the colorants used in cultural properties. Our approach, which we have deployed and are continuing to deploy in a wide range of practical scenarios, involves collecting and analyzing five mutually complementary sets of experimental data: infrared photographs, microscopy images, visible/near-infrared reflection spectra, X-ray fluorescence spectra, and three-dimensional fluorescence spectra. For each type of data, we described the experimental procedures and analytical instruments used to perform the measurements, discussed the interpretation of the data and the information that may be derived from it, and presented practical case studies involving pigments and dyes used in cultural properties. Colorants may be pigments, dyes, or mixtures of both, and no single experimental approach could possibly suffice to lend full insight into both possibilities. At the very least, one requires a minimum of two types of analysis: X-ray spectroscopy, to investigate the elemental composition of pigments, and fluorescence spectroscopy (three-dimensional fluorescence spectra), to analyze dyes; moreover, the procedures and instruments used for these analyses must be amenable to on-site measurements at museums or other remote sites at which cultural properties are curated. As we demonstrated, the insights delivered by these two types of spectroscopy are richly complemented by other types of data, which—despite being somewhat classical in nature—furnish valuable insights and lend comprehensive understanding: infrared

photographs, microscopy images depicting enlarged views of fine-grained structures, and visible/near-infrared reflection spectra.

Finally, in analyzing colorants in cultural properties our discussion emphasized the importance of combining experimental data with other types of knowledge—including awareness of the physical and chemical properties of colorants, understanding of the historical factors responsible for the introduction of various types of colorants, and an appreciation of the techniques by which cultural properties were created—to foster the broad multidisciplinary mindset required for a proper interpretation of measured data.



An on-site analysis in progress at the Museum of Fine Arts, Boston.

Acknowledgements

Our research has been made possible by the generous assistance of many individuals who have provided opportunities to work closely with precious cultural properties. It is thanks to these individuals that we have been able to perform non-destructive analyses to study pigments, dyes, and other materials used to colorize famous works of art, and we are grateful for their assistance. In particular, we acknowledge the following institutions for offering access to the cultural properties discussed in the text.

- Hiroshima Museum of Art
 - “Daubigny’s Garden”, oil painting (Figures 2,4)
- Naha City Museum of History
 - Japanese National Treasure: Naha-15-, Bingata-07-07, Bingata garment (Figures 5,7,16)
 - Japanese National Treasure: Naha-20-, Bingata-12-12, Bingata garment (Figures 8,14)
 - Japanese National Treasure: Naha-26-, Bingata-18-18, Bingata garment (Figure 20)
 - Japanese National Treasure: Naha-80-, Earthenware-01, Square-base candle stand with green Ryoku-Yu glaze (Figures 17,23)
 - Japanese National Treasure: Naha-77-, Lacquerware-09, Tsuikin-decorated black-lacquered wood tray (Figure 22)
- Hagi Uragami Museum
 - Katsushika Hokusai, "Thirty-Six Views of Mount Fuji" (Front views of Mount Fuji) (Figure 13)
- Asakura Shrine, Fukui Prefecture
 - Images of Rashomon (Figure 21)

References

- 1) Nobuyuki Kamba: Technical note on ultraviolet photography, ultraviolet fluorescent photography, and infrared photography. *Sokei Art Academy Restoration Laboratory Report*, Vol. 5, pp. 56-62 (1987) (in Japanese).
- 2) Susumu Shimoyama: Part 1, Scientific investigation. In: All about van Gogh's "Daubigny's Garden": An illustrated guide. Hiroshima Museum of Art and Kibi International University, pp. 16-26 (2008) (in Japanese).
- 3) Susumu Shimoyama: Discovering the black cat hidden beneath van Gogh's "Daubigny's Garden." *ISOTOPE NEWS*, Apr. No. 660, pp. 11-16 (Japan Radioisotope Association, 2008) (in Japanese).

- 4) Hiroshi Yoshida: Essay 3, "Daubigny's Garden": Partial reproduction and van Gogh's oil-painting technique. In: All about van Gogh's "Daubigny's Garden": An illustrated guide. Hiroshima Museum of Art and Kibi International University, pp. 37-41 (2008) (in Japanese).
- 5) Susumu Shimoyama and Hideo Matsui: Studying the historical trajectory leading to the use of Prussian blue in *ukiyo-e* prints via non-destructive analysis based on fluorescent X-rays from radio isotopes, Reports of the Reimei Research Report Meeting, Japan Atomic Energy Research Institute, pp. 440-454 (2004) (in Japanese).
- 6) Susumu Shimoyama, Hideo Matsui, Yasuko Shimoyama: Non-Destructive Identification of Blue Colorants in Ukiyo-e Prints by Visible-Near Infrared Reflection Spectrum Obtained with a Portable Spectrophotometer Using Fiber Optics, *Bunseki Kagaku*, **55**, pp. 121-126 (2006) (in Japanese).
- 7) Hideo Matsui, Susumu Shimoyama, Yasuko Shimoyama: Studying the historical trajectory leading to the use of Berlin blue and the *Thirty-Six Views of Mount Fuji* as a pioneering work establishing the genre of landscape images in *ukiyo-e* prints via non-destructive methods for identifying blue colorants in *Nishiki-e*. *Hokusai Kenkyu*, No. 37, pp. v-liv/pp. 5-54 (Tokyo Bijutsu, 2005) (in Japanese).
- 8) Bokashi, in *Ukiyo-e Daijiten* (Ukiyo-e Encyclopedia), edited by the International Ukiyo-e Society p. 445 (Tokyodo Shuppan, 2008) (in Japanese).
- 9) Omohan, in *Ukiyo-e Daijiten* (Ukiyo-e Encyclopedia), edited by the International Ukiyo-e Society p. 110 (Tokyodo Shuppan, 2008) (in Japanese).
- 10) Susumu Shimoyama and Yasuko Noda: Simple portable X-ray fluorescence analyzer using a low-level radioactive isotope (^{241}Am) as an X-ray source: Applications to non-destructive analysis of inorganic colorants in traditional Japanese *ema* artworks. *Bunseki Kagaku*, **49**, pp. 1015-1021 (2000) (in Japanese).
- 11) Izumi Nakai, editor, the Japan Society for Analytical Chemistry X-Ray Analysis Discussion Group, supervisor: X-ray fluorescence analysis in practice, 2nd edition, pp. 257-260, Appendix D: Table of characteristic X-ray energies and absorption edges, etc. (Asakura Shoten, 2016) (in Japanese).
- 12) Izumi Nakai, editor, the Japan Society for Analytical Chemistry X-Ray Analysis Discussion Group, supervisor: X-ray fluorescence analysis in practice, 2nd edition, Chapter 1: Foundations of fluorescent X-rays, Section 1.5.6, Fluorescent X-ray Spectra, pp. 7-9 (Asakura Shoten, 2016) (in Japanese).
- 13) R. J. Gettens and G. L. Stout, *Painting Materials: A Short Encyclopedia*, Revised Edition, translated into Japanese by Tsuneyuki Morita, pp. 153-156 (Bijutsu Shuppan-Sha, 1999).
- 14) B. Muhlethaler, J. Thissen: 'Smalt', Artists' Pigments, A Handbook of Their History and Characteristics, Vol. 2, pp. 113-130 (© 1993 National Gallery of Art).
- 15) Susumu Shimoyama and Yasuko Shimoyama: Japanese National Treasures: Materials related to the Ryukyu Kingdom Sho Dynasty, Report on Non-destructive analytical investigations of craft works (including dyed fabrics), "Lacquerware-09: Tsuikin-decorated black-lacquered wood tray", *Naha City Museum of History Bulletin*, No. 1, pp. 31 (Naha City Museum of History, 2009) (in Japanese).
- 16) Kazumi Murose, *Lacquer Culture: Inheriting a Tradition of Japanese Beauty*, <Creating Lacquer Colors>, Kadokawa Selections 343, pp. 33-37 (Kadokawa Shoten, 2002) (in Japanese).
- 17) Noboru Terada, Kadoaki Oda, Yasushi Oyabu, Tetsu Asami, editors: *Lacquer: Science and practical techniques*, History of Colored Lacquer Pigments, pp. 200-211 (Riko Shuppansha, 2002) (in Japanese).
- 18) Holbein Industrial Technology Division, editors: *Handbook of Painting Materials: Indigo-shaded Blue*, pp. 59-60 (Chuokouronbijyutsu, 1991) (in Japanese).
- 19) Susumu Shimoyama, Koji Oshita, Yasuko Shimoyama: 2010/2011 Report on non-destructive investigation of colorants in dyed fabric artifacts related to Uzagaku costumes, Gold-inlaid green-lacquered wood tray, *Shurijo Koen Kanri Center Annual Report on Investigative Research and Programs to Promote Outreach and Awareness*, No. 3, p. 84 (Okinawa Churashima Foundation, 2013) (in Japanese).
- 20) Susumu Shimoyama, Koji Oshita, Yasuko Shimoyama: 2015 Report on investigation of colorants in Sho-dynasty artifacts, Square-base candle stand with green Ryoku-Yu glaze, *Shurijo Koen Kanri Center Annual Report on Investigative Research and Programs to Promote Outreach and Awareness*, No. 7, pp. 66-68 (Okinawa Churashima Foundation, 2015) (in Japanese).
- 21) Masataro Onishi, *Pottery Glazes: Theory and Practical Recipes* (Revised Edition), pp. 180-183 (Rikogakusha, 2007) (in Japanese).
- 22) Susumu Shimoyama and Yasuko Noda: Extinction effects of iron (II) and copper (II) ions on measured fluorescent intensities in three-dimensional fluorescence spectra of dyed fabric specimens, *Bunseki Kagaku*, **47**, pp. 295-301 (1998) (in Japanese).

- 23) Juuji Yoshimura, editor: Kagaku Jiten (Chemical Dictionary), Second Edition, Susumu Shimoyama, author: Natural dyes: Primary natural dyes used in traditional Japanese crafts, pp. 936-937 (Morikita Shuppan, 2009) (in Japanese).
- 24) Susumu Shimoyama and Yasuko Noda: Non-destructive identification of dyes used in ancient dyed-fabric artifacts from three-dimensional fluorescence spectra, *Bunseki Kagaku*, **41**, pp. 243-249 (1992) (in Japanese).
- 25) Susumu Shimoyama and Yasuko Noda: NON-DESTRUCTIVE ANALYSIS OF UKIYO-E PRINTS Determination of Plant Dyestuffs used for Traditional Japanese Woodblock Prints, Employing a Three-Dimensional Fluorescence Spectrum Technique and Quartz Fiber Optics, *Dyes in History and Archaeology*, **15**, pp. 27-42 (1997).
- 26) Susumu Shimoyama, Yasuko Noda, Shinya Katsuhara: Using optical fibers to measure three-dimensional fluorescence spectra for non-destructive identification of colorants used in ancient Japanese *ukiyo-e* prints, *Bunseki Kagaku*, **47**, pp. 93-100 (1998) (in Japanese). Winner of the 1998 Analytical Chemistry Research Article award from the Japan Society for Analytical Chemistry.
- 27) Susumu Shimoyama: Spaulding Collection, Museum of Fine Arts, Boston, Collaborative investigation of colorants in *ukiyo-e* woodblock prints, Colorants used in works by Torii Kiyonaga (Initial Report), *Research on Information Science for Cultural Properties*, **5**, pp. 43-53 (Kibi International University Research Advancement Center for Cultural Property, 2008) (in Japanese).
- 28) Juuji Yoshimura, editor: Kagaku Jiten (Chemical Dictionary), Second Edition, Susumu Shimoyama, author: Three-dimensional fluorescence spectra, pp. 553-554 (Morikita Shuppan, 2009) (in Japanese).
- 29) Seiji Hasegawa, author, Japan Forensic Science Appraisal Center, supervisor: Scientific Forensics, pp. 24-33 (Natsume-sha, 2004) (in Japanese).
- 30) Analysis Laboratory of Forensic Science, supervisor: Scientific Forensics, pp. 31-46 (Shufunotomoshia, 2010) (in Japanese).
- 31) Katsumi Kurosaka, editor: Newly revised and expanded supplement, Kokushi Taikei (Annals of Japanese History), Popular Edition, Engishiki, Volume 2, Engishiki, Scroll 14, Nuidonoryou, Zassenyoudo, pp. 400-405 (Yoshikawa Koubunkan, 1990) (in Japanese).
- 32) Susumu Shimoyama and Yasuko Noda: Using optical fibers to measure three-dimensional fluorescence spectra to identify dyes used to colorize fabrics for brocades in ancient China, *Bunseki Kagaku*, **46**, pp. 571-578(1997) (in Japanese).
- 33) Susumu Shimoyama, Yasuko Noda: Using optical fibers to measure three-dimensional fluorescence spectra to identify dyes used for ancient green dyed fabrics, *Bunseki Kagaku*, **46**, pp. 791-799 (1997) (in Japanese).
- 34) Susumu Shimoyama and Yasuko Noda: NON-DESTRUCTIVE ANALYSIS OF DYES IN A CHINESE BROCADE: Determination of Plant Dyes in a 16th/17th-Century Textile by a Three-Dimensional Fluorescence Spectrum Technique with Fibre Optics, *Dyes in History and Archaeology*, **15**, pp. 70-84 (1997).
- 35) Yasuko Noda and Susumu Shimoyama: Non-Destructive Analysis Ukiyo-e, Traditional Japanese Woodblock Prints, Using a Portable X-ray Fluorescence Spectrometer, *Dyes in History and Archaeology*, **18**, pp. 73-86 (1999).
- 36) Yasuharu Nishikawa and Keizou Hiraki: Practical instrumental analysis series: Fluorescence and phosphorescence methods: Extinction due to paramagnetic ions, pp. 37-38 (Kyoritsu Shuppan, 1987) (in Japanese).
- 37) Susumu Shimoyama: Natural colorants: A new trend in coloring agents for developing makeup cosmetics. *COSMETIC STAGE*, Vol. 2, No. 4, pp. 1-7 (Technical Information Institute, 2008) (in Japanese).
- 38) Kazuhisa Ono, Naoko Watanabe, Shigenori Kumagai, Susumu Shimoyama, Yasuko Noda: Cosmetic ingredient, Open patent, Public application number 8-059427 (Shiseido, Den Material, 1996) (in Japanese).
- 39) Naoki Moritaka, Masao Oinuma, Susumu Shimoyama: Natural colorant-processed powder and cosmetic ingredient based thereon, Open patent, Public application number 2000-044828 (Den Material, Chifure, 2000) (in Japanese).