


Phosphorescence Excitation Spectrum of Benzophenone at Liq.N₂ Temperature

INTRODUCTION

By using the low temperature measurement accessory, fluorescence/phosphorometric analysis at the liquid nitrogen temperature (-196°C) is possible. The fine structure of a sample, that is not shown at a room temperature, can be analyzed. F-7000 can measure the phosphorescence spectrum and phosphorescence lifetime and the information on the triplet excited state of the sample can be obtained. Even by the high resolution analysis with a narrow slit width, the peaks can be measured depending on the difference in the vibrational states of benzophenone with a sufficient S/N.

SAMPLE		ACCESSARY
SAMPLE NAME	: Benzophenone : (C ₆ H ₅) ₂ CO mol.wt. 182.22 (Wako, Osaka, Japan)	Low temperature measurement accessory (P/N : 5J0-0112)
SOLVENT	: Ethanol	
CONCENTRATION	: 1.0 × 10 ⁻³ mol/L	
INSTRUMENT CONDITIONS		PEAKS (nm)
INSTRUMENT	: F-7000	1 : 327
DATA MODE	: Phosphorescence	2 : 339
CHOPPING FREQUENCY	: 40 Hz	3 : 354
EM WAVELENGTH	: 444 nm	4 : 369
EX BANDPASS	: 2.5 nm	
EM BANDPASS	: 20 nm	
	SCAN SPEED : 60 nm/min RESPONSE : 0.5 s PHOTOMULTIPLIER : R3788 PHOTOMULTIPLIER VOL.: 250 V	



KEY WORDS

Material·Processing Material Related, Industrial Chemical, Benzophenone, Diphenyl Ketone, Organic Solvent, Low Temperature, Liquid Nitrogen, 77K, Phosphorescence Excitation Spectrum, Liq.N₂ Temperature, LN₂, FL, F-7000

Fluorophotometer (FL)


Sheet No. FL080018-01

Fluorescence Spectrum of Benzophenone at Liq.N₂ Temperature

INTRODUCTION

By using the low temperature measurement accessory, fluorescence/phosphorometric analysis at the liquid nitrogen temperature (-196°C) is possible. The fine structure of a sample, that is not shown at a room temperature, can be analyzed. F-7000 can measure the phosphorescence spectrum and phosphorescence lifetime and the information on the triplet excited state of the sample can be obtained. Even by the high resolution analysis with a narrow slit width, the peaks can be measured depending on the difference in the vibrational states of benzophenone with a sufficient S/N.

SAMPLE	ACCESSARY
SAMPLE NAME : Benzophenone : (C ₆ H ₅) ₂ CO mol.wt. 182.22 (Wako, Osaka, Japan) SOLVENT : Ethanol CONCENTRATION : 1.0 × 10 ⁻³ mol/L	Low temperature measurement accessory (P/N : 5J0-0112)
INSTRUMENT CONDITIONS	PEAKS (nm)
INSTRUMENT : F-7000 EX WAVELENGTH : 282 nm RESPONSE : 0.5 s EX BANDPASS : 10 nm PHOTOMULTIPLIER : R3788 EM BANDPASS : 2.5 nm PHOTOMULTIPLIER VOL. : 250 V SCAN SPEED : 60 nm/min	1 : 414 2 : 443 3 : 477 4 : 517



[Corrected spectrum]

KEY WORDS

Material·Processing Material Related, Industrial Chemical, Benzophenone, Diphenyl Ketone, Organic Solvent, Low Temperature, Liquid Nitrogen, 77K, Fluorescence Spectrum, Liq.N₂ Temperature, LN₂, FL, F-7000

Fluorophotometer (FL)


Sheet No. FL080018-02

Phosphorescence Spectrum of Benzophenone at Liq.N₂ Temperature

INTRODUCTION

By using the low temperature measurement accessory, fluorescence/phosphorometric analysis at the liquid nitrogen temperature (-196°C) is possible. The fine structure of a sample, that is not shown at a room temperature, can be analyzed. F-7000 can measure the phosphorescence spectrum and phosphorescence lifetime and the information on the triplet excited state of the sample can be obtained. Even by the high resolution analysis with a narrow slit width, the peaks can be measured depending on the difference in the vibrational states of benzophenone with a sufficient S/N.

SAMPLE		ACCESSARY
SAMPLE NAME	: Benzophenone : (C ₆ H ₅) ₂ CO mol.wt. 182.22 (Wako, Osaka, Japan)	Low temperature measurement accessory (P/N : 5J0-0112)
SOLVENT	: Ethanol	
CONCENTRATION	: 1.0 × 10 ⁻³ mol/L	
INSTRUMENT CONDITIONS		PEAKS (nm)
INSTRUMENT	: F-7000	1 : 414
DATA MODE	: Phosphorescence	2 : 443
CHOPPING FREQUENCY	: 40 Hz	3 : 477
EX WAVELENGTH	: 282 nm	4 : 517
EX BANDPASS	: 20 nm	5 : 562
EM BANDPASS	: 2.5 nm	



[Corrected spectrum]

KEY WORDS

Material·Processing Material Related, Industrial Chemical, Benzophenone, Diphenyl Ketone, Organic Solvent, Low Temperature, Liquid Nitrogen, 77K, Phosphorescence Spectrum, Liq.N₂ Temperature, LN₂, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL080018-03

3D Phosphorescence Spectrum of Benzophenone at Liq.N₂ Temperature

INTRODUCTION

By using the low temperature measurement accessory, fluorescence/phosphorometric analysis at the liquid nitrogen temperature (-196°C) is possible. The fine structure of a sample, that is not shown at a room temperature, can be analyzed. F-7000 can measure the phosphorescence spectrum and phosphorescence lifetime and the information on the triplet excited state of the sample can be obtained. Even by the high resolution analysis with a narrow slit width, the peaks can be measured depending on the difference in the vibrational states of benzophenone with a sufficient S/N.

SAMPLE

SAMPLE NAME : Benzophenone
 : (C₆H₅)₂CO mol.wt. 182.22
 (Wako, Osaka, Japan)

SOLVENT : Ethanol


CONCENTRATION : 1.0 × 10⁻³ mol/L

ACCESARRY

Low temperature measurement accessory
 (P/N : 5J0-0112)

INSTRUMENT CONDITIONS

INSTRUMENT	: F-7000				
DATA MODE	: Phosphorescence	SCAN SPEED	: 240 nm/min	FULLSCALE	: 300
CHOPPING FREQUENCY	: 40 Hz	RESPONSE	: Auto	DIVISION NUMBER	: 5
EX BANDPASS	: 5 nm	PHOTOMULTIPLIER	: R3788		
EM BANDPASS	: 5 nm	PHOTOMULTIPLIER VOL.:	: 400 V		



[Corrected spectrum]

KEY WORDS

Material·Processing Material Related, Industrial Chemical, Benzophenone, Diphenyl Ketone, Organic Solvent, Low Temperature, Liquid Nitrogen, 77K, Three-dimensional Phosphorescence Spectrum, 3D, Liq.N₂ Temperature, LN₂, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL080018-04

Phosphorescence Lifetime of Benzophenone

INTRODUCTION

By using the low temperature measurement accessory, fluorescence/phosphorometric analysis at the liquid nitrogen temperature (-196°C) is possible. The fine structure of a sample, that is not shown at a room temperature, can be analyzed. F-7000 can measure the phosphorescence spectrum and phosphorescence lifetime and the information on the triplet excited state of the sample can be obtained. Even by the high resolution analysis with a narrow slit width, the peaks can be measured depending on the difference in the vibrational states of benzophenone with a sufficient S/N.

SAMPLE

SAMPLE NAME : Benzophenone
 : (C₆H₅)₂CO mol.wt. 182.22
 (Wako, Osaka, Japan)

SOLVENT : Ethanol


CONCENTRATION : 1.0 × 10⁻³ mol/L

INSTRUMENT CONDITIONS

INSTRUMENT	: F-7000		
DATA MODE	: Phosphorescence Lifetime	EX BANDPASS	: 10 nm
CHOPPING FREQUENCY	: 40 Hz	EM BANDPASS	: 10 nm
EX WAVELENGTH	: 282 nm	RESPONSE	: 0.002 s
EM WAVELENGTH	: 444 nm	PHOTOMULTIPLIER	: R3788
		PHOTOMULTIPLIER VOL.	: 250 V

ACCESARRY

Low temperature measurement accessory
 (P/N : 5J0-0112)



[τ = 5.539 ms]

KEY WORDS

Material·Processing Material Related, Industrial Chemical, Benzophenone, Diphenyl Ketone, Organic Solvent, Phosphorescence Lifetime, FL, F-7000

Fluorophotometer (FL)

Sheet No. FL080018-05