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- 1) Y.Watanabe and K.Imai : Pre-Column Labelling for High-Performance Liquid Chromatography of Amino Acids with 7- Fluoro-4-nitrobenzo-2-oxa-1,3-diazole and Its Application to Protein Hydrolysates. J.Chromatogr., 239, 723-732(1982)
- 2) C.Aoyama, T.Santa, M.Tsunoda, T.Fukushima, C.Kitada and K.Imai : A fully automated amino acid analyzer using NBD-F as a fluorescent derivatization reagent. Biomed.Chromatogr., 18, 630-636 (2004)
- 3) M.Tomiya, T.Fukushima, J.Kawai, C.Aoyama, S.Mitsuhashi, T.Santa, K.Imai, T.Toyo ' oka : Alternation of plasma and cerebrospinal fluid glutamate levels in rats treated with N-methyl-D-aspartate receptor antagonist, ketamine. Biomed.Chromatogr., 20, 628-633(2006)
- 4) P.Prados, H.Matsunaga, T.Mori, T.Santa, T.Fukushima, H.Homma, C.Kasai, K.Imai : Changes of plasma L-arginine levels in spontaneously hypertensive rats under induced hypotension. Biomed.Chromatogr., 13, 27-32(1999)
- 5) Y.Watanabe and K.Imai : Sensitive Detection of Amino acids in Human Serum and Dried Blood Disc of 3mm Diameter for Diagnosis of Inborn Errors of Metabolism. J.Chromatogr., 309, 279-286 (1984)
- 6) K.Hashimoto, T.Fukushima, E.Shimizu, N.Komatsu, H.Watanabe, N.Shinoda, M.Nakazato, C.Kumakiri, S.Okada, H.Hasegawa K.Imai and M.Iyo : Decreased Serum Levels of D-Serine in Patients With Schizophrenia. Arch.Gen.Psychiatry, 60, 572-576(2003)
- 7) C.Toriumi and K.Imai : An Identification Method for Altered Proteins in Tissues Utilizing Fluorescence Derivatization, Liquid Chromatography, Tandem Mass Spectrometry, and a Database-Searching Algorithm. Anal.Chem., 75, 3725-3730(2003)
- 8) M.Masuda, C.Toriumi, T.Santa and K.Imai : Fluorogenic Derivatization Reagents Suitable for Isolation and Identification of Cysteine-Containing Proteins Utilizing High-Performance Liquid Chromatography-Tandem Mass Spectrometry. Anal. Chem., 76, 728-735(2004)
- 9) M.Masuda, H.Saimaru, N.Takamura and K.Imai :

An improved method for proteomics studies in C. elegans by fluorogenic derivatization, HPLC isolation, enzymatic digestion and liquid chromatography-tandem mass spectrometric identification. Biomed.Chromatogr., 19, 556-560(2005)

- 10) T.Toyo'oka, S.Uchiyama, Y.Saito and K.Imai : Simultaneous Determination of Thiols and Disulfides by High-Performance Liquid Chromatography with Fluorescence Detection. Anal.Chim.Acta, 205, 29-41(1988)
- 11) T.Toyo'oka, F.Furukawa, T.Suzuki, Y.Saito, M.Takahashi, Y.Hayashi, S.Uzu and K.Imai : Determination of Thiols and Disulfides in Normal Rat Tissues and Hamster Pancreas Treated with N-Nitrosobis(2-oxopropyl)amine using 4-(Amino-sulfonyl)-7-fluoro-2,1,3-benzoxadiazole and Ammonium 7-Fluoro-2,1,3-benzoxadiazole-4-sulfonate. Biomed.Chromatogr., 3, 166-172(1989)
- 12) K.Moriya, K.Nakagawa, T.Santa, Y.Shintani, H.Fujie, H.Miyoshi, T.Tsutsumi, T.Miyazawa, K.Ishibashi, T.Horie, K.Imai, T.Todoroki, S.Kimura and K.Koike : Oxidative Stress in the Absence of Inflammation in a Mouse Model for Hepatitis C Virus-associated Hepatocarcinogenesis. Cancer Res., 61, 4365-4370(2001)
- 13) K.Imai, T.Fukushima, K.Hagiwara and T.Santa : Occurrence of D-Aspartic Acid in Rat Brain Pineal Gland. Biomed.Chromatogr., 9, 106-109(1995)
- 14) K.Sakai, H.Homma, J.Lee, T.Fukushima, T.Santa, K.Tashiro, T.Iwatsubo and K.Imai : Emergence of D-Aspartic Acid in the Differentiating Neurons of the Rat Central Nervous System. Brain Res., 808, 65-71(1998)
- 15) K.Sakai, H.Homma, J.A.Lee, T.Fukushima, T.Santa, K.Tashiro, T.Iwatsubo and K.Imai : Localization of D-Aspartic Acid in Elongate Spermatids in Rat Testis. Arch.Biochem.Biophys., 351, 96-105, 1998.
- 16) K.Hashimoto, T.Fukushima, E.Shimizu, N.Komatsu, H.Watanabe, N.Shinoda, M.Nakazato, C.Kumakiri, S.Okada, H.Hasegawa K.Imai and M.Iyo : Decreased Serum Levels of D-Serine in Patients With Schizophrenia. Arch.Gen.Psychiatry, 60, 572-576(2003)
- 17) H.Hasegawa, T.Fukushima, J-A Lee, K.Tsukamoto, K.Moriya, Y.Ono and K.Imai : Determination of serum D-lactic and L-lactic acids in normal subjects and Diabetic patients by column-switching HPLC with pre-column fluorescence derivatization. Anal. Bioanal. Chem., 377, 886-891(2003)

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Hair mercury analysis with cold vapor - gold amalgamation - atomic absorption spechometry

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【Summary】The challenges of comparing human serum proteins between the patient and healthy individuals lead to discover novel disease markers. However, to identify these markers require a platform, which is automated and efficient for the comprehensive analysis of human serum. In this study, we established fully automated system of on-line

2DnanoLC/LIT-TOFMS, which the data processing protocol called Information Based Acquisition(IBA) was introduced to support the analysis of complex protein mixtures. Here, we evaluated the performance of this system, whether it achieves to increase the identification number of unique proteins automatically from human serum.

Y. Sasakura(日立ハイテクノロジーズ), 他 : Identification of the Interacting Protein of Heme-regulated Phosphodiesterase from Escherichia coli Using Protein Microarray

【Summary】In our previous study, we have developed the novel protein immobilization method for the protein-protein interaction

analyses on the solid surface. Anti-Tag antibody(Ab) was in advance, immobilized on the solid surface and tag-fused target proteins

were captured by antigen-antibody interaction. Using this method, the already known interaction of the heme-regulated phosphodiesterase from Escherichia coli(Ec DOS) was successfully detected on the solid surface. In this study, we adopted this novel immobilization method to identify the protein interacting with Ec DOS. Ec DOS was immobilized on the solid surface via anti-Tag Ab, and E. coli lysate was applied in it. Interacting pro-

teins were captured on the surface, eluted from it, and analyzed by Nano-LC/LIT-TOF MS. We succeeded to identify GntR, which is the transcriptional regulator of Gnt operon, as the novel protein interacting with Ec DOS, and demonstrated their physiological relationships in vivo. These results suggested that this method would be useful to clarify the unknown protein-protein interactions.

K. Deguchi(北海道大学), 他 : Structural Analysis of N-glycopeptides by Using NanoHPLC/ESI-Linear IT-TOF MS and MSn Spectral Matching

4. 2006 China-Japan-Korea Symposium on Environmental Analytical Chemistry (2006.9/14 ~ 18 中国)

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8. Human Proteome Organisation, 5th Annual World Congress(2006.10/28~11/1 USA)

M. Watanabe(日立ハイテクノロジーズ), 他: The Identification of Serum Proteins using Information Based Acquisition Techniques on Nano-LC/LIT-TOFMS

【Summary】 This study reports the evaluation of Information Based Acquisition(IBA) method on NanoFrontier LD for the identification of serum proteins.

1. The repetitive measurement of the human

serum by NanoFrontier LD with IBA increased the number of identified proteins.

2. A model biomarker added to the human serum at the patient level was identified by the repetitive measurement with IBA.

Y. Sasakura(日立ハイテクノロジーズ), 他: Sample preparation methods using immobilized proteins prior to the Mass spectrometry analyses

【Summary】 The concentration of the target proteins and the rapid proteolytic digestion would become necessary for the focused proteome analyses using Mass spectrometry (MS) In this study, we adapted the protein immobilization techniques to develop these sample preparation methods.

First, a method to capture the interacting protein was constructed. Target protein was immobilized on the solid surface via anti-Tag Ab, and clued lysate was applied in it. Interacting proteins were captured on the surface, eluted from it, and analyzed by Nano-LC/LIT-TOF MS. We succeeded to identify the novel

protein interacting with model protein, and demonstrated their physiological relationships in vivo. These results suggested that this method would be useful to clarify the unknown protein-protein interactions.

Second, protein digestion system using immobilized enzymes was developed. We constructed the reaction unit for the micro-scale sample convection on the enzyme-immobilized solid surface. Proteins were digested using this unit, and successfully identified by MS. Glycopeptidase F(NGF) an enzyme that cleaves N-glycans from glycoproteins, was also immobilized and used to remove the gly-

cochains from human immunoglobulin G (IgG) Trypsin and NGF were immobilized on the single solid surface and used to obtain peptides and glycochains from IgG, simultaneously. Glycochains were labeled with fluorescence reagent and analyzed by HPLC. Several

peaks corresponding to the glycochains of IgG were detected. These results suggested that the simultaneous digestion system, by immobilized multiple enzymes(trypsin and NGF)would be effective for the rapid structural analysis of proteins and glycochains.

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
Fig. 1 Total ion chromatogram and extracted ion chromatograms of the glycerophospholipids standard mixture.

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11. THE 2006 MEETING OF THE SOCIETY FOR GLYCOBIOLOGY(2006.11/15~19 USA)

H. Ito(日立ハイテクノロジーズ), 他: Structural Analysis of O-Glycopeptides Employing Negative- and Positive-Ion MSⁿ Spectra Acquired by nanoHPLC/ESI-Linear Ion Trap Time-of-Flight Mass Spectrometer

【Summary】Recently, we reported that MSⁿ spectral matching simply based on correlation coefficient calculations between positive/negative-ion MSⁿ spectra is useful and applicable to the structural assignment of PA N-glycans and N-glycans binding to peptides without releasing them, enzymatically or chemically. In this presentation, this approach is extended to the structural assignment of Oglycan of glycopeptides. This direct assignment method of O-glycans is particularly suitable for O-glycopeptides which is still missing a versatile enzyme like PNGase F for N-glycopeptides. Experiments were performed by using a Hitachi NanoFrontierL system consisting of a capillary HPLC and an ESI-Linear IT-TOF MS. The samples used were synthetic O-glycopeptides binding a sialyl N-acetylglucosamine(3'SLN) and a sialyl Lewis x(sLex). The samples were dissolved in water/acetoni-

trile. Flow rate was 200 nL/min. MSⁿ spectra were acquired varying CID gain(energy) in 0.6-2.6. Positive-ion MS² spectra caused neutral losses of a fucose and a sialic acid, and therefore, they were useless for the structural assignment of O-glycans. In the negative-ion mode, the CID MS² spectra of O-glycopeptides showed a relatively abundant glycoside-bond cleavage between the core Nacetylglucosamine(GlcNAc) and serine that yields deprotonated C3-type fragment ions of O-glycan. The structure of sLex oligosaccharide was simply assigned by comparing the CID MS³ spectrum derived from the C3-type fragment ion with the CID MS² spectra of the sLex and sLea standards(i.e., negative-ion MSⁿ spectral matching). The amino acid sequence of the peptide including the glycosylation site was determined from the MS² spectrum in the positive-ion mode.

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