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Analysis of capsaicinoids in foods using LaChromUltra System

Jaeho HA* Hye-Young SEO You-Shin SHIM Hyun-Jin NAM

Introduction

Capsaicin has been recognized for many years as the pungent principle of capsicum. In a few studies, it has become clear that other structurally similar components, capsaicinoids, also contribute to the pungency of red pepper⁽¹⁻³⁾. Among them, capsaicin (CAP) is the most abundant heat principle of capsicum and its heat pungency is almost twice as much as that of homocapsaicin (HCAP), homodihydrocapsaicin (HDHC), or nordihydrocapsaicin (NDHC)⁽⁴⁻⁵⁾.

Because the chemical properties and solubility of capsaicin are virtually the same as those of the accompanying capsaicinoids, the natural capsaicin was considered to be an individual substance for a long time. The first separation of the natural substance into capsaicin and dihydrocapsaicin was carried out by paper chromatography⁽⁶⁾. To date, about

20 natural capsaicinoids have been identified in the fruits of red pepper, although even this number can hardly be considered as final (**Table 1**).

The relative contents of the components depend on the species and sort of red pepper may vary within rather wide limits. For example, in the Japanese sort of pepper known as “Santaka” the ratio of capsaicin to dihydrocapsaicin is 70:30, where in the Cayenne species, this ratio is 47:53.

The obtaining of natural capsaicin on a commercial scale is usually based on the infusion of crushed dry pepper fruits with 80-95% ethyl alcohol. Some of the large scale production can be achieved using

Jaeho HA

Table 1. Characteristics of capsaicinoids

Capsaicinoids	Abbrev.	MF	MW	Scoville heat units	Structure
Capsaicin	CAP	C ₁₈ H ₂₇ NO ₃	305.4	16,000,000	
Dihydrocapsaicin	DHC	C ₁₈ H ₂₉ NO ₃	307.4	15,000,000	
Nordihydrocapsaicin	NDHC	C ₁₇ H ₂₇ NO ₃	293.2	9,100,000	
Homodihydrocapsaicin	HDHC	C ₁₉ H ₃₁ NO ₃	321.5	8,600,000	
Homocapsaicin	HC	C ₁₉ H ₂₉ NO ₃	319.4	8,600,000	
Nonivamide	PAVA	C ₁₇ H ₂₇ NO ₃	293.4		

(Reference : Wikipedia, 2009)

* Food Analysis Center, Korea Food Research Institute 516, Baekhyeon, Bundang, Seongnam, Gyeonggi, 463-746, Republic of Korea

diethyl ether, acetone, chloroform and the various concentration of alcohol from 50 to 85%.

Capsaicinoids are characterized by a high biological activity and their pharmacological, neurological and dietetic effectiveness is well known⁽⁷⁾. It has been reported that capsaicin is a highly selective agonist for the transient receptor potential vanilloid receptor 1 (TRPV1; also known as VR1). Toxicity studies with pure *trans*-capsaicin was performed by a few scientist, indicating that pure *trans*-capsaicin does not alter the duration of cardiac action potentials, and causes only very minimal organ toxicities⁽⁸⁾. Other researchers have conducted the experiments of the effect of capsaicin on the herpes simplex virus (HSV), showing that when the onset of treatment was delayed until after intravaginal virus challenge, primary genital skin disease severity was significantly reduced⁽⁹⁾. In animal models, capsaicin was shown to be rather effective against HSV. However, it was noxious when applied to the mucous membranes and caused an unpleasant burning sensation⁽¹⁰⁾.

Kang *et al.*⁽¹¹⁾ have made a report that capsaicin can suppress obesity-induced inflammation by modulating adipokine release from and macrophage behavior in obese mice adipose tissues. Moreover, capsaicin suppressed not only macrophage migration induced by the adipose tissue-conditioned medium, but also macrophage activation to release proinflammatory mediators. Capsaicin may be a useful phytochemical for attenuating obesity-induced inflammation and obesity-related complications. In the view point of analytical method for capsaicinoids in foods, many trials have been made to extract and to analyze the capsaicinoids from food matrix, even though relatively multifarious, the quite successful

results have been achieved to date. The analytical method of capsaicinoids in foods has been developed in many ways, including spectrophotometric and chromatographic technologies such as gas chromatography, liquid chromatography coupled with mass spectrometry (**Table 2**).

Analysis of capsaicinoids in such foods as chilli, paprika, chilli oil is considered to be relatively simple due to simple method of sample preparation, including solvent extraction, filtration by syringe-type membrane filter, and analysis by chromatographic techniques. In the meantime, some of the foods such as *Gochujang* (one of the traditional foods in Korea), *Kimchi* and chilli sauce have to follow the multi steps, including extraction, purification, filtration and analysis.

Recently analytical method of capsaicinoids in *Gochujang* by high performance liquid chromatography (HPLC) has been endorsed as type II and by gas chromatography as type IV in the 30th section of the CODEX Alimentarius Committee on Methods of Analysis and Sampling, Balatonalmadi, Hungary, 9-13, March 2009. Since then the needs for capsaicinoids analysis has been rapidly increased in Korea as well as in many Asian countries. Moreover, Korea government tries to establish the criteria of hot taste label based on concentration of capsaicinoids in foods, instead of scoville heat units. This kind of implementation has been bringing out the necessity of analysis of capsaicinoids in foods to be more effective and rapid than conventional method.

Meanwhile, manufacturing *Gochujang* requires many pretreatment processes to remove such impurities as interfere with analysis before the content of capsaicin is determined, since the raw materials of

Gochujang include soybean, wheat flour, glucose, and sugar. It has been reported that the concentration of capsaicinoids in *Gochujang* varieties ranges from 10 mg/kg to 100 mg/kg; varieties of mild taste contain from 10 mg/kg to 30 mg/kg and strong taste are characterized by a content higher than 70 mg/kg, reaching about to 100 mg/kg⁽¹²⁻¹³⁾. Therefore, the effective analytical method of capsaicinoids in foods such as *Gochujang*, *Kimchi*, and snacks using chromatographic technology has been a challenge due to their complex matrix and interferences in the raw materials.

Analytical method of capsaicinoids has been studied and developed for the chilli and chilli pepper since these have frequently been used as a food ingredient and raw material. Most of difficulties in the analysis are related with the method of sample preparation to remove the interferences, separation of individual capsaicinoids and analysis time.

Ha *et al.*⁽¹⁴⁻¹⁵⁾ have performed the gas chromatographic analysis of capsaicin in *Gochujang* and foods. For the gas chromatographic analysis of capsaicinoids in foods, the sample preparation not only follow the multi-step to remove the interferences in foods but also requires a large consumption of solvent during sample preparation, resulting in environmentally unfriendly, time consuming, and high cost (**Figure 1**).

With some advantages of high performance liquid chromatography (HPLC) over chemical analysis, its application to the analysis of capsaicinoids has

frequently been proposed⁽¹⁵⁻¹⁷⁾. Even though sample preparation for HPLC analysis may be simple and convenient, the interfering substances extracted with capsaicinoids hinder the chromatogram and can easily be absorbed into the column-packing materials, increasing pressure in the system and shortening the column life.

Some of the researches have performed to detect the extremely low concentration of capsaicinoids in self-defense weapon using liquid chromatography-mass spectrometry and -tandem mass spectrometry⁽²⁾.

Recently ultra high performance liquid chromatography (u-HPLC) method has been adopted in many areas of food analysis due to its rapid analysis and remarkably excellent separation. U-HPLC method has been known to be economical and environmentally friendly due to extremely rapid analysis. In concomitant with the fast analysis, the consumption of solvent for mobile phase can be reduced up to 5 to 10-fold in savings.

Therefore, our study is focused to develop the analytical method for the capsaicin content in foods, especially *Gochujang*, using u-HPLC (LaChromUltra System, Hitachi-High Technologies Corporation, Japan) and to validate the new method according to single laboratory validation, including linearity, limit of detection (LOD) and limit of quantitation (LOQ), accuracy, precision, repeatability, and applicability.

Table 2. Analytical instruments for the quantification of capsaicinoids in foods

Instrument	Stationary phase	Mobile phase	Detector	Foods	Reference
Spectrophotometry	N.A.	N.A.	UV/VIS at 590 nm	Hot pepper	I.Perucka <i>et al.</i> Food Chemistry 71, 287-291, 2000
Gas chromatography	medium polar capillary column (0.25 mm×30 m)	Helium	FID	Gochujang	Ha, J <i>et al.</i> , J.AOAC Int. 92(2), 387-391, 2008
Liquid chromatography	C18 (4.6 × 15 cm)	Acetonitrile/ 1% Acetic acid in water (6:4)	UV 280 nm FLD280/325	Rep pepper oleoresins Hot pepper	AOAC 995.03, 2005


Figure 1. Chromatograms of capsaicin in standard (A) and *Gochujang* (B) by GC method.

Experimental

Material

Gochujang were purchased from the super market in Seongnam city, Korea, and stored at below 4°C in an air-tight container prior to analysis. All other samples for the applicability were also purchased from the super market in Seongnam city, Korea and should be mixed thoroughly so as to be a homogeneous mixture prior to analysis.

Sample Preparation

Accurately weigh ca 5 g *Gochujang* into 250 mL boiling flask. Add 50 mL 95% ethanol and several glass beads, and attach flask to reflux condenser. Gently reflux it for 5 hours and let it cooled. Filter 1 mL solution through a 0.20 µm syringe filter into small glass vial. Use it for LC analysis⁽⁹⁾.

Chromatographic conditions

Analysis of capsaicinoids was performed using LaChromUltra L-2000 U Series apparatus (Hitachi-High Technologies Corporation, Japan), including an eluant reservoir, HPLC pump, an autoinjection system of 2 µL injection at a fixed volume. Liquid chromatographic separation of the analytes was achieved using LaChromUltra C₁₈ (2 µm, 2 mm i.d. × 50 mm L, Hitachi-High Technologies Corporation, Japan). The column was equilibrated at a flow rate of 0.6 mL with a mobile phase consisting of acetonitrile-distilled water (40:60, v/v) containing

1% (v/v) acetic acid at 40-°C. L-2485U Fluorescence-Detector (Hitachi-High Technologies Corporation, Japan) under the condition of *Ex.* 280 nm, *Em.* 325 nm was used and the autosampler injection volume was set at 2 µL.

Reagents

Capsaicin (CAP) and *dihydrocapsaicin* (DHC) were purchased from Sigma (St. Louis, MO, USA). Ethanol (ACS reagent grade), acetonitrile (chromatographic grade), acetic acid (ACS reagent grade) were used.

Calibration Graph

Calibration graph for u-HPLC were based on peak area prepared by injecting 2 µL of the 0.1 0.5, 1.0 2.5 and 10.0 µg/mL solutions prepared through the dilution of capsaicin and dihydrocapsaicin stock solutions (1000 µg/mL) with 95% ethanol.

Method validation

The method validation has been carried out based on FDA guidance for bioanalytical validation (2001) and ICH Q2B Guideline using *Gochujang*.

Results and Discussions

Comparison of elution time

by conventional HPLC with u-HPLC

The comparison of chromatograms of capsaicinoids standard by conventional high performance

liquid chromatography (c-HPLC) and ultra high performance liquid chromatography (u-HPLC) was represented in **Figure 2**. The typical c-HPLC requires a considerable amount of analysis time (around 25 minutes) and solvent consumption for successful chromatographic analysis. However, u-HPLC requires a short analysis time (around 6 minutes) for the achievement of the same chromatographic analysis as c-HPLC. Therefore, it can be considered that u-HPLC method has been economical and environmentally friendly due to extremely rapid analysis. Because of the fast analysis, the consumption of solvent for mobile phase can be reduced up to 5-fold in saving.

Linearity

To determine the linearity, 5 different concentrations of capsaicin standards were used in a working range from 0.1 to 10.0 µg/mL for ultra high performance liquid chromatograph (u-HPLC)

method. Each solution was injected 3 times. Regression analysis revealed a good relation (correlation coefficient (r)=0.9997 for capsaicin and 0.9999 for dihydrocapsaicin). The graph showed a negligible intercept, which was calculated by the least-square method's regression equation (**Figure 3**).

Limit of detection (LOD)

and limit of quantitation (LOQ)

The limits of detection (LOD) and limits of quantitation (LOQ) values were estimated at an SD/b ratio of 3 and 10, where SD and b stand for the standard deviation of the intercept and slope of the regression line, respectively. LOD was 0.054 µg/mL for capsaicin and 0.053 µg/mL for dihydrocapsaicin. LOQ was 0.163 µg/mL for capsaicin and 0.160 µg/mL for dihydrocapsaicin (**Table 3**).

Recovery and accuracy

In order to study the accuracy of the methods, recovery experiments were carried out with the standard spiking a sample at three levels. The results were given in **Table 4**. The average recoveries obtained were quantitative (91.1-94.8% for capsaicin, 91.4-97.0% for dihydrocapsaicin), indicating good accuracy of u-HPLC methods.

Applicability

To apply the analysis method to various *Gochujang* products, chilli oil, *Kimchi* and snack containing a different amount of capsaicin, were analyzed with u-HPLC methods. The results showed that u-HPLC method can be applicable to the wide variety of foods containing chilli power and obtain the good separation without hindrance of interferences in foods as shown in **Table 5** and **Figure 4**.


Figure 2. Comparison of chromatograms of capsaicinoids standard by c-HPLC (upper) and u-HPLC (lower). Key : 1. Capsaicin, 2. Dihydrocapsaicin.


Figure 3. Calibration graph of capsaicin and dihydrocapsaicin by u-HPLC method.

Table 3. Calibration and sensitivity data of capsaicinoids

Component	Linear range (µg/mL)	R ²	LOD (µg/mL)	LOQ (µg/mL)
Capsaicin	0.10 - 10.0	0.9995	0.054	0.163
Dihydrocapsaicin	0.10 - 10.0	0.9999	0.053	0.160

Table 4. Intra- and Inter-day precision and accuracy data of the developed method for the determination of capsaicinoids in *Gochujang*^a

Components	Spiked amount (μg/g)	Intra-day (CV%)	Inter-day (CV%)	Recovery (%)
Cap	19.7	3.07	13.12	94.0
	39.4	4.51	4.53	94.8
	98.4	2.42	3.27	91.1
DHC	8.8	6.51	9.79	97.0
	17.5	5.81	5.18	93.4
	43.8	5.00	6.63	91.4

^a : sample weight 5.0 g, concentration of capsaicin and dihydrocapsaicin was 29.3 and 15.1 mg/kg, respectively.

Table 5. Applicability of the u-HPLC methods to various foods

Foods	Capsaicin	Dihydrocapsaicin
<i>Gochujang 1</i>	16.3±0.06	15.7±0.03
<i>Gochujang 2</i>	35.8±0.05	22.6±0.02
<i>Gochujang 3</i>	95.0±0.06	44.4±0.07
<i>Gochujang 4</i>	28.5±0.9	14.6±0.6
<i>Chilli oil</i>	290±0.5	146±0.5
<i>Kimchi</i>	2.29±0.15	1.30±0.06
<i>Snack</i>	4.50±0.02	3.64±0.06


Figure 4. Separation of capsaicinoids in various foods using u-HPLC.

Conclusions

For the analysis of capsaicinoids in foods, the u-HPLC method using LaChromUltra System (Hitachi-High Technologies Corporation, Japan) was validated in terms of precision, accuracy, and linearity according to FDA and ICH guidelines, and good results obtained within acceptable criteria.

Benefits of LaChromUltra System for capsaicinoids analysis in foods are reducing the analysis time to about 1/4, reducing the volume of solvent to about 1/10, remaining high resolution, good separation and easy handling.

It seems that LaChromUltra System can offer significant improvements in speed, sensitivity, and resolution compared with conventional HPLC, this bodes well for future applications in the analysis of capsaicinoids as well as in the analysis of food components.

Acknowledgement

We especially thank Hitachi-High Technologies Corporation for its financial support and also appreciate the technical advice received from Analytical Systems Design Department, Naka Division.

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Evaluation of Solar Cell Components using UV-VIS-NIR Spectrophotometer

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3.2 MuGFET (multi-gate field effect transistor) の観察

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H. Sakamoto(日立ハイテクノロジーズ)他: Determination of total arsenic in environmental sample with decomposition by UV-Irradiation using hydride generation atomic absorption spectrometry

【Summary】 It is well known that the sensitivity of arsenic measurement by hydride generation atomic absorption spectrometry (HG-AAS) depends on its chemical form. The measurement of the concentration of total arsenic often results in considerable errors in a sample containing various arsenic compounds. Therefore, the arsenic compounds need to be reduced to a homogeneous form. The optimum conditions for ultraviolet irradiation as a pretreatment have been investigated in this study. The main aspects examined were reagent concentrations and time of exposure to ultraviolet radiation. The sample

was mixed with $K_2S_2O_8$ and NaOH, and ultraviolet light was irradiated using a UV reactor (Model WA-UVR, KYORITSU CHEMICAL-CHECK Lab., Corp.). This pretreatment method was applied to the analysis of total arsenic in a standard reference material (The Japan Society for Analytical Chemistry, JSAC0302). Analytical results were in good agreement with the certified values. The recovery of each of the five arsenic compounds tested (methylarsonic acid, dimethylarsinic acid, trimethylarsine oxide, arsenobetaine and arsenocholine) was about 100%.

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