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Arsenic speciation with gradient hydride generation interfacing liquid chromatography and atomic absorption spectrometry

Yong Tian* Ming-Li Chen* Xu-Wei Chen* Jian-Hua Wang* Yoshihiro Hirano** Hideyuki Sakamoto** Ikumei Setsu**

This work was published in *Journal of Analytical Atomic Spectrometry*, 2010, 25, 48-54.

Introduction:

Arsenic in various sample matrixes causes a multitude of serious health problems.¹ Arsenic in drinking water in South Asia has become one of the biggest environmental problems of the century.² Generally, arsenic in environmental waters is almost exclusively inorganic,³ however, methylated arsenic species were frequently identified from biological samples and alga (Hijiki). Therefore, total amount of arsenic is no longer sufficient for providing comprehensive information about their toxicity or ease of removal among the various arsenic oxidation states, as both properties depend strongly on the chemical form of a specific arsenic species. Thus, the accurate identification of arsenic species has become one of the most important issues when evaluating the toxic effects of arsenic.⁴⁻⁶ Although a number of analytical protocols have been dedicating to arsenic speciation, however, there are no universal methodologies available, not to say providing standard procedures for such purposes. At this juncture, it is highly desirable to develop suitable and reliable approaches for arsenic speciation.

So far, extensive efforts have been devoted to the hyphenation of LC separation with ICP-MS.⁷⁻¹¹ LC provides powerful separation for isolating the individual arsenic species, while the high sensitivity of ICP-MS and its flow through feature facilitate accurate and fast detection. However, in many circumstances, this approach is not suitable for routine analysis attributed to the high running cost of the instrumentation. For this purpose, studies have been dedicated to coupling LC with hydride generation and some other atomic spectrometric detectors, e.g., atomic absorption spectrometry (AAS)¹²⁻¹⁵ and atomic fluorescence spectrometry (AFS),^{16,17} for arsenic

speciation.

The hyphenation of hydride generation with atomic absorption spectrometry is among the most frequently used techniques for the separation and quantification of arsenic species.¹⁸ There are, however, some intrinsic shortcomings of this technique when compared to those with detection by ICP-MS and AFS, e.g., the sensitivity of AAS can hardly compete with the latter systems. In order to be competitive, a significant improvement on the sensitivity is required with AAS as detector.¹²⁻¹⁵ When considering the inherent detection power of AAS limited by its intrinsic design, one of the most feasible approaches to improve the sensitivity is to ensure a complete capture of the analyte of interest at each single step of reaction. In arsenic speciation by using hydride generation, the conversion efficiencies of arsenic species into corresponding arsenic hydrides generally vary significantly. Although fixed experimental condition might be well suitable for certain arsenic species, while the conversion efficiencies for some others could be seriously deteriorated.¹²⁻¹⁴ A pH-gradient hydride generation approach has been proposed for the improvement of hydride generation efficiency.¹⁹ On the other hand, as an indispensable component in the hyphenation of a hydride generation unit with an atomic spectrometric detector, the gas-liquid separator is an extremely important part when desiring to improve the sensitivity and stability. Recently, this issue has caught pertinent attentions.^{17,20-22} Generally, a sufficient chamber volume in the gas-liquid separator is required to provide a buffering for the hydride and the carrier gas in order to depress the noise signal. On the other hand, however, the chamber volume should be limited to the extent possible to minimize the dispersion of the hydride before it entering the atomization area.²³ In this respect, it still deserves more efforts to be

Jian-Hua Wang

* Research Center for Analytical Sciences, Northeastern University, Box 332, Shenyang 110004, China

** Hitachi High-Technologies Corporation, Hitachinaka, Ibaraki, 312-8504, Japan

focused on this issue to further improve the performance of gas-liquid separation.

This work reports a gradient hydride generation protocol to improve the conversion efficiency of As (V), As(III), MMA, DMA and TMAO into their corresponding hydrides before AAS detection. In addition, a novel design of gas-liquid separator with limited dead volume was used providing minimum dispersion of the hydrides during their transportation from the separator to the detection point and thus facilitates high detection sensitivity. Arsenic speciation in Hijiki samples was performed by this procedure.

Materials and Methods

Instrumentation

AL-2130 pump and a L-2300 column oven (Hitachi High-Technologies) with Wakopak Navi C30-5 columns (F6×250 mm and F6×150 mm) were used for chromatographic separation of arsenic species. A six-port injection valve with a 100 μ L sample loop was used for sample introduction. A Z-2000 atomic absorption spectrophotometer (Hitachi High-Technologies) with a quartz atomizer was used for arsenic detection.

A Minipuls evolution peristaltic pump (Gilson) was used for delivering HCl and NaBH₄ solutions in the hydride generation process. A 60 °C water bath was used to improve the hydride generation efficiency of the arsenic species, followed by a cooling water bath to control the temperature of the reaction mixture before entering the gas-liquid separator. Four types of gas-liquid separator (**Figure 1**) were employed

and their separation efficiencies evaluated. The hydrides after separation were swept by argon to flow into the T-type quartz atomizer for arsenic quantification.

Figure 2 illustrated the flow manifold of the liquid chromatographic separation followed by gradient hydride generation-AAS detection for arsenic speciation.

Reagents

The chemicals used were at least of analytical-reagent grade with 18.2 M Ω cm de-ionized water.

Arsenic stock solutions of 1000 mg/L were prepared for sodium arsenite, disodium hydrogen arsenate heptahydrate, monomethylarsonic acid (MMA), dimethylarsinic acid (DMA) and trimethylarsine oxide (TMAO).

The mobile phase contains sodium 1-butanedisulfonate, malonic acid, tetramethylammonium hydroxide, methanol and ammonium tartrate. NaBH₄, HCl and L-cysteine were used for hydride generation.

Sampling and sample pretreatment

Hijiki digestion with nitric acid for total arsenic: 0.5 gram of Hijiki was randomly taken from the homogenized sample bulk. The sample was transferred into a 10 mL quartz glass crucible with 5 mL HNO₃ (65 %) and the mixture was heated gently on a hot plate by controlling the temperature not to exceed 100°C for about 3 hours until the solution had nearly dried. After cooling, the mixture was filtered through a 0.22 μ m membrane filter and the filtrate


Figure 1. The designs of gas-liquid separators used in the investigations.


Figure 2. The flow manifold of the LC separation-gradient hydride generation-AAS system. LC: liquid chromatography; FAAS: quartz atomizer flame atomic absorption spectrometer; P: peristaltic pump; SV: selector valve; RC: reaction coil; GLS: gas-liquid separator; Ar: argon; W: waste.
Valve position I: 6 mol/L HCl for the hydride generation of As(V), As(III) and MMA.
Valve position II: 0.6 mol/L HCl for the hydride generation of DMA and TMAO.

was diluted to 50.0 mL.

Ultrasonic-assisted water extraction for arsenic species: 5 gram of homogenized Hijiki was taken into a 100 mL vessel with 50 mL water. The mixture was sonicated for 30 min and the leachate was collected by filtration which was used for arsenic speciation.

Total arsenic quantification

Electrothermal atomic absorption spectrometry (ETAAS) was used for quantifying the total arsenic with 10 μ L of $\text{Mg}(\text{NO}_3)_2$ solution (1000 mg/L) as matrix modifier. In practice, 50-fold and 500-fold dilutions with HNO_3 (1 %, v/v) were made when analyzing total arsenic in the wet acid digests with HNO_3 and the individual arsenic species in the ultrasonic-assisted water leachate.

LC separation-gradient hydride generation-AAS system for arsenic speciation

Figure 2 illustrated the flow manifold. The system was initiated as soon as 100 μ L sample solution was injected into LC. The effluent from the LC system was directed to meet streams of NaBH_4 and HCl solutions to initiate hydride generation, which was significantly accelerated when the reaction mixture flowing through the reaction coil in a heated water bath at 60 $^\circ\text{C}$. After cooling down in a water bath, the reaction mixture was directed into the gas-liquid separator where the hydrides were isolated and further transported into the quartz atomizer for arsenic quantification. The gradient hydride generation was performed as such that a flow of 6 mol/L HCl was used for the hydride generation of As(V), As(III) and MMA, and when the signal of MMA has been

recorded, the selector valve was switched to introduce a stream of 0.6 mol/L HCl as hydride generation medium for the other arsenic species, i.e., DMA and TMAO. A fixed NaBH_4 concentration of 1.0 % (m/v) was used.

Results and discussion

The chromatographic separation process

The effect of the mobile phase constitutes on the absorbance and the retention time for each arsenic species were investigated, including the concentrations of sodium 1-butanesulfonate, tetramethylammonium hydroxide, methanol, malonic acid, ammonium tartrate and the acidity of the mobile phase.

The hydride generation process

Gas-liquid separation: The gas-liquid separation efficiency plays a vital role in the evaluation of hydride generation process. Four types of gas-liquid separators were evaluated (Figure 1), which provide 9 possible scenarios for designing the LC-GHG-AAS flow manifold. Figure 3 illustrated three different flow manifold designs. Generally, the separators II, III, IV receive the reaction mixture from the hydride generation process swept by an argon flow, and the argon stream can be introduced either right after the confluence of the fluidic zones or just before their entering into the separator (Figure 3B, 3C). In addition, when separator I was used in the flow manifold B and C, the bottom argon inlet was sealed. The separation of hydrides by separators II and III were not complete, and significant peak broadening was observed when using III and IV. Also, obvious fluctuations on the blank were recorded which deteriorate

the precision and detection limit. On the other hand, separator I provides an opportunity to introduce the argon stream into the separator via its bottom gas inlet. In practice, the liquid-seal above the out-let significantly reduced the dead volume for the gas phase which minimized the dispersion of the separated

hydrides and thus avoided signal tailing, i.e., much steeper signal peaks were recorded which facilitated high sensitivity.²³ In addition, the introduction of argon via the bottom of the separator well confined the fluctuation of the flow and created a much smoother blank signal which improved the precision.


Figure 3. The flow manifolds employed for evaluating the performance of the four types of gas-liquid separators.

FAAS: quartz atomizer flame atomic absorption spectrometer; P: peristaltic pump; RC: reaction coil; GLS: gas-liquid separator; Ar: argon; W: waste.

Table 1. The performance (characterized by peak height of As(V), peak half-width and relative standard deviation) of the system by using different designs of gas-liquid separators in various flow manifolds.

	Separator I			Separator II			Separator III			Separator IV		
	h	$Y_{1/2}$ (s)	RSD (%)	h	$Y_{1/2}$ (s)	RSD (%)	h	$Y_{1/2}$ (s)	RSD (%)	h	$Y_{1/2}$ (s)	RSD (%)
Flow manifold A	1	6.9	2.5	-	-	-	-	-	-	-	-	-
Flow manifold B	0.96	6.3	3.4	0.93	6.4	3.4	0.60	11.0	4.6	0.50	17.9	4.1
Flow manifold C	0.65	7.2	5.4	0.77	7.0	5.2	0.55	11.6	7.5	0.46	17.0	3.6

Note: HCl: 6 mol/L (with 2% L-cysteine), NaBH_4 : 1% (in 0.5% NaOH); A same flow rate of 1.1 mL/min was used for both HCl and NaBH_4 solutions. The 100 μL sample solution contains 100 $\mu\text{g/L}$ of As(V) at a sampling flow rate of 1.1 mL/min. The water bath was set at 60 $^{\circ}\text{C}$. h: normalized peak height; $Y_{1/2}$: peak half-width.

The results for As(V) in **Table 1** showed favorable analytical performance by incorporating separator I into the flow manifold of **Figure 3A**.

Gradient hydride generation: General conditions for hydride generation were first scrutinized with 0.5%–3% NaBH₄ (m/v, 0.5% NaOH) and 0.2–6.0 mol/L HCl. Within this concentration range of HCl, substantial increments of the absorbance signals were recorded for the 5 arsenic species when increasing the NaBH₄ concentration up to 1%. Thereafter, higher concentration of NaBH₄ not only caused declines of the signal, but also generated more noise for the blank.

Figure 4 illustrates the dependences of hydride generation on HCl concentration. For As(V), As(III) and MMA, continual increase of the signal was observed when enhancing the acid concentration within the whole range of 0.2–6 mol/L. For DMA and TMAO, however, a rapid increase of the signal was first encountered by increasing the acid concentration up to 0.6 mol/L where maximum signals were achieved, while a further increase leads to decrease of the signal, especially for TMAO whose signal sharply dropped to a negligible level. These observations could be attributed to the protonation of a methyl-arsine, e.g., (CH₃)₂ AsH from DMA and (CH₃)₃ As from TMAO, which can readily occur in an acid medium and the produced water-soluble

methyl-arsonium species in equilibrium with gaseous methyl-arsine reduces the hydride generation efficiency,^{24, 25} corresponds to a rapid drop of the signal arising from DMA and TMAO. The much more significant drop of the signal from TMAO is explained by the favorable protonation of (CH₃)₃ As than that of the (CH₃)₂ AsH due to its higher basic nature.

Figure 4 indicated that each arsenic species requires a specific acid concentration for optimal efficiencies of hydride generation.^{26, 27} For 6 M HCl at which maximum signals were recorded for As(V), As(III) and MMA, the signals from DMA and TMAO corresponds to only 20% and 10% or less of the signal as obtained from a same amount of As(V). Thus, it is unable to find a suitable acid concentration at which an identical or similar response could be achieved for the 5 arsenic species. In practice, a compromise acid concentration is generally selected which might be favorable for the hydride generation of some arsenic species, but inevitably sacrifices the hydride conversion efficiency for other species.^{12–14, 17, 28–30} In this respect, gradient hydride generation should be favorable to improve the efficiency of hydride generation for all the arsenic species. In the present case, gradient hydride generation is performed such that a 6 M of HCl was used for the hydride generation of As(V), As(III) and MMA, when the MMA signal has been recorded, a


Figure 4. The dependences of hydride generation efficiencies on the concentration of HCl. 100 µg/L of each arsenic species was introduced. 2% L-cysteine (m/v) was included in the HCl solutions. A flow rate of 1.1 mL/min was used for both HCl and NaBH₄ solutions.


Figure 5. The chromatograms recorded for arsenic speciation by LC separation-hydride generation-AAS detection. A and B were obtained by conventional hydride generation with a fixed HCl concentration for the 5 arsenic species, while C was achieved by a gradient hydride generation, i.e., a 6.0 M HCl was used for the hydride conversion of As(V), As(III) and MMA, while a 0.6 M HCl for those of DMA and TMAO.

stream of 0.6 mol/L of HCl as hydride generation medium was introduced for the other two arsenic species, i.e., DMA and TMAO. **Figure 5** clearly showed that the retention time interval between MMA and DMA is long enough for performing the gradient hydride generation. During this process, a fixed concentration of NaBH₄ at 1.0 % (m/v) was used for the 5 arsenic species.

Figure 5 illustrates the chromatograms recorded by gradient hydride generation employing a gradient of acid concentration, i.e., 6 M HCl for As(V), As(III) and MMA, while 0.6 M HCl for DMA and TMAO. As a comparison, those results achieved by a conventional hydride generation protocol with compromised HCl concentrations, i.e., 0.6 M and 6 M HCl for all the 5 arsenic species, were included. Much improved signals for the 5 arsenic species were recorded, ensuring both favorable hydride generation efficiency and detection sensitivity for these arsenic species. It is also worth mentioning that a maximum of 13-fold improvement on the detection sensitivity of As(V) was achieved in comparison with the previous reports.¹²⁻¹⁵

L-cysteine has been widely used for reducing arsenic species prior to hydride generation. In this study, 2 % L-cysteine (m/v) in HCl solution further improves the hydride generation efficiency and the detection sensitivity for the 5 arsenic species up to 16 %, 19 %, 24 %, 13 %, 32 % for As(V), As(III), MMA, DMA, TMAO, respectively. These results confirmed the previous observations.³¹ A most likely explanation might be that in aqueous medium the formation of borane complexes of X-BH₃⁻ (X = SR), which are more active than BH₄⁻ in the hydride transfer to analytical substrates, is possible in the case of hydrolysis of BH₄⁻ in the presence of thiols in L-cysteine.³²

The temperature dependence of hydride generation efficiency: We have also investigated the temperature dependence of hydride generation efficiency of the 5 arsenic species within 20-100 °C. A gradual increase of the hydride conversion efficiency for As(V) and As(III) was observed in the whole range studied. For MMA, DMA and TMAO, however, an increment was achieved when enhancing the temperature from 20-60 °C, while thereafter a decline of the signal was recorded. The hydride generation was performed at a compromised temperature of 60 °C.

The dimension of the reaction coil: The diameter and length of the reaction coil should be optimized to

achieve maximum signal. 1.0 mm i.d. PTFE tubing well facilitates the conversion of arsenic species into hydrides when fixing the other parameters. When evaluating the effect of the reaction coil length, a slight improvement of the signal was observed by increasing its length within a certain range. However, an excessive long reaction coil created unexpected flow impedance, causing unstable flow rates for the acid and NaBH₄ zones during hydride generation which produced an exaggerated noise on the blank signal. Finally, a reaction coil length of 3500 mm was used.

Analytical performance

By using a 100 µL sample volume, a complete analytical run could be finished within 13.5 min. The limits of detection for the 5 arsenic species, i.e., As(V), As(III), MMA, DMA and TMAO, were 0.9, 1.4, 1.4, 1.6 and 1.5 µg/L (As), respectively. In addition, RSD of less than 3 % and 6 % were obtained at 100 µg/L and 20 µg/L.

Table 2 compares the characteristic performance data of the present system with some of the relevant procedures based on LC separation hyphenated with detection by mass spectrometry and atomic spectrometry. A significant improvement on the overall sensitivities for the 5 arsenic species was achieved by the present system incorporating gradient hydride generation when compared to the reported similar procedures with detection by AAS¹²⁻¹⁵ and ICPAES.^{29, 30} Although the limits of detection of this procedure cannot compete with those achieved by ICPMS^{11, 33}, it is sufficient for routine arsenic speciation in environmental water and seafood or seaweed samples.

Arsenic speciation in real samples

Total arsenic concentrations in the wet acid digest with HNO₃ and the ultrasonic-assisted water leachate of Hijiki were quantified by ETAAS after appropriate dilution (**Table 3**). The ultrasonic-assisted water extraction efficiencies of total arsenic in Hijiki 1 and Hijiki 2 were ca. 76.1 % and 82.2 % respectively. A longer extraction time has been used but no significant improvement on the extraction efficiencies was obtained. Thus, the low extraction efficiency might be attributed to the presence of some unextractable arsenic species in the Hijiki samples.

The speciation of the 5 arsenic species in Hijiki was performed by employing the present LC separation -

Table 2. Comparisons of the performance data of the present procedure with some of the arsenic speciation protocols using LC separation coupled to ICPMS, HG-AFS, HG-ICPAES and HG-AAS.

Method	LODs ($\mu\text{g/L}$)					RSD (%)	Sample volume (μL)	Ref.
	As(V)	As(III)	MMA	DMA	TMAO			
LC-ICPMS	0.1	0.02	0.04	0.06	-	10	200	11
	0.3	0.1	0.2	0.1	-	2.7	50	33
LC-HG-AFS	0.2	0.1	0.2	0.1	-	10.3	200	16
	8.5	1.0	2.2	2.3	-	2.2	20	17
LC-HG-ICPAES	20	5	10	20	-	4.8	100	29
	15	10	30	65	-	5.9	200	30
LC-HG-AAS	2.6	2.4	2.4	2.3	1.1	6.4	100	12
	8.4	1.5	3.3	4.3	-	3.7	100	13
	4.2	2.0	2.9	2.4	-	8.7	100	14
	12.0	7.8	-	-	-	12.1	200	15
	0.9	1.4	1.4	1.6	1.5	2.8	100	This work

Table 3. The results of arsenic speciation in Hijiki samples and the corresponding spiking recoveries.

Species	Found Hijiki 1 ($\mu\text{g/g}$)	Found Hijiki 2 ($\mu\text{g/g}$)	Spiked ($\mu\text{g/L}$)	Spiking Recov. Hijiki 1 (%)	Spiking Recov. Hijiki 2 (%)
As(V)	53.6 ± 9.2	60.8 ± 13.2	50	107.4	111.2
As(III)	n.d.*	n.d.	50	97.9	97.1
MMA	n.d.*	n.d.	50	106.9	104.9
DMA	2.7 ± 0.4	2.4 ± 0.5	50	112.5	105.7
TMAO	0.34 ± 0.12	0.15 ± 0.03	50	104.6	102.2
Sum of the arsenic species determined	56.6 ± 9.7	63.4 ± 13.7			
Total arsenic (water leachates)	63.9 ± 3.5	63.8 ± 1.9			
Total arsenic (wet acid digest)	84.0 ± 8.2	77.6 ± 4.2			

Mean value \pm 95 % confidence interval (number of replicates = 3)

* n.d.: not detected.


Figure 6. The chromatograms for arsenic speciation in Hijiki. The ultrasonic-assisted water leachates were 5-fold diluted.

gradient hydride generation-AAS system. After appropriate dilutions of the ultrasonic-assisted water leachates, large amounts of As(V) were found from both Hijiki samples. The chromatograms were illustrated in Figure 6 and the speciation results were given in Table 3. For Hijiki 2, a reasonable agreement was achieved between the sum of the three individual arsenic species identified (As(V), DMA and TMAO), i.e., $63.4 \mu\text{g/g}$, and the total arsenic concentration quantified by ETAAS, i.e., $63.8 \mu\text{g/g}$. However, a discrepancy was observed for Hijiki 1, which attributed probably to the presence of some other arsenic species, e.g., arsenosugars, except for those investigated in the present work.³⁴ In addition, favorable recoveries were obtained for both samples by spiking appropriate amount of the arsenic species.

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Analysis of environmental samples by atomic absorption spectrometry with pretreatment using polyamino-polycarboxylic acid type chelating resin

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H.Nakagawa(日立ハイテクノロジーズ)他: Post-column labeling using HPLC for food analysis

【Summary】 Post-column labeling using high performance liquid chromatography (HPLC) is powerful method for analyses of samples containing contaminants, such as foodstuff, or components without a specific UV absorbance or fluorescence. The analytical system is shown in **Figure 1**. Ultra-HPLC (UHPLC) combines high-speed analysis with high separation efficiency through the use of fine silica particles and tubes. It was

believed that the post-column method, in which another pump is connected after separation by the column, is incompatible with high-separation by UHPLC. We realized a post-column method using UHPLC for organic acid analysis at the first in the world. Here we improved post-column labeling by UHPLC for amino acid analysis using OPA as a labeling reagent.


Figure 1. Analytical system of post-column method

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6. 58th ASMS Conference on Mass Spectrometry and Allied Topics(2010/5/22～29 USA)

M. Watanabe(日立ハイテクノロジーズ)他: Peptide identification using electroncapture dissociation combined with subsequent radical-ion/molecule reactions on a chromatographic time scale

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