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## Combination of Dynamic Hollow Fiber Liquid-phase Microextraction with a Hitachi L-2000 HPLC for the Determination of UV Filters in Cosmetic Products

Hongyun Yang<sup>1,2</sup>, Haifang Li<sup>2</sup>, Masahito Ito<sup>3</sup>, Jin-Ming Lin<sup>2</sup>, Guangsheng Guo<sup>1</sup> and Mingyu Ding<sup>2</sup>

### 1. Introduction

Sunscreen cosmetic products are widely used in the world to prevent radiations emitted by solar radiation, and especially UV radiation (UVR) which lead to sunburns and carcinogenesis<sup>[1, 2]</sup>. Sunscreen cosmetic products generally are prepared as creams, lotions, emulsions or sprays. A variety of UV filters are added into the cosmetics as strong absorbers, reflectors or scatters of UVR<sup>[3-6]</sup>. It has been documented that organic UV filters can result in some dermatological reactions during exposure to the sunlight<sup>[7]</sup>, and become environmental pollutants when they are discharged into the living waste water<sup>[8-10]</sup>. Lists of those approved UV filters with their maximum allowed concentrations in commercial products have been set by regulatory authorities in USA, Europe and Japan<sup>[11-13]</sup>. Therefore, it is necessary and urgent to establish convenient quantification methods for determination of UV filters in cosmetic products.

In recent years, various pretreatment and analytical methods have been used to determine UV filters<sup>[14, 15]</sup>. Based on the literature, liquid-liquid extraction (LLE)<sup>[16, 17]</sup> and solid-phase extraction (SPE)<sup>[18-20]</sup> as the most commonly sample pretreatment methods have been used to extract UV filters. Other extraction techniques such as cloud-point extraction<sup>[21]</sup>, supercritical fluid extraction<sup>[22]</sup>, microwave assisted extraction<sup>[23]</sup> have also been proposed for concentration of UV filters. The extraction pretreatments always combine with high performance liquid chromatography (HPLC)<sup>[21-27]</sup>, liquid or gas chromatography coupled to mass spectrometry (GC-MS)<sup>[16, 21, 28-30]</sup> or size-exclusion liquid chromatography<sup>[31]</sup>.

With recent development of microextraction

Prof. Jin-Ming Lin

Dr. Haifang Li

techniques, solid phase microextraction (SPME)<sup>[32]</sup> and liquid phase microextraction<sup>[33]</sup> in terms of high speed, high sensitivity and convenient have been recommended to be more convenient techniques than LLE and SPE. The most widely used extraction materials of SPME are fragile and expensive. Furthermore, carry-over of sample is sometimes difficult to be eliminated for extraction and injection performed in different steps. Single-drop LPME (SDME) methods<sup>[24]</sup> based on a droplet of a few microliters of organic solvent at the tip of a syringe needle have become the most extended approaches. However, the droplet easily dislodges from the needle tip of the syringe during extraction, especially in the case of vigorous stirring.

Hollow fiber liquid-phase microextraction (HF-LPME) technique based on the principle of membrane liquid microextraction<sup>[34]</sup>, is employed as an alternative liquid microextraction. The micropores of hollow fiber can confine organic liquid solvent inside of the hollow fiber, which prevent solvent loss during extraction. The large molecules and extraneous materials in sample can be prevented due to the selective transmission of the micropore membrane. Furthermore, micropore structure increases a large specific surface and supports a more amount of extraction phase volume. HF-LPME is an efficient sampling pretreatment method with many advantages such as simplify, high efficiency and low cost. The HF-LPME can be classified into static HF-LPME and dynamic HF-LPME. The static HF-LPME is performed by jointing a hollow fiber on a microsyringe

needle and immersing into a stirred sample solution<sup>[35]</sup>. Dynamic HF-LPME was developed by continuous injection of sample solution into the extractor to well contact with the acceptor phase, which contribute to high extraction efficiency and shorten extraction time. Ouyang et al<sup>[36]</sup> has introduced the kinetic calibration for automated HF-LPME.

In the present study, we developed a simple dynamic HF-LPME device for pretreatment of five UV filters in cosmetic products, which followed by HPLC analysis. The cannular extractor was assembled to realize continuous-flow microextraction. The sheath structure of the extractor made the sufficient contact between the acceptor phase and aqueous sample and led to high extraction efficiency. The main parameters of dynamic HF-LPME were investigated and optimized. Under the optimized HF-LPME conditions, the proposed dynamic HF-LPME presented good extraction efficiency for five UV filters in cosmetic products.

## 2. Experimental

### 2.1 Reagents and materials

2-Hydroxy-4-methoxybenzophenone-5-sulphonic acid (BZ4, >96%), oxybenzone (BZ3, >98%), 2-ethylhexyl-4-dimethyl aminobenzoate (ODP, >98%), 2-ethylhexyl-4-methoxy-cinnamate (OMC, >99%), 2-ethylhexyl salicylate (ES, >99%) were purchased from Sigma-Aldrich (Steinheim, Germany) and were used as standards. The chemical properties of the standards are listed in Table 1.

The solvents of toluene, n-hexane, propylene glycol monomethyl ether acetate (PMA), and dichloromethane were of HPLC grade and were obtained from Sigma-Aldrich (Phillipsburg, NJ, USA). HPLC-grade methanol from J.T. Baker (Phillipsburg, NJ, USA) and de-ionized water purified by a Milli-Q ultrapure purification system (Millipore, Bedford, MA, USA) were used as mobile phase in HF-LPME. Acetone and other reagents were purchased from Beijing Chemical Works at analytical grade.

The Q3/2 Accurel polypropylene (pp) hollow fiber membrane (600  $\mu\text{m}$  i.d., 200  $\mu\text{m}$  wall thickness, and 0.2  $\mu\text{m}$  pore size) was purchased from Membrane GmbH (Wuppertal, Germany). Disposable flow control valve lines for the visiprep<sup>TM</sup>-DL (1.0 mm i.d., 500  $\mu\text{m}$  wall thickness, and 50 mm length) from Supelco (Bellefonte, USA) were used as external tube. The 10  $\mu\text{L}$  syringe (Shanghai, China) was attached with the hollow fiber to fill in and recover

| Name <sup>a</sup>                           | Structure | Log $K_{ow}$ <sup>b</sup> | $pK_a$ <sup>b</sup> |
|---------------------------------------------|-----------|---------------------------|---------------------|
| Benzophenone-4 (BZ4)                        |           | 0.89                      | -0.7                |
| Benzophenone-3 (BZ3)                        |           | 3.79                      | 7.56                |
| 2-Ethylhexyl-4-dimethyl aminobenzoate (ODP) |           | 6.15                      | 2.39                |
| 2-Ethylhexyl-4-methoxy Cinnamate (OMC)      |           | 5.80                      | -                   |
| 2-Ethylhexyl salicylate (ES)                |           | 5.97                      | 8.13                |

<sup>a</sup> Name: International Nomenclature for Cosmetic Ingredient (INCI) elaborated by CTFA and COLIPA.

<sup>b</sup> Values obtained from SciFinder Scholar Database, <http://www.cas.org/products/sfacad/>.

the acceptor phase. The 1.0 mL syringe was used to continuously introduce the sample into the extractor. Two micro syringe pumps (Baoding, China) were used to inject the solutions during extraction.

### 2.2 Chromatography

A Hitachi L-2000 HPLC system, supplied by Hitachi High-Technologies (Tokyo, Japan), consisted of a wavelength UV detector (Hitachi L-2400) with a standard flow cell, a quaternary pump with online degasser (Hitachi L-2130), an autosampler (Hitachi L-2200) and a column oven (Hitachi L-2300). The separation was performed on a reversed phase  $C_{18}$  HPLC column (Lichrospher ODS, 5.0  $\mu\text{m}$ , 4.6 mm  $\times$  250 mm, Merck Millipore, China). The mobile phase was a mixture of methanol-pure water (80:20, v/v). The five filters were separated within 25 min at flow rate of 0.8 mL  $\text{min}^{-1}$ . The column oven was controlled at 40  $^{\circ}\text{C}$  and UV detector was set at 310 nm.

### 2.3 Preparation of standard solutions and cosmetic samples

The stock solution was prepared in methanol weekly. The concentrations of BZ4, BZ3, ODP and OMC were 1000 mg  $\text{L}^{-1}$ , while ES was 5000 mg  $\text{L}^{-1}$ , respectively. And the stock solution was stored in the dark at -15  $^{\circ}\text{C}$  to minimize photodegradation of the

UV filters. The working standard solution (0.1–500 mg  $\text{L}^{-1}$ ) was prepared daily with appropriate dilution of the stock solution with methanol.

All of the cosmetic samples were purchased from local markets. Sample of about 0.010 g was dissolved in 2.0 mL n-hexane completely through ultrasonic at ambient temperature for 5 min, then centrifugated at 5000 rpm for 3.0 min. The supernatant was filtered through a 0.22  $\mu\text{m}$  membrane and nitrogen dried at 60  $^{\circ}\text{C}$ , and then diluted to 100 mL with de-ionized water. The prepared sample solution was stocked for the HF-LPME procedure.

### 2.4 Dynamic HF-LPME device

The extraction device used in this study consisted of external tuber, hollow fiber and microsyringe pump. Fig. 1 shows schematic illustration of the HF-LPME extraction system. Q3/2 Accurel polypropylene hollow fiber was cut into 3.7 cm hollow fiber segments. Those hollow fiber segments ultrasonically cleaned in acetone for 15 min and dried in the air. A disposable flow control valve line for the visiprep<sup>TM</sup>-DL as external tube was installed on the sample injection syringe. Then aqueous sample of 0.6 mL was loaded into the syringe. The hollow fiber attached onto the microsyringe needle was inserted into the visiprep<sup>TM</sup>-DL external tube, and then 10  $\mu\text{L}$  acceptor phase was filled into the hollow fiber. The sample was continuously injected into the extractor


Figure 1 Schematic illustration of HF-LPME method. The cannular extractor was assembled by mounting a hollow fiber inside an external tube with a tee-connector.

by the pump. During the extraction, the analytes in the aqueous sample were largely extracted into the organic solvent by diffusion. The analyze-enriched acceptor phase was directly collected into the micro-syringe after the extraction. Finally, 5.0  $\mu\text{L}$  acceptor phases was used for HPLC analysis. Different acceptor phases and different flow rates were investigated. The extraction was preformed at ambient temperature (25  $^{\circ}\text{C}$ ).

## 3. Results and Discussion

This dynamic HF-LPME relying on membrane liquid microextraction is an excellent technique for extraction. The organic acceptor phase was immobilized in the micropores of the hollow fiber to form liquid membrane. The targets were extracted from the aqueous sample matrix into the organic liquid membrane on the wall of hollow fiber, and then diffused into the lumen of the hollow fiber. The process of mass transfer is shown in Fig. 2.

In this study, in order to optimize the dynamic HF-LPME for the determination of five UV filters in cosmetic products, several parameters controlling extraction efficiency including types of extraction solvent, sample volume, sample flow rate, pH values and ionic strength were assessed. A series of 0.6 mL of sample solution spiked with 5.0 mg  $\text{L}^{-1}$  of ES and 1.0 mg  $\text{L}^{-1}$  of the other four UV filters were extracted in triplicate.



Figure 2 Schematic representation of HF-LPME method. The analyte diffuses across a microporous membrane while in contact with an organic acceptor.  $K_p$  is defined as the partition coefficient of the analyte between the aqueous and organic phase.

### 3.1 Effect of selection of acceptor phase


Organic solvent serving as acceptor phase was an essential consideration for HF-LPME. The organic solvent used in HF-LPME must satisfy several criteria. Firstly, the organic solvent should have good affinity for hollow fiber on which thin liquid membrane can be formed. Secondly, the organic solvent should possess a high distribution coefficient for the analytes and be immiscible with water. Thirdly, the organic solvent should provide appropriate extraction selectivity and high extraction recovery. Finally, the organic solvent should be compatible with HPLC analysis. Given on the above consideration, n-hexane, toluene, and 1-octanol, which have been the most commonly used extraction solvents in LPME, were estimated in our study. **Fig. 3** shows the performance of the organic solvents for extraction of five UV filters under identical conditions. The results indicated that toluene was the most suitable as the acceptor phase with the highest extraction efficiency. However, the extraction efficiencies of n-hexane, dichloromethane, 1-octanol and PMA were much lower than those of toluene. The good compatibility of the hollow fiber and the low solubility in water were contributed to a high enrichment factor for UV filters. Furthermore, the hollow fiber could become transparent when it is filled with toluene, which made the operation easier. Thus, toluene was chosen for the further experiments.

**Figure 3** Effects of different extraction solvent on HF-LPME method. Extraction conditions: 0.6 mL of sample containing 25 mg L<sup>-1</sup> of ES, 5.0 mg L<sup>-1</sup> of other four target analytes, extracted by toluene at 3.0 mL h<sup>-1</sup>, ambient temperature.

### 3.2 Effect of sampling volume

HF-LPME is a process dependent on equilibrium rather than exhaustive extraction. So extraction efficiency is normally controlled by the distribution coefficient of the analyte between the aqueous phase and organic acceptor phase. The adsorption capacity of organic phase was also an important parameter for 10 µL organic acceptor phase. Sample solution containing 25 mg L<sup>-1</sup> of ES, 5.0 mg L<sup>-1</sup> of the other four UV filters was extracted. The different sample volumes in the range of 0.5–2.0 mL were performed. As shown in **Fig. 4**, the extraction efficiency increases with sample volume increasing and reaches the maximum up to 1 mL volume. When the sample volume higher than 1 mL, the extraction efficiency was not further increased. This is due to that the extracted amount of UV filters was increased with the sampled volume increasing until the extraction reached balance. On the other hand, acceptor phase slightly percolated through the micropores of hollow fiber during the sample continuous injection, which may lead to a little decrease of extraction efficiency. The adsorption capacities of 10 µL toluene were reached to 4.0 µg of BZ4, 4.5 µg of BZ3, 4.7 µg of ODP and OMC, and 23 µg of ES, respectively.

**Figure 4** Effect of extraction sample volume on the HF-LPME method. Extraction conditions: sample containing 25 mg L<sup>-1</sup> of ES, 5.0 mg L<sup>-1</sup> of other four target analytes, extracted by toluene at 3.0 mL h<sup>-1</sup>, ambient temperature.

### 3.3 Effect of flow rate

Generally speaking, the extraction efficiency could be enhanced with the increasing extraction time before the equilibrium reached. The extraction time could be controlled by the injection flow rate

of the sample in this study. Therefore, the flow rate of the sample solution played a significant effect on the extraction of the dynamic HF-LPME. The flow rate in the range of 0.5–7.0 mL h<sup>-1</sup> was considered and the result showed that the extraction of five UV filters have no significant fluctuation at the range of low rate (<2.0 mL h<sup>-1</sup>). However, the extraction was curtailed when the flow rate was higher than 2.0 mL h<sup>-1</sup>. This result might be due to the fact that the extraction equilibrium could not be reached balance under the very high flow rate condition. Finally, the flow rate of 2.0 mL h<sup>-1</sup> was chosen for subsequent experiments.

### 3.4 Effect of salt concentration

The addition of salt to aqueous sample decreases the solubility of analytes in the sample [34] and increases their partition into the organic solvent through changing solvent environment. In this study, the effect of ionic strength on the extraction was evaluated by addition of different concentration of NaCl (0–20 %, w/v) solutions to the sample. The best extraction efficiencies were achieved without addition of the NaCl and the extraction efficiencies of five UV filters decreased dramatically as more NaCl was added. The results indicated that the addition of NaCl had a significant negative effect for extraction. The extraction efficiency of BZ4 decreased slightly with the ionic strength increase compared with other four UV filters. At concentration up to 20 %, BZ3 can not be extracted. This phenomenon could be assumed that apart from the “salting-out”, the physical properties of the organic membrane were changed with the presence of salt, which reduced the diffusion rates of the analytes into the acceptor phase. The viscosity of the sample was also changed at high concentration of salt, which affected the diffusion rate. Thus, the best conditions for the HF-LPME were performed without the addition of NaCl.

### 3.5 Effect of pH of sample

The extraction efficiency and selectivity of the HF-LPME can be optimized by the pH values of the sample solution. The effect of sample pH values on the extractability of the target analytes by HF-LPME was investigated from 1.0 to 8.0. The result indicates that the extraction efficiencies of BZ3, BZ4 and ODP were decreased significantly flowing with pH values increase. The pH values had positive effects on the extraction efficiency of OMC and ES. This

phenomenon was in accordance to the pK<sub>a</sub> and could be explained by the molecular properties of the five UV filters. According to the chemical structure of five UV filters (**Table 1**), most of them have the hydrophilic groups and the UV filters were protonated at high pH values, which make their partition much more into the aqueous phase and reduce their transportation into the organic solvent. Thus, UV filters display a high affinity for organic solvent in low pH values and good extraction efficiency was achieved in acidic pH environment. The chemical structure of OMC indicates that only hydrophobic groups existed can not protonated and the molecules of OMC and ES could hydrolyze easily in acidic pH environment, leading to low response. Then it was observed that the better extraction efficiency for five UV filters were accomplished at pH value of 2 after taking all factors into consideration. Therefore, pH value of 2 was maintained throughout all the work.

### 3.6 Method evaluation

As shown before, the best extraction efficiency was obtained under the following condition: 0.6 mL of sample solution at no salt addition in acidic environment, acceptor phase of toluene, sample injection flow rate at 2 mL h<sup>-1</sup>. Under the optimized conditions, limits of detection (LODs), linearity and relative standard deviations (RSDs) were investigated. The calibration curves were constructed by plotting the mean peak areas against respective concentration (mg L<sup>-1</sup>) as listed in the **Table 2**. Good correlation coefficients (R<sup>2</sup>) with excellent linearity for target analytes were higher than 0.9960.

**Table 2** Analytical characteristics of the method

| UV filters | Time (min) | Range (mg L <sup>-1</sup> ) | R <sup>2</sup> | LOD (ng L <sup>-1</sup> ) <sup>a</sup> | RSD (%) <sup>b</sup> | F <sup>c</sup> |
|------------|------------|-----------------------------|----------------|----------------------------------------|----------------------|----------------|
| BZ4        | 2.31±0.01  | 0.5-500                     | 0.9978         | 1.00                                   | 4.5                  | 24             |
| BZ3        | 10.67±0.04 | 0.5-500                     | 0.9960         | 20.00                                  | 5.1                  | 34             |
| ODP        | 14.99±0.01 | 0.5-500                     | 0.9960         | 15.00                                  | 5.2                  | 57             |
| OMC        | 17.73±0.03 | 0.5-500                     | 0.9970         | 10.00                                  | 4.8                  | 55             |
| ES         | 20.92±0.02 | 2.5-2500                    | 0.9965         | 100.00                                 | 5.0                  | 53             |

<sup>a</sup> S/N=3; <sup>b</sup> N=3; <sup>c</sup> F=Enrichment factor

The LODs (S/N=3) of the five UV filters were found in the range from 1 to 100 µg L<sup>-1</sup>. A blank sample was carried out every time to confirm that no carryover occurred when determining LODs. The RSDs (n=3) for the five UV filters were lower than 5.2 %. All the results were calculated and summarized in **Table 2**. Compared to the traditional SPE


method, the proposed method exhibited better selectivity with comparative sensitivity and recovery. The reason was that the undesirable large molecules can not cross through the micropores of the hollow fiber. To investigate the extraction recovery, the method of standard additions was applied to all sunscreen cosmetic by spiking with 5 mg L<sup>-1</sup> ES and 1 mg L<sup>-1</sup> of other analytes, respectively. The relative recoveries of the five UV filters ranged between 80 % and 102 %.

### 3.7 Extraction of UV filters in cosmetic products

In order to test the applicability of the HF-LPME method to real sample, the proposed procedure was applied to analysis the cosmetics purchased from local super market. After the pretreatment of original cosmetic samples, all the cosmetic products were extracted under the optimum HF-LPME conditions and the existent amount of the UV filters in the cosmetics were determined and listed in **Table 3**. The OMC was contained in all the test cosmetic products. BZ3 and ES were presented in most of the test cosmetic products. Whereas, ODP was not found in some analyzed samples. The determined results were in the accordance to the service condition of the UV filters. Hence, HF-LPME method enables the precise determination of UV filters in cosmetic products. **Fig. 5** shows a chromatogram of a cream sample with and without spiking of 5 mg L<sup>-1</sup> ES and 1 mg L<sup>-1</sup> of the other four filters after HF-LPME extraction.

### 4. Conclusions

In the present study, a simple and effective dynamic HP-LPME method was introduced for determination of UV filters in cosmetic samples combined with HPLC analysis. High extraction efficiency was obtained by continuous injection of sample into the extractor to accelerate the extraction. This dynamic HF-LPME method with the advantages of rapid speed, low cost, and convenient operation were successfully applied to the investigation of UV filters in cosmetic products.

**Figure 5** Extraction of a blank cream and a cream spiked with 5.0 mg L<sup>-1</sup> ES and 1.0 mg L<sup>-1</sup> of other four target analytes by dynamic HF-LPME with HPLC-UV. (a) Cream sample analyzed by HF-LPME-HPLC-UV. (b) Spiked cream sample analyzed by HF-LPME-HPLC-UV. (1) BZ4, (2) BZ3, (3) ODP, (4) OMC, (5) ES.

**Table 3** Five UV filters found in real cosmetic samples after applying the proposed HF-LPME method

| Cosmetic products       | UV filters in cosmetics (g 100g <sup>-1</sup> ) |      |      |      |      |
|-------------------------|-------------------------------------------------|------|------|------|------|
|                         | BZ4                                             | BZ3  | ODP  | OMC  | ES   |
| Cream <sup>1</sup>      | N.D. <sup>a</sup>                               | 2.02 | N.D  | 2.36 | 4.35 |
| Cream <sup>2</sup>      | 2.56                                            | N.D  | 4.32 | 3.28 | N.D  |
| Cream <sup>3</sup>      | 2.79                                            | N.D  | N.D  | 1.58 | 3.12 |
| Cream <sup>4</sup>      | 0.60                                            | N.D  | N.D  | 4.42 | 2.03 |
| Cream <sup>5</sup>      | 0.54                                            | N.D  | N.D  | 4.20 | N.D  |
| Lotion <sup>1</sup>     | N.D                                             | 2.35 | N.D  | 2.46 | 2.82 |
| Lotion <sup>2</sup>     | N.D                                             | 3.42 | N.D  | 2.42 | N.D  |
| Lotion <sup>3</sup>     | 1.42                                            | 2.96 | N.D  | 3.62 | 4.18 |
| Lotion <sup>4</sup>     | 2.32                                            | 4.51 | N.D  | 3.88 | 4.06 |
| Lipstick <sup>1</sup>   | N.D                                             | N.D  | 0.8  | 4.38 | N.D  |
| Lipstick <sup>2</sup>   | 1.48                                            | N.D  | N.D  | 4.39 | N.D  |
| Foundation <sup>1</sup> | N.D                                             | 4.12 | N.D  | N.D  | 3.23 |
| Foundation <sup>2</sup> | 2.82                                            | 4.34 | N.D  | 4.04 | N.D  |
| Home made <sup>1</sup>  | 2.43                                            | 3.36 | 4.03 | 4.16 | 7.08 |
| Home made <sup>2</sup>  | 1.30                                            | 1.68 | 1.26 | 1.54 | 5.70 |

<sup>a</sup> N.D. not detected

### 5. Acknowledgements

This work was supported by the Research Fund for the Doctoral Program of Higher Education (No. 20110002110052).

### References

- [1] van der Rhee H J, de Vries E, Coebergh J W W., *Eur J Cancer*, 2006, 42: 2222-2232
- [2] Armstrong B K, Kricker A., *J Photoch Photobio B*, 2001, 63: 8-18
- [3] Serpone N, Dondi D, Albini A. *Inorg Chim Acta*, 2007, 360: 794-802

- [4] Chatelain E, Gabard B. *Photochem Photobiol*, 2001, 74: 401-406
- [5] Sayre R M, Dowdy J C, Gerwig A J, Shields W J, Lloyd R. V. *Photochem Photobiol*, 2005, 81 (2): 452-456
- [6] Imanaka N, Masui T, Hirai H, Adachi G. *Chem Mater*, 2003, 15: 2289-2291
- [7] Norval M, Cullen A P, de Gruijl F R, Longstreth J, Takizawa Y, Lucas R M, Noonan F P, van der Leun J C. *Photoch Photobio Sci*, 2007, 6: 232-251
- [8] Diaz-Cruz M S, Barcelo D. *TrAC-Trend Anal Chem*, 2009, 28(6): 708-717
- [9] Diaz-Cruz M S, Llorca M, Barcelo D. *TrAC-Trend Anal Chem*, 2008, 27 (10): 873-887
- [10] Giokas D L, Salvador A, Chisvert A. *TrAC-Trend Anal Chem*, 2007, 26: 360-374
- [11] European Directive 76/768/EEC and its successive amendments, basic act 31976L0768
- [12] Japanese, SCI, Japanese Standard of Cosmetic Ingredients, Yakuji Nippo Ltd., Tokyo, 1985
- [13] FDA, Department of Health and Human Services, 21CFR Parts 310, 325, 700 and 740, RIN 0910-AA01, Sunscreen Drug Products for over-the-counter Human Use Final Monograph, Federal Register, Rules and Regulations, 1999
- [14] Salvador A, Chisvert A. *Anal Chim Acta*, 2005, 537: 1-14
- [15] Salvador A, Chisvert A. *Anal Chim Acta*, 2005, 537: 15-24
- [16] Tarazona I, Chisvert A, Leon Z, Salvador A. *J Chromatogr A*, 2010, 1217: 4771-4778
- [17] Rodil R, Schrader S, Moeder M. *J Chromatogr A*, 2009, 1216: 4887-4894
- [18] Oliveira H M, Segundo M A, Lima J L F C, Miro M, Cerda V. *J Chromatogr A*, 2010, 1217: 3575-3582
- [19] Giokas D L, Sakkas V A, Albanis T A. *J Chromatogr A*, 2004, 1026: 289-293
- [20] Negreira N, Rodriguez I, Ramil M, Rubi E, Cela R. *Anal Chim Acta*, 2009, 654: 162-170
- [21] Giokas D L, Sakkas V A, Albanis T A. *J Chromatogr A*, 2005, 1077: 19-27
- [22] Scalia, S. *J Chromatogr A*, 2000, 870: 199-205
- [23] Shih Y, Cheng F C. *J Chromatogr A*, 2000, 876: 243-246
- [24] Vidal L, Chisvert A, Canals A, Salvador A. *Talanta*, 2010, 81: 549-555
- [25] Gaspar L R, Goncalves P M B. *Talanta*, 2010, 82: 1490-1494
- [26] Schakel D J, Kalsbeek D, Boer K. *J Chromatogr A*, 2004, 1049: 127-130
- [27] Smyrniotakis C G, Archontaki H A. *J Chromatogr A*, 2004, 1031: 319-324
- [28] Moeder M, Schrader S, Winkler U, Rodil R. *J Chromatogr A*, 2010, 1217: 2925-2932
- [29] Negreira N, Rodriguez I, Rubi E, Cela R. *J Chromatogr A*, 2009, 1216: 5895-5902
- [30] Rodil R, Moeder M. *J Chromatogr A*, 2008, 1179: 81-88

- [31] Rastogi S C, Jensen G H. *J Chromatogr A*, 1998, 828: 311-316
- [32] Negreira N, Rodriguez I, Ramil M, Rubi E, Cela R. *Anal Chim Acta*, 2009, 638: 36-44
- [33] Sarafraz-Yazdi A, Amiri A. *TrAC-Trend Anal Chem*, 2010, 29(1): 1-14
- [34] Hylton K, Mitra S. *J Chromatogr A*, 2007, 1152: 199-214
- [35] Zhang J, Su T, Lee H K. *J Chromatogr A*, 2006, 1121: 10-15
- [36] Ouyang G, Pawliszyn J. *Anal Chem*, 2006, 78: 5783-5788

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Quantitative analysis of the rat brain dopamine using mass spectrometry.

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
Fig. 1 BF-STEM image of the single-crystal silicon  
(a) Before Cs reduction, (b) After Cs reduction  
Specimen: Single-crystal Si<011>,  
Acceleration voltage: 30 kV


Fig. 2 BF-STEM image around the gate oxide  
Specimen: 3X nm NAND flash memory  
Specimen thickness: 30 nm,  
Acceleration voltage: 30 kV

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井上 陽子(日立ハイテクノロジーズ)他: Amino acid analysis of peptide based on the regulation method

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宝泉 雄介(日立ハイテクノロジーズ)他: Analyses of Biopharmaceuticals using Amino Acid Analyzer

【Summary】Quantitative analysis of amino acids is required in several fields, including biomedical research, bio-engineering, amino acid metabolism and food science. Moore and Stein developed the first automated amino acid analyzer (AAA), combining cation-exchange chromatographic separation of amino acids with postcolumn ninhydrin detection<sup>[1]</sup>. This method is excellent in stability and repeatability. Recently, the production of biopharmaceuticals, such as immunoglobulins

and erythropoietin, has increased greatly. As the majority of these biopharmaceuticals are proteins produced by the cell culture, the control of cell culture conditions and confirmation of amino acid components of product, is very important. In this presentation, we will present several AAA applications, amino acid monitoring of Cell culture medium, Amino acid components analysis of monoclonal antibody often used as a model sample for biopharmaceutical

片山 裕美子(日立ハイテクノロジーズ)他: Food Safety and Environmental Safety for Human Health (Introduction for Hitachi ZA3000 series Polarized Zeeman ASS)

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