

Food Analysis for U-5100 and UH5300

U-5100 Retio beam spectrophotometer/UH5300 Double -Beam Spectrophotometer

Feature



- Simple operation and high-throughput measurement are achieved by a 6-cell turret as standard.
- A wide range of accessories is available to support various application such as short and long path length.
- Environmental friendly Xe flash lamp is a low-power consumption, long-life light source※1.

※1: It is assumed ten years under the measurement conditions of 1 wavelength, 300 measurements / days, 240 days / year (U-5100)



U-5100





Measurement image



UH5300

Measurement



- It is possible to measure absorbance of target component with easy step when measurement kits are used.※2

- L-ascorbic acid in drink
- Nitrite ion in ham
- Ethanol in soy sauce

※2: Please contact the measurement kit maker for the effects of interference substances. For more information, please refer to the our HP "S.I.navi."