

# Tabletop Microscope® TM4000Plus/TM4000


Tabletop Microscope has developed to the next stage.  
The All NEW Tabletop Microscope® TM4000 debuts.



- Support obtaining data to making reports with simple operation
- Newly developed optical system
- Environmentally-friendly pumping system allows you to install almost anywhere.



**Start observing images in just 3 minutes.**  
**Enables you to obtain the data and make reports quickly.**




Camera image helps you find the region of interest and support observation with MAP function.




Newly developed optical system enables enhanced image quality

TM3030Plus



5 μm

TM4000



5 μm

Sample : cement  
Image signal : Secondary electrons  
Accelerating voltage 5 kV