

HIRANUMA Automatic Coulometric Karl Fischer Titrator AQ-2200A Series

HIRANUMA Automatic Volumetric Karl Fischer Titrator AQV-2200A Series

The HIRANUMA Karl Fischer Titrator 2200A series are easy-to-use instruments with touch-panel controls and clearly legible displays. They include an extensive selection of titration cells and vaporizers and facilitate highly efficient analyses.

Main features

(1) Simultaneous measurements


The addition of the unit allows 2 types of water-content measurements simultaneously.

(2) Settings to restrict user privileges

The ability to restrict operation of the instrument to specific users improves security and prevents inadvertent modifications of measurement settings.

(3) Titration cells

By single-chamber electrolysis cells and small-capacity cells that use only a small amount of reagent, these instruments cut down on waste solution and reduce running costs.



The Automatic Coulometric Karl Fischer Titrator AQ-2200A



The Automatic Volumetric Karl Fischer Titrator AQV-2200A