

New! Contact Lens Luminous Transmittance Analyzer

Integrating Sphere Method per ISO 18369-3



- **Based on the High Performance U-3900 UV-Vis Spectrophotometer which provides highly accurate measurements**
- **In-solution measurement of Soft, Hard and Intraocular contact lenses**
- **Specially designed accessories facilitate measurement**
- **60-mm integrating sphere detector for total luminous transmittance measurement**
- **Optional eyeglass measurement accessories for a complete prescription optics evaluation**

To learn more about Hitachi's Contact Lens Measurement System and other Hitachi products, please visit www.hitachi-hta.com/U-3900 or contact us.

Hitachi—The Reliability Leader in Spectroscopy

Hitachi High Technologies America, Inc.
Life Sciences Division
1375 N. 28th Ave
Dallas, TX 75261

Toll Free: (800) 548-9001
E-mail: sales-LS@hitachi-hta.com
Web site: www.hitachi-hta.com/UVspectroscopy