


U-4100 UV-Vis-NIR Spectrophotometer Continuously Variable Absolute Reflectance/Transmittance Accessory



Main features/specifications:

- Sample size: 30x30 mm to 90x90 mm, thickness < 5 mm
- Absolute Variable Angle Reflectance Measurement
- Transmission measurement angle: 0° to 60°
- Reflectance measurement angle: 5° to 70°
- Measurement wavelength range: 340 to 2000 nm
- Ordering information: P/N 1J0-0062

Hitachi—The Reliability Leader in Spectroscopy

Hitachi High Technologies America, Inc.
Life Sciences Division
1375 North 28th Avenue, Dallas, TX 75261

Toll Free: (800) 548-9001
E-mail: sales-LS@hitachi-hta.com
Web site: www.hitachi-hta.com/LShome